TRIGLICÉRIDES MONORREAGENTE

REF **K316**

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FINALIDADE

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in Adults (Adult Treatment Panel III). JAMA 2001; 285:2486-97.

7 - Expert Panel on Integrated Guidelines for cardiovascular H, Risk Reducon in C, Adolescents, Naonal Heart L, Blood I. Expert panel on integrated guidelines for cardiovascular health and risk reducon in children and adolescents: summary report. Pediatrics. 2011; 128 Suppl 5: S213&56.

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


INSTRUCTIONS FOR USE

FUNCTION

Method for determination of Triglycerides in biological samples of serum or plasma (Heparin or EDTA). Enzymatic colorimetric test, only for *in vitro* diagnostic use.

PRINCIPLE OF ACTION

Methodology: Enzymatic Colorimetric (Trinder)

O H₂O₂, 4 Aminoantipyrine and p-Chlorophenol in the presence of Peroxidase result in a cherry colored compound whose color intensity is proportional to the concentration of Triglycerides.

REAGENTS

Number 1 - Enzyme Reagent - Store between 2 and 8°C. Contains: 4-Chlorine Phenol < 5 mmol/L, Lipoprotein Lipase < 5000 U/L, Glycerol Kinase < 3000 U/L, Peroxidase < 5000 U/L, Glycerol-3-Phosphate Oxidase < 5000 U/L, 4-Amineantipyrine < 1 mmol/L, ATP < 5 mmol/L, buffers, activators, stabilizer, surfactant and preservative.

Number 2 - Standard - Store between 2 and 8°C. Contains: Balstilled Glycerin 0,098 g/L (equivalent to 100,0 mg/dL Triglycerides) and preservative.

PRESENTATION

Presentation	Reagent N° 1	Reagent N° 2
K316-2-VET	2 x 100 mL	1 x 3 mL

EQUIPMENTS AND OPERATIONAL INPUTS

Spectrophotometer or colorimeter, water bath at 37°C, watch or stopwatch, pipettes, test tubes, Biocontrol N and Biocontrol P Bioclin. They can be found at markets specialized on Laboratories of Clinical Analysis.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be between 2 to 8°C. The transport at temperatures up to 30°C should not exceed 5 days. Protect from light and avoid moisture. **Do not freeze.**

SPECIAL CARE

1- For *in vitro* diagnostic use only.

2- Strictly follow the methodology proposed to obtain exact results.
3- Water used in material cleaning must be recent and free of contaminants.

4- Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers that may alter the results significantly.

5- Reagents N° 1 and N° 2 contains Sodium Azide, and should be handled with care.

6- The Sodium Hypochlorite is a contaminating agent that can significantly alter the results, so the materials used to perform the tests must be properly washed and free from this type of residue.

7- Water level in water bath should be above the level of reagents in test tubes.

8- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

9- To obtain information related to biosafety or in case of accidents with the product, consult the SDS (Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa

10- Do not use the product in case of damaged packaging.

11- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Serum obtained free of hemolysis (to avoid falsely elevated results) or plasma collected with EDTA or Heparin. The compound is stable for up to 3 days between 2 to 8°C, and 30 days at -20°C. The blood can be collected without fasting or with a fasting period of 12 to 14 hours. Lipemic samples must first be diluted with Sodium Chloride 0,85%, in proportion 1:2.

PROCESS DESCRIPTION

The calibration stability of the Triglycerides Monoreagent kit installed on refrigerated equipment is at least 30 days. This stability may vary depending on the conditions of the test, equipment and environment. Therefore, it is suggested to follow the product performance using control serum.

TECHNIQUE

The Bioclin recommends, as control serum, Biocontrol N and P Bioclin kits.

Mark 3 test tubes: B (Blank), A (Sample), P (Standard), and proceed as follows:

	Blank	Standard	Sample
Sample	---	---	10 µL
Reagent N°1	1.0 mL	1.0 mL	1.0 mL
Reagent N°2	---	10 µL	---

Homogenize well and put in a water bath at 37°C for 10 minutes. Read the absorbance of Sample and Standard at 500 nm (490 - 540 nm), hitting the zero with the Blank. The color is stable for 30 minutes.

CALCULATIONS

$$\text{Triglycerides (mg/dL)} = \frac{\text{Sample Absorbance} \times 100}{\text{Standard Absorbance}}$$

As the reaction follows the Beer-Lambert Law, the Calibration Factor can be used.

$$\text{Calibration Factor} = \frac{\text{Standard Concentration}}{\text{Standard Absorbance}} \quad (100 \text{ mg/dL})$$

$$\text{Triglycerides (mg/dL)} = \text{Sample Absorbance} \times \text{Calibration Factor}$$

Results are expressed in mg/dL.

VETERINARY USE

TO OBTAIN THE INSTRUCTIONS FOR USE IN PRINTED FORMAT, AT NO ADDITIONAL COST, CONTACT CUSTOMER ADVISORY SERVICE:

SAC: +55 (31) 3439 5454 / 0800 031 5454 / sac@bioclin.com.br

INTERFERENCES

Icteric samples with Bilirubin above of 5 mg/dL, interfere with the methodology, leading to falsely diminished. Corticosteroid use and lack of fasting in dogs and cats can cause falsely increased results.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present an analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

TRACEABILITY

The kit's standard is traceable to the reference material NIST (National Institute of Standards and Technology) SRM 1951.

REFERENCE VALUES

The reference values were determined in serum samples of healthy dogs, cats, cattle and horses, males and females, adults and without predilection of breed.

Specie	Reference Value
Dog	22.67 - 103 mg/dL
Cat	20.5 - 112.5 mg/dL
Cattle	6 - 22 mg/dL
Horse	5.045 - 44.19 mg/dL

These values should be used as guidelines and each laboratory should establish its own reference value, according to the population served. The results provided by this kit must be interpreted by the veterinary medical professional, and are not the only criterion for determining the diagnosis and/or treatment of the patient.

To convert the values from mg/dL to mmol/L (SI), multiply by 0.0113.

PRODUCT PERFORMANCE**ACCURACY****Comparison of Methods**

Triglycerides Monoreagent kit was compared with another method for measurement of Triglycerides commercially available. 42 tests were performed and the results were evaluated. The linear equation obtained was $Y = 0.9694X + 3.5955$ and coefficient of correlation 0.9978. With these results we can conclude that the kit shows good methodological specificity.

PRECISION**Repeatability**

The repeatability was calculated from 40 successive determinations, using 3 samples with different concentrations, obtaining the following results:

Repeatability	Sample 1	Sample 2	Sample 3
Average (mg/dL)	102.47	173.31	139.40
SD (mg/dL)	0.63	1.08	0.75
CV (%)	0.61	0.62	0.53

Reproducibility

The reproducibility was calculated from 40 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

Reproducibility	Sample 1	Sample 2	Sample 3
Average (mg/dL)	103.00	176.67	139.87
SD (mg/dL)	0.48	1.36	0.54
CV (%)	0.46	0.78	0.38

SENSITIVITY

The sensitivity was calculated from 40 determinations of a sample free of Triglycerides. The average, that indicates the method detection limit, was 1.34 mg/dL, with standard deviation of 0.41 mg/dL. The sensitivity, that indicates the method quantification limit, corresponds the average plus 3 times the standard deviation, and is equal to 2.58 mg/dL.

LINEARITY

The reaction is linear up to the concentration of 900 mg/dL. For samples with values above 900 mg/dL dilute the sample with Sodium Chloride 0.85%, repeat the dosage and multiply the result by dilution factor.

DIAGNOSTIC SIGNIFICANCE

Triglycerides are important sources of energy for the cell in the form of lipoproteins. Lipemia in samples is frequent and due to the accumulation of chylmicrons in samples from dogs and cats without fasting. The same can occur with the use of corticosteroids due to insulin inhibition. Inflammatory processes of infectious origin, hyperadrenocorticism, pancreatitis and negative energy balance (such as the hyperglycemic syndrome in horses) can lead to an increase in serum triglycerides. The reduction of serum triglycerides has no clinical relevance.

BIBLIOGRAPHIC REFERENCES

- 1 - BUCOLO, G.; DAVID, H., Clin. Chem., 1973, 476.
- 2 - TONKS, D. B., Quality Control in Clinical Laboratories, 1983.
- 3 - MC GOWAN, M. W.; ARTISS, J. D.; STRANDBERGH, D. R.; ZAK, B., Clin. Chem., 1983, 29-538.
- 4 - TRINDER, P.; Ann. Clin. Biochem, 1969, 6-24.
- 5- Reunião Conjunta – Laudos Laboratoriais. SBC/DA, SBAC, SBPC/ML, SBBM; 15 de Outubro de 2013.
- 6- Executive Summary of the Third report of the National Cholesterol Education Program (NCEP) Expert Panel on Detection, Evaluation and Treatment of High Blood Cholesterol

in Adults (Adult Treatment Panel III). JAMA 2001: 285:2486-97.

7 - Expert Panel on Integrated Guidelines for cardiovascular H, Risk Reducon in C, Adolescents, Naonal Heart L, Blood I. Expert panel on integrated guidelines for cardiovascular health and risk reducon in children and adoloscents: summary report. Pediatrics. 2011; 128 Suppl 5: S213&56.

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QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

**QUIBASA QUÍMICA BÁSICA Ltda**

Rua Teles de Menezes, 92 - Santa Branca

CEP 31565-130 - Belo Horizonte - MG - Brasil

Phone: +55 31 3439.5454 | E-mail: bioclin@bioclin.com.br

CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service

Phone: 0800 0315454 | E-mail: sac@bioclin.com.br

Registration Number: Product exempt from registration with the Ministry of Agriculture, Livestock and Supply.

Technical Responsible: Dr. Camila Eckstein (CRMV/MG 20.611)

Review: June/2025

UNIVERSAL SYMBOLOGY

CATALOG NUMBER

LOT NUMBER

MANUFACTURING DATE

VALIDITY DATE
(last day of the month)

TEMPERATURE LIMIT
(store)

CONTENT IS SUFFICIENT
FOR <N> TEST

SEE INSTRUCTIONS
FOR USE

IN VITRO DIAGNOSTIC
PRODUCT

KEEP AWAY
FROM SUNLIGHT

DO NOT REUSE

CAUTION

MADE BY

CONTROL

POSITIVE CONTROL

NEGATIVE CONTROL

BIOLOGICAL RISK

FLAMMABLE

CORROSIVE

TOXIC

DO NOT USE IF
PACKAGE IS
DAMAGED

PRODUCT
STERILIZED

DANGER