**BIO GEÑE**  
Bioclin · QUIBASA

**HPV ALTO RISCO  
MULTIPLEX PCR**

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1. CLINICAL AND LABORATORY STANDARDS INSTITUTE (CLSI). Collection, transport, preparation and storage of specimens for molecular methods; approved guideline. CLSI document MM13-A. Pennsylvania, USA: Clinical and Laboratory Standards Institute, 2005.
2. Schmitz, M; Scheungraber, C; Herrmann, J. et al. Quantitative multiplex PCR assay for the detection of the seven clinically most relevant high-risk HPV types. J Clin Virol. 2009 Apr;44(4):302-7. Epub 2009 Feb 14. PubMed PMID: 19223232.
3. Seaman, W.T; Andrews, E; Couch, M. Detection and quantitation of HPV in genital and oral tissues and fluids by real time PCR. Virol J. 2010 Aug 19; 7:194.

4. Seth, R; Rippin, J; Guo, L; Jenkins, D. Detection and quantitation of HPV gene expression using real-time PCR. Methods Mol Med. 2005; 119:61-72.

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**Bioclin · QUIBASA**

## TABLE OF CONTENTS

Intended Use .....	34
Principle of Action .....	34
Presentation .....	34
Reagents .....	35
Equipment and Operational Supplies .....	35
Storage and Transport Conditions .....	36
Special Care .....	36
Samples .....	37
Procedure .....	38
A . DNA Extraction .....	38
B . Reagent Preparation .....	38
C . PCR Preparation .....	38
D . Thermocycler Settings for Real-Time PCR .....	39
E . Result Validation .....	40
F . Result Interpretation .....	41
Process Limitations .....	42
Product Performance/Quality Control .....	42
Method Comparison and Methodological Specificity .....	42
Repeatability .....	42
Reproducibility .....	43
Clinical Sensitivity .....	43
Analytical Sensitivity .....	44
Linearity .....	44
Cross-Reactivity .....	44
Diagnostic Significance .....	44
References .....	45
Customer Service .....	45
Universal Symbology .....	46

## INTENDED USE

Qualitative detection of high-risk Human Papillomavirus (HPV) DNA genotypes HPV 16, HPV 18, and HPVs (31, 33, 35, 39, 45, 51, 52, 56, 58, 59, 66, and 68) in biological samples through real-time polymerase chain reaction (PCR).

## PRINCIPLE OF ACTION

The **Bio Gene HPV Alto Risco Multiplex PCR** kit is an in vitro assay based on the qualitative detection of high-risk oncogenic HPV genotype DNA through real-time PCR. The kit discriminates HPV genotypes 16 and 18, as well as detects HPVs genotypes (31, 33, 35, 39, 45, 51, 52, 56, 58, 59, 66, and 68).

The real-time PCR method is used to amplify the pathogen's DNA. A real-time PCR thermal cycler is used to amplify and detect specific fluorescent probes for each target. The instrument's software detects the presence of the pathogen.

## PRESENTATION

Reagent	Presentation
	Bio Gene HPV Alto Risco Multiplex PCR
	100 Tests
<b>R1</b>	1 x 110 µL
<b>R2</b>	1 x 1 mL
<b>R3</b>	1 x 500 µL*
<b>R4</b>	1 x 1,5 mL
<b>R5</b>	1 x 1,5 mL

\*Lyophilized reagent. The volume described above corresponds to the final volume after resuspension of the reagent, as described in the PROCEDURE section, subsection B (Reagent Preparation).

## REAGENTS

- R1. PCR Solution:** Primers, probes, TRIS-HCl.  
**R2. Taq Mix:** Polymerase, dNTPs, MgCl<sub>2</sub>, stabilizers.  
**R3. Positive Control:** Plasmid, TRIS-HCl, EDTA.  
**R4. Diluent:** TRIS-HCl, EDTA.  
**R5. Negative Control:** DNase/RNase-free water.

\* R1 contains primers and probes for high-risk HPV genotypes and for the endogenous control.

The kit may include Sample Collection Accessories in the following configurations:

### 1<sup>st</sup> Bio Gene Collection Kit:

- 1) Brush: for endocervical sample collection
- 2) Spatula: for ectocervical sample collection
- 3) Tube with fixative/transport solution
- 4) Vaginal speculum: instrument designed to separate the walls of the vaginal canal to allow gynecological examination

### 2<sup>nd</sup> Bio Gene Collection Kit:

- 1) Brush: for endocervical sample collection
- 2) Vaginal speculum: instrument designed to separate the walls of the vaginal canal to allow gynecological examination
- 3) Tube with fixative/transport solution

### 3<sup>rd</sup> Bio Gene Collection Kit:

- 1) Brush: for endocervical sample collection
- 2) Tube with fixative/transport solution

## EQUIPMENT AND OPERATIONAL SUPPLIES

### Materials included in the kit:

- Reagents listed in the Reagents section.
- Instruction manual.

### Materials needed but not included in the kit:

- 1- Programmable optical fluorescence detection system (real-time PCR thermal cycler);
- 2- Laminar flow hood (for PCR);
- 3- Biological safety cabinet;
- 4- Centrifuge;
- 5- Vortex mixer;
- 6- Nucleic acid extraction kit;
- 7- Micropipettes (0.5–10 µL, 100–1000 µL);

- 8-** Sterile, DNase/RNase-free filter pipette tips;
- 9-** 0.5 mL and 1.5 mL centrifuge tubes;
- 10-** PCR tubes or plates;
- 11-** Optical adhesives;
- 12-** Sample storage tubes;
- 13-** Tube and microtube racks;
- 14-** Personal protective equipment (lab coat, disposable mask, cap, powder-free disposable gloves, beard covers, safety glasses).

### **STORAGE AND TRANSPORT CONDITIONS**

The storage temperature is -20 °C. Transportation can be carried out between 2 and 8 °C for up to 5 days. Freezing and thawing of reagents should be avoided to prevent degradation. Keep the product away from light and protect it from moisture.

### **SPECIAL CARE**

- 1-** This product is for in vitro diagnostic use only.
- 2-** Strictly follow the proposed methodology to obtain accurate results.
- 3-** Handle and dispose of all biological samples, reagents, and materials used in the assay as if they were capable of transmitting infectious agents. Avoid direct contact with biological samples and reagents. Avoid spills or aerosols. Waste should be handled and disposed of in accordance with appropriate safety measures.
- 4-** Molecular biology procedures - such as nucleic acid extraction, reverse transcription, amplification, and detection - require qualified personnel to avoid the risk of erroneous results, particularly due to degradation of nucleic acids in the samples or contamination by amplification products.
- 5-** Separate areas are required for extraction/reaction preparation and for amplification/detection. Never introduce amplification products into the area designated for extraction or preparation of amplification reactions.
- 6-** All samples and reagents must be handled under a laminar flow hood. Pipettes must be used with filter tips. Tips must be sterile and DNase/RNase-free.
- 7-** Avoid repeatedly freezing and thawing the reagents.
- 8-** Store DNA samples at -20 °C if not used immediately.
- 9-** Do not use the kit after the expiration date.
- 10-** We recommend following local, state, and federal environmental protection regulations to ensure that the disposal of reagents and biological materials is carried out in accordance with applicable laws.

**11-** For information regarding biosafety or in the event of accidents involving the product, consult the Safety Data Sheets (SDS) available at [www.bioclin.com.br](http://www.bioclin.com.br) or upon request from Quibasa Customer Advisory Service (SAC).

**12-** Do not use the product if the packaging is damaged.

**13-** It is essential that all instruments and equipment used are properly calibrated and subjected to regular maintenance.

## **SAMPLES**

- The Bio Gene HPV Alto Risco Multiplex PCR kit must be used with DNA samples extracted from cervical, urethral, and anorectal swabs. Other sample types may be used according to medical or laboratory recommendations, following appropriate validation.

- Samples must be collected by trained professionals in accordance with the laboratory's recommendations for molecular testing.

- Prior to sample collection, it is recommended that the patient follow these instructions:

- Do not use vaginal creams or douches for 3 days before collection.

- Abstain from sexual intercourse for 24 hours before collection.

- Preferably collect the sample outside the menstrual period.

- Do not apply acetic acid (vinegar) and/or iodine solution to the vagina before collection.

- Cervical scraping samples or biopsies must be transported and stored in appropriate media for subsequent extraction and real-time PCR.

- These may include media intended exclusively for molecular diagnostics or media suitable for both cytological and molecular diagnostics, in which the same sample (liquid biopsy) will be used.

- Fresh biopsy samples in liquid medium should be processed within 24 hours if kept at room temperature. If stored frozen at -20 °C, samples should be processed within 30 days. These measures ensure the preservation and integrity of the genetic material obtained from collected samples.

- Avoid fixing cervical scraping samples or biopsies in formalin, as this compromises DNA extraction and real-time PCR performance.

- It is essential to follow the manufacturer's instructions regarding the stability and storage temperature of samples in liquid medium.

- Proper handling ensures accurate and reliable laboratory analysis results. For long-term storage, freeze samples at -70 °C<sup>1</sup>.

- More than five freeze-thaw cycles should be avoided to maintain DNA integrity.

## PROCEDURE

### A. DNA Extraction

Nucleic acids (DNA) from the samples must be extracted following the instructions for use of the selected extraction kit, ensuring that the integrity, quality, and purity of the genetic material are appropriate for PCR. A minimum input volume of 200  $\mu\text{L}$  of biological sample must be used for extraction. Sample volumes below 200  $\mu\text{L}$  may be used provided that the procedure has been previously validated.

### B. Reagent Preparation

Reagent R3 contains a DNA template and must be handled (resuspended) in an appropriate area to prevent contamination of the other reagents.

- 1- Centrifuge (spin down) the **Positive Control (R3)** before opening the microtube.
- 2- Resuspend the **Positive Control (R3)** with 500  $\mu\text{L}$  of **Diluent (R4)**.

*\*Passive reference dye (ROX) must not be added to the **Taq Mix (R2)**, as ROX is used in the detection system for the HPV 18 target.*

### C. PCR Preparation

- 1- Prepare in advance the microtubes/wells to be used, according to the number of samples and controls to be analyzed.
- 2- Prepare the total volume of the final PCR solution according to the number of reactions to be performed.

Reagents	1 Reaction	25 Reactions	50 Reactions	100 Reactions
Taq Mix (R2)	10 $\mu\text{L}$	250 $\mu\text{L}$	500 $\mu\text{L}$	1mL
PCR Solution (R1)	1 $\mu\text{L}$	25 $\mu\text{L}$	50 $\mu\text{L}$	100 $\mu\text{L}$

For preparing a different number of reactions, multiply the reagent volumes for 1 reaction by the required number of reactions.

- 3- Pipette 11  $\mu\text{L}$  of the final PCR solution into the designated tubes or wells.
- 4- Add 9  $\mu\text{L}$  of **extracted sample DNA** ( $\geq 20 \text{ ng}/\mu\text{L}$ ), or 9  $\mu\text{L}$  of the **Positive Control (R3)**, or 9  $\mu\text{L}$  of the **Negative Control (R5)**.
- 5- Homogenize carefully to avoid bubble formation.
- 6- Note that the total reaction volume is 20  $\mu\text{L}$ , and each PCR run must include the appropriate controls (Negative Control and Positive Control).
- 7- Homogenize carefully to avoid bubble formation.
- 8- Transfer the tubes/plate to the thermal cycler.

#### ***D. Thermocycler Settings for Real-Time PCR***

Refer to the real-time PCR instrument's operating manual for experiment programming.

##### **1- Set the experiment type to:**

Qualitative Test.

##### **2- Configure the fluorescent detectors (probes) as follows:**

Target	Detector	Quencher
HPV 16	VIC/HEX	NFQ-MGB
HPV 18	ROX	
HPVs (31, 33, 35, 39, 45, 51, 52, 56, 58, 59, 66, 68)	FAM	
Endogenous Control	Cy5	

#### **NOTE:**

- Instruments that use ROX as a passive reference must be programmed with the passive reference option set to "none."
- Extracted samples must be labeled with the VIC/HEX, ROX, FAM, and Cy5 detectors.

### 3- Set the Cycling Conditions

Step	Temperature	Time	Cycles
1	95°C	3 minutes	1
2	95°C	15 seconds	6
	60°C	60 seconds	
3	95°C	15 seconds	39
	60°C	60 seconds	

Set “Data Collection” to “stage 3, step 2 (60°C - 0:60)”.

### E. Result Validation

#### 1- Controls

Control	Acceptable CT Range	Amplification/ Detection
Positive	VIC/HEX $\leq$ 24 ROX $\leq$ 24 FAM $\leq$ 24	Valid
Negative	Not detected	Valid
Endogenous	Cy5 $\leq$ 35	Valid

**NOTE:** CT values for the Endogenous Control may vary depending on process conditions, such as DNA extraction efficiency, sample concentration, and thermal cycler settings. Samples with a high number of DNA copies may, in some cases, inhibit amplification of the Endogenous Control.

- A low CT value indicates a high concentration of biological material in the sample and that the extraction process was efficient.
- A high CT value indicates a low concentration of biological material in the sample and/or an inefficient extraction.

Therefore, these conditions should be evaluated when CT values fall outside the acceptable range and, when appropriate, the results may still be validated.

If the above requirements are not met, the assay is considered invalid and the test must be repeated.

### ***F. Result Interpretation***

<b>HPV Alto Risco Multiplex</b>				<b>Result</b>
<b>VIC/HEX</b>	<b>ROX</b>	<b>FAM</b>	<b>CY5</b>	
<b>CT DETECTED</b>	CT NOT DETECTED	CT NOT DETECTED	<b>CT DETECTED</b>	<b>Positive for HPV 16</b> Negative for HPV 18 Negative for HPVs (31, 33, 35, 39, 45, 51, 52, 56, 58, 59, 66, 68)
CT NOT DETECTED	<b>CT DETECTED</b>	CT NOT DETECTED	<b>CT DETECTED</b>	<b>Positive for HPV 18</b> Negative for HPV 16 Negative for HPVs (31, 33, 35, 39, 45, 51, 52, 56, 58, 59, 66, 68)
CT NOT DETECTED	CT NOT DETECTED	<b>CT DETECTED</b>	<b>CT DETECTED</b>	<b>Positive for HPVs (31, 33, 35, 39, 45, 51, 52, 56, 58, 59, 66, 68)</b> Negative for HPV 16 and HPV 18
CT NOT DETECTED	CT NOT DETECTED	CT NOT DETECTED	<b>CT DETECTED</b>	<b>Negative for all high-risk HPV genotypes tested</b>

**CT detected:** presence of an amplification curve with a sigmoidal pattern.

**CT not detected:** absence of an amplification curve.

\*The results apply only to the high-risk genotypes detected by this kit, as described in the Intended Use section.

**Note:** Amplification (CT detected) may occur for more than one genotype; in such cases, this indicates a co-infection. The result is valid provided that all controls fall within their acceptable CT ranges and the amplification curve of the tested sample resembles the curves obtained for the Positive Control (R3), showing a sigmoidal pattern.

Failure to detect pathogen DNA does not exclude the possibility of infection when the pathogen load is below the detection limit of this kit. The results provided by this kit must be interpreted by the responsible healthcare professional and should not be used as the sole criterion for determining the patient's diagnosis and/or treatment. Results must be evaluated in conjunction with the patient's clinical data and additional laboratory tests.

### **Process Limitations**

Cross-contamination during sample collection, processing, transport, or storage may lead to false results.

## **PRODUCT PERFORMANCE/QUALITY CONTROL**

### **Accuracy**

#### **METHOD COMPARISON AND METHODOLOGICAL SPECIFICITY**

The **Bio Gene HPV Alto Risco Multiplex PCR** kit was compared with another kit designed for quantifying HPV DNA. A total of 100 analyses were performed, and the results were evaluated. Sample results from both products were compared using agreement between results as the acceptance criterion. Based on these findings, the kit demonstrates good methodological specificity.

### **Precision**

#### **REPEATABILITY**

Ten measurements were performed in a single day using two samples. The product met the established acceptance criteria, demonstrating 100% agreement in the results obtained for both samples.

## REPRODUCIBILITY

Ten measurements were performed over three consecutive days using one sample. The product met the established acceptance criteria, demonstrating 100% agreement in the results obtained for the analyzed sample.

## CLINICAL SENSITIVITY

The clinical sample panel shown below demonstrates that the clinical sensitivity and specificity of the **Bio Gene HPV Alto Risco Multiplex PCR** kit are greater than 99.9%.

Clinical Sensitivity			
Target	True Positive	False Negative	Sensitivity
HPV 16 (VIC/HEX)	20	0	100%
HPV 18 (ROX)	18	0	100%
Other high-risk HPVs subtypes (12 subtypes) (FAM)	12	0	100%

Clinical Specificity			
Target	True Positive	False Negative	Specificity
HPV 16 (VIC/HEX)	27	0	100%
HPV 18 (ROX)	13	0	100%
Other high-risk HPVs subtypes (12 subtypes) (FAM)	10	0	100%

## ANALYTICAL SENSITIVITY

The analytical sensitivity determined in the performance studies was:  
11.11 copies/ $\mu$ L (100 copies/reaction) for HPV 16;  
11.13 copies/ $\mu$ L (100.17 copies/reaction) for HPV 18;  
11.16 copies/ $\mu$ L (100.44 copies/reaction) for the other 12 HPV types.

**NOTE:** The analytical sensitivity of the product may be affected by factors such as the efficiency of the nucleic acid extraction kit/method used and the performance capability of the real-time thermal cycler.

## LINEARITY

The linearity observed was approximately 5.0 log ( $1 \times 10^5$  copies/ $\mu$ L).

## CROSS-REACTIVITY

No cross-reactivity was observed with any pathogens that infect the urinary and/or reproductive systems, nor with sexually transmitted pathogens. The species included in the study were: *Ureaplasma urealyticum*, *Chlamydia trachomatis*, *Mycoplasma hominis*, *Neisseria gonorrhoeae*, HIV, HCV, and HBV.

## DIAGNOSTIC SIGNIFICANCE

Human papillomavirus (HPV), belonging to the Papillomaviridae family, infects epithelial cells of the skin and mucosa, causing lesions that result from abnormal cellular growth. There are more than 150 different HPV subtypes, approximately 35 of which infect the anogenital region of both females and males, potentially causing conditions ranging from warts to low- and high-grade dysplastic lesions.

HPV transmission occurs through direct contact with the infected area; in the case of genital HPV, transmission occurs through sexual contact. HPV types are classified as low-, intermediate-, or high-risk according to the types of lesions they cause.

HPV genotypes 16, 18, 31, 33, 35, 39, 45, 51, 52, 56, 58, 59, 66, and 68 are categorized as high-risk viruses because they may be associated with tumors. Among these, genotypes 16, 18, 31, and 45 account for 85% of cervical cancer cases.

Therefore, early identification of these viruses is essential to prevent the development of malignant lesions. Real-time PCR detection is highly valuable due to its sensitivity, precision, and rapid results.

## REFERENCES

1. CLINICAL AND LABORATORY STANDARDS INSTITUTE (CLSI). Collection, transport, preparation and storage of specimens for molecular methods; approved guideline. CLSI document MM13-A. Pennsylvania, USA: Clinical and Laboratory Standards Institute, 2005.
2. Schmitz, M; Scheungraber, C; Herrmann, J. et al. Quantitative multiplex PCR assay for the detection of the seven clinically most relevant high-risk HPV types. J Clin Virol. 2009 Apr;44(4):302-7. Epub 2009 Feb 14. PubMed PMID: 19223232.
3. Seaman, W.T; Andrews, E; Couch, M. Detection and quantitation of HPV in genital and oral tissues and fluids by real time PCR. Virol J. 2010 Aug 19; 7:194.
4. Seth, R; Rippin, J; Guo, L; Jenkins, D. Detection and quantitation of HPV gene expression using real-time PCR. Methods Mol Med. 2005; 119:61-72.

### **QUIBASA QUÍMICA BÁSICA Ltda**

Rua Teles de Menezes, 92 – Santa Branca  
CEP 31565-130 – Belo Horizonte – MG – Brazil  
Phone: +55 (31) 3439-5454  
Email: [sac@bioclin.com.br](mailto:sac@bioclin.com.br)  
CNPJ: 19.400.787/0001-07 – Brazilian Industry

## **CUSTOMER SERVICE**

Customer Advisory Service  
Phone: 0800 031 5454  
E-mail: [sac@bioclin.com.br](mailto:sac@bioclin.com.br)

Registration number of the BIO GENE HPV ALTO RISCO MULTIPLEX PCR kit with ANVISA: 10269360457.

## UNIVERSAL SYMBOLOGY

	CATALOG NUMBER		MADE BY
	LOT NUMBER		CONTROL
	MANUFACTURING DATE		POSITIVE CONTROL
	VALIDITY DATE (last day of the month)		NEGATIVE CONTROL
	TEMPERATURE LIMIT (store)		BIOLOGICAL RISK
	CONTENT IS SUFFICIENT FOR <N> TEST		FLAMMABLE
	SEE INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC PRODUCT		TOXIC
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGED IS DAMAGED
	DO NOT REUSE		PRODUCT STERILIZED
	CAUTION		DANGER

**Bioclin · QUIBASA**

 Rua Teles de Menezes, 92 . Belo Horizonte . MG . Brasil . CEP: 31565-130  
Tel +55 31 3439 5454 . Fax +55 31 3439 5455 . [www.bioclin.com.br](http://www.bioclin.com.br)  
FARM. RESP. Silvio Wandalsen Arndt - CRF MG 7422  
C.N.P.J.: 19.400.787/0001-07 - Indústria Brasileira