LDL DIRETO

REF K313

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FINALIDADE

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
INSTRUCTIONS FOR USE

FUNCTION

Method for determination of LDL Cholesterol in biological samples of serum or plasma (EDTA or Heparin). Enzymatic colorimetric test. For *in vitro* diagnostic only.

PRINCIPLE OF ACTION

Methodology: Enzymatic colorimetric (end point)

The intensity of the formed color is directly proportional to the LDL Cholesterol concentration in the sample.
CE = Cholesterol Esterase
CO = Cholesterol Oxidase
4 - AA = 4 - Aminoantipyrine
DSBmT = N,N - Bis (4 - Sulfobutyl) - m - Toluidine Disodium

REAGENTS

Number 1 - Enzymatic Reagent - Store between 2 and 8°C. Contains: Buffer, Cholesterol Esterase 1,5 KU/L, Cholesterol Oxidase 1,5 KU/L, Peroxidase 1,3 KU/L, 4 - Aminoantipyrine 0,1%, Ascorbic Acid Oxidase 3 KU/L, detergent 1 and preservative.
Number 2 - Color Reagent - Store between 2 and 8°C. Contains: Buffer, N,N - Bis (4 - Sulfobutyl) - m - Toluidine Disodium (DSBmT) 1,0 mM, detergent 2 and preservative.
Number 3 - Calibrator* - Store between 2 and 8°C. Contains: Lyophilized serum with LDL value set and preservative.
Potentially Infectious.
Warning: The concentration varies with the lot - see bottle label.

PRESENTATION

Reagent Nº 1	Reagent Nº 2	Reagent Nº 3*
1 x 30 mL	1 x 10 mL	1 x 2 mL
2 x 30 mL	2 x 10 mL	1 x 2 mL
4 x 30 mL	4 x 10 mL	1 x 2 mL
2 x 60 mL	2 x 20 mL	1 x 2 mL
1 x 30 mL	1 x 10 mL	1 x 5 mL
2 x 30 mL	2 x 10 mL	1 x 5 mL
4 x 30 mL	4 x 10 mL	1 x 5 mL
2 x 60 mL	2 x 20 mL	1 x 5 mL

*The R3 Calibrator is not present in all presentations.

EQUIPMENTS AND OPERATIONAL INPUTS

Spectrophotometer or colorimeter, water bath at 37°C, watch and stopwatch, pipettes, test tubes, Biocontrol N and Biocontrol P Bioclin. They can be found at market specialized on Clinical Analysis Laboratories.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be between 2 and 8°C The transport at temperatures up to 30°C should not exceed 5 days. Protect from light and avoid moisture. **Do not freeze.**

SPECIAL CARE

- For *in vitro* diagnostic use only.**
- Strictly follow the methodology proposed to obtain exact results.
- Water used in material cleaning must be recent and free of contaminants.
- Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers that may alter the results significantly.
- To avoid possible analytical errors, it is advisable that the distilled water used for reconstitution of the calibrator is of adequate quality. For this, it is recommended that water used in clinical laboratories follow the specifications of the Clinical Laboratory Standards Institute (CLSI).
- The water level in water bath must be superior to the reagents level in the test tubes.
- Handle the reagents carefully because they contain Sodium Azide, irritating to skin and membrane mucous.
- Determine the factor periodically and in each product batch.
- Protect the reagent from direct light.
- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.
- To obtain information related to biosafety or in case of accidents with the product, consult the SDS (Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.
- Do not use the product in case of damaged packaging.
- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLE

Serum obtained free of hemolysis or plasma harvested with Heparin or EDTA. LDL cholesterol is stable in serum for 5 days if stored at 2 and 8°C and 3 months at -20°C ¹¹.
Note: Samples can only be frozen once.
Anticoagulants containing Citrate can not be used.

PROCESS DESCRIPTION

The calibration stability of the LDL Direct kit installed on refrigerated equipment is at least 12 days. This stability may vary depending on the conditions of the test, equipment and environment. Therefore, it is suggested to follow the product performance using control serum.

WORKING REAGENT PREPARATION

Before reconstitute the Calibrator should be noted the volume indicated on the bottle and proceed as described below.
Presentation 1, 2, 4 and 9: Open carefully the bottle and add exactly **2 mL of distilled water**, as volume described in the label of the bottle.
Presentation 27, 28, 30 and 35: Open carefully the bottle and add exactly **5 mL of distilled water**, as volume described in the label of the bottle.
Mix. Leave the bottle resting for 30 minutes.

TO OBTAIN THE INSTRUCTIONS FOR USE IN PRINTED FORMAT, AT NO ADDITIONAL COST, CONTACT CUSTOMER ADVISORY SERVICE:
SAC: +55 (31) 3439 5454 / 0800 031 5454 / sac@bioclin.com.br

Every 10 minutes, agitate the bottle carefully using circular motion to prevent material sticking to the wall of the bottle. Avoid formation of foam.
After reconstitution, the calibrator is stable for one week at 2 to 8°C or for 30 days at -20°C.

TECHNIQUE

Bioclin recommends, as control serum, Biocontrol N and P Bioclin kits.

Mark 3 test tubes with letters B (Blank), A (Sample), C (Calibrator) and proceed as follows:

	Blank	Calibrator	Sample
Reagent 1	300 µL	300 µL	300 µL
Distilled Water	3 µL	---	---
Calibrator	---	3 µL	---
Sample	---	---	3 µL

Agitate and after 5 minutes incubation at 37°C, measure the absorbance at 546 nm (540 - 600 nm).

Reagent 2	100 µL
-----------	--------

Agitate and after 5 minutes of incubation at 37°C, measure the absorbance in 546 nm (540 - 600 nm).

CALCULATIONS

$$\frac{(A_2 - A_1)_{\text{Sample}}}{(A_2 - A_1)_{\text{Standard}}} \times n$$

n = Calibrator Concentration

INTERFERENTS

The essay equilibrium is affected by a few substances such as Bilirubin above 20 mg/dL, Hemoglobin above 500 mg/dL, Ascorbic Acid above 50 mg/dL and Gamma Globulins above 5000 mg/dL.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

TRACEABILITY

The kit's calibrator is traceable to the reference material NIST (National Institute of Standards and Technology) SRM 1951.

REFERENCE VALUES ⁷

The reference values in mg/dL, for this method, were obtained through the determination of LDL Cholesterol in healthy populations of male and female.
Reference values below are for individuals with 12 hours fasting and follow Adult Treatment Panel III (ATP III).

Adults (≥ 20 years)	Optimum	< 100 mg/dL
	Desirable	100 - 129 mg/dL
	Limit	130 - 159 mg/dL
	Elevated	160 - 189 mg/dL
	Very Elevated	≥ 190 mg/dL

Childrens and Adolescents (2 to 19 years)	Desirable	<100 mg/dL
	Limit	100 - 129 mg/dL
	Elevated	≥ 130 mg/dL

To convert the mg/dL values to mmol/L (SI) multiply by 0.026. These values should be used as guidelines, each laboratory should establish its range of reference values, according to the population served. The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

PRODUCT PERFORMANCE
ACCURACY

Comparison of Methods
The LDL Direct kit was compared with other methods for LDL Cholesterol dosage commercially available, 42 analyses were performed and the results were evaluated. The linear equation obtained was Y = 1.0265X - 1.1202 and the correlation coefficient 0.9987. With these results, we can conclude that the kit shows good methodological specificity.

PRECISION
Repeatability

The repeatability was calculated from 40 successive determinations, using 3 samples with different concentrations, obtaining the following results:

Repeatability	Sample 1	Sample 2	Sample 3
Average Concentration (mg/dL)	13.28	43.30	28.98
Standard Deviation (mg/dL)	0.45	0.61	0.70
Coefficient of Variation (%)	3.41	1.40	2.41

Reproducibility

The reproducibility was calculated from 40 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

Reproducibility	Sample 1	Sample 2	Sample 3
Average Concentration (mg/dL)	13.82	43.43	28.87
Standard Deviation (mg/dL)	0.51	0.62	0.55
Coefficient of Variation (%)	3.66	1.43	1.89

SENSITIVITY

The sensitivity was calculated from 40 determinations of a sample free of LDL Cholesterol. The average found was of 0.234 mg/dL with standard deviation of 0.049 mg/dL. The sensitivity, which indicates the method detection limit, corresponds the average plus 3 times the standard deviation and is equal to 0.382 mg/dL.

LINEARITY

The test is linear in a concentration track that varies from 5.6 to 998 mg/dL.

SPECIFICITY

The specificity of the kit is guaranteed by the selective action of reagent N° 1 detergent that allows the consumption of the cholesterol presented only in non LDL particles.

DIAGNOSTIC SIGNIFICANCE

Cholesterol, an insoluble molecule, circles associated to lipoproteins (HDL, LDL and VLDL). LDL (Low Density Lipoprotein) is obtained from VLDL (Very Low Density Lipoprotein) through the hydrolysis by different lipolytic enzymes. LDL, who transports approximately 60% of total plasmatic cholesterol, are in its majority absorbed by specific receptors of extra hepatic and hepatic tissues. There is a confirmed relation between cardiac coronary diseases incidences and the LDL Cholesterol. LDL are atherogenic lipoproteins: the increase of LDL Cholesterol is one of the main causes of atherosclerosis apparition and evolution, mainly the coronary. Therefore, the elevated LDL Cholesterol treatment is the primary target of the cholesterol reduction therapy. An increase in LDL Cholesterol can be seen in different pathology states, inclusive hyperpolyproteinemia due to hepatic or kidney dysfunctions, hypothyroidism and diabetes.

BIBLIOGRAPHIC REFERENCES

- 1 - RIFAI, N., BACHORIK, P.S., ALBERS, J.J., Tietz Fundamentals of Clinical Chemistry, 5ª Ed., (2001), 463.
- 2 - NAITO, H.K., Clinical Chemistry: Concepts and Application, (2003), 603.
- 3 - TIETZ, N.W., Clinical Guide to Laboratory Tests, 3ª Ed., (1995), 404.
- 4 - VASSAULT, A., et al, Ann. Biol. Clin., (1986), 44, 686.
- 5 - VASSAULT, A., et al, Ann. Biol. Clin., (1999), 57, 685.
- 6 - YOUNG, D.S., Effects of Preanalytical Variables on Clinical Laboratory Tests, 2ª Ed., (1997).
- 7 - Executive Summary of the Third report of the National Cholesterol Education Program (NCEP) Expert Panel on Detection, Evaluation and Treatment of High Blood Cholesterol in Adults (Adult Treatment Panel III). JAMA 2001; 285:2486-97.
- 8 - YOUNG, D.S., Effects of Drugs on Clinical Laboratory Tests, 4ª Ed., (1995).
- 9 - ESTEBAN-SALAN, M., et al, Clin. Chem., (2000), 46, 1121.
- 10 - Reunião Conjunta – Laudos Laboratoriais. SBC/DA, SBAC, SBPC/ML, SBBM; 15 de Outubro de 2013.
- 11 - WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.
- 12 - QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda
Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone: +55 31 3439.5454 | E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service
Phone: 0800 0315454 | E-mail: sac@bioclin.com.br

ANVISA registration for LDL Direct Family kit: 10269360372

Review: March/2025(1)

REF

CATALOG NUMBER

LOT

LOT NUMBER

MANUFACTURING DATE

VALIDITY DATE
(last day of the month)

TEMPERATURE LIMIT
(store)

CONTENT IS SUFFICIENT
FOR <N> TEST

SEE INSTRUCTIONS
FOR USE

IN VITRO DIAGNOSTIC
PRODUCT

KEEP AWAY
FROM SUNLIGHT

DO NOT REUSE

CAUTION

MADE BY

CONTROL

POSITIVE CONTROL

NEGATIVE CONTROL

BIOLOGICAL RISK

FLAMMABLE

CORROSIVE

TOXIC

DO NOT USE IF
PACKAGE IS
DAMAGED

PRODUCT
STERILIZED

DANGER