LEISHMANIOSE IgG VET FAST

REF VET004

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FINALIDADE

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INSTRUCTIONS FOR USE

FUNCTION

Fast immunochromatographic test for the qualitative determination of anti-*Leishmania* antibodies in serum, plasma or whole blood canine samples.

PRINCIPLE OF ACTION

Methodology: Immunochromatography.

Leishmaniasis IgG VET FAST is an *in vitro* diagnostic test for the rapid and qualitative detection of anti-*Leishmania* antibodies in serum, plasma or whole blood samples. The cassette consists of a nitrocellulose membrane impregnated with the *Leishmania* recombinant antigen in test line. The conjugate consists of an antibody anti-dog IgG conjugated to colloidal gold particles. When the sample is added to the cassette, it migrates by capillary action until it reaches the conjugate. In this region, the anti-*Leishmania* IgG antibodies present in the sample interact with the conjugate, forming the complex: dog IgG – conjugate. The sample continues to migrate through the immunochromatographic strip until it reaches the test line, where the dog IgG antibody also interacts with the immobilized antigen of *Leishmania* in the test line, forming the complex: *Leishmania* antigen – dog IgG – conjugate. This shows the formation of a pink-red line on the test line; indicating positive result.

Sample excess continues to migrate through the test strip until it reaches the control line, which is composed of impregnated anti-conjugate antibody. Thus, the formation of a second line (control line) is observed; indicating that the test yielded a valid result.

REAGENTS

Number 1 (R1): Cassette - Store between 15 and 30°C. **Do not freeze.** It contains:

- A) 01 plastic cassette with 01 test strip containing
 - Control region (C): Antibody anti-conjugated.
 - Test region (T): Immobilized *Leishmania* recombinant antigen.
 - Conjugate: Antibody anti-dog IgG, colloidal gold conjugate.

B) 01 sachet of silica

Number 2 (R2): Diluent - Store between 15 and 30°C. **Do not freeze.** Contains: Buffer solution and preservative.

Collection Accessories:

- Disposable Plastic Pipette
- Collection tube with EDTA

PRESENTATIONS

Reagent	Presentation	Number of Cassettes for Package	Diluent
Cassette + Diluent	VET004-1	10	1 x 2.5 mL
	VET004-2	25	1 x 5.0 mL
	VET004-3*	10	10 x 0.5 mL
Diluent	VET004-4*	25	25 x 0.5 mL
	VET004-5	5	1 X 1.25 mL

* Individualized presentations.

EQUIPMENTS AND OPERATIONAL INPUTS

Clock or stopwatch. They are in the specialized market of articles for Clinical Analysis Laboratories.

TRANSPORTATION AND STORAGE CONDITIONS

The storage and transport temperature of the product must be between 15 and 30°C. Keep away from direct sunlight.

Do not freeze.

The cassette is sensitive to air humidity. Therefore, it should always be kept in its inviolate primary packaging (aluminized sachet) until the time of use.

SPECIAL CARE

1- Product for *in vitro* diagnosis, for veterinary use only.

2- - Strictly follow the methodology proposed for obtaining results.

3- Cassettes (test devices) do not contain contaminants; but since the samples are biological, the professional must strictly follow the safety rules and routines. The use of disposable gloves and other personal protective equipment is essential.

4- Use clean, dry containers to package biological samples.

5- Open the aluminum envelop to remove the cassette only at the time of use.

6- The kit cannot be frozen or exposed to high temperatures, since this will cause deterioration of its components.

7- The reagents, as well as the samples, should reach room temperature before starting the test.

8- Do not use reagents from different lots or kits.

9- Do not use kit beyond its expiration date. The expiry date is stated on the product label.

10- The cassette (test device) and the plastic pipette are disposable and cannot be reused.

11- The EDTA collector tube accessory must be used exclusively with this product.

12- Blood samples with clots, microclots and formation of fibrin mesh can cause clogging of the filter mesh, making it impossible for the sample to pass correctly through the test membrane.

13- The kit materials that had contact with the patient's sample must have a special destination. Plastic pipettes, cassettes (test devices) and diluent bottles must be discarded in biological waste.

14- Secondary packaging (box), Instruction for use, aluminum envelope and silica sachet can be disposed of in trash common.

15- We recommend to apply the local, state and federal environmental protection so that the disposal of reagents and biological material is done in accordance with the current legislation.

16- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or through request by the Customer Service Department (SAC) of the Bioclin.

SAMPLES

Perform the test with serum, plasma or whole blood.

Serum: Allow blood to clot and then centrifuge to separate serum.

Plasma: Collect whole blood in a tube containing anticoagulants such as heparin, citrate, or EDTA. Centrifuge the blood and separate the plasma.

VETERINARY USE

Whole blood: Collect the blood sample in a tube containing EDTA, heparin, or citrate using standard laboratory procedures.

Whole blood samples should be analyzed within 3 days, can be stored at 2-8°C and never be frozen.

Serum and plasma samples should be analyzed immediately after collection, but if necessary, they can be stored at 2-8°C for up to 7 days or at -20°C for longer periods.

PROCESS DESCRIPTION

Stability After Opening

The stability of the Leishmaniasis IgG VET FAST cassettes (test devices), after opening the aluminized sachet, is 2 hours. This stability may vary according to environmental conditions after opening the aluminized sachet.

Technique

1- The sample should be at room temperature before starting the test.

2- Remove the cassette (test device) from the protective packaging and place it on a clean, dry and level surface

3- Collect the sample using the disposable pipette that accompanies the kit.


4- Press the pipette bulb and touch it to the sample. Relieve the pressure on the bulb so that the sample can be aspirated. Aspirate up to the trace to obtain a volume of 10 μ L.

5- To dispense the sample, press the bulb again, touching the lower part of the pipette to the test membrane in the sample application well.

Note: Alternatively, an automatic pipette can be used to collect and add the sample.

6- Hold the Diluent bottle (Reagent N° 2) vertically and apply 2 drops (70 μ L) in sample well.

7- Wait for the lines to form. Interpret the results between 15 and 20 minutes. Do not interpret after 20 minutes.

INTERPRETATION OF RESULTS

Reactive Test: Formation of two lines: one in the control region (C) and one in the test region (T). Do not interpret after 20 minutes.	
Not Reactive Test: Line formation only in control region (C). Do not interpret after 20 minutes.	
Invalid Test: The absence of line formation in the control region (C), with or without line formation in the test region (T), indicates procedure error or deterioration of the cassette. In this case, repeat the test using a new cassette.	 

According to the Ministry of Health (2011), dogs with a reactive result in the screening carried out by a rapid test (immunochromatographic) must be sent for confirmatory examination by immunoenzymatic assay (ELISA).

TEST LIMITATIONS

This product is intended for screening for Leishmaniasis in dogs. Interpretation should not be based on a single essay. Other confirmatory tests must be included before a sample is considered positive. Therefore, for the diagnosis of the disease, the professional must take into account not only the result of this test; but also, clinical signs, animal health conditions, and other diagnostic results that may be available. Serological tools are not recommended to monitor the effectiveness of dogs undergoing treatment, considering that the amount of antibodies produced by the animal's body is not proportional to the parasite load. For this purpose, it is recommended to measure the parasite load using molecular tests in bone marrow samples or lymph node aspirates.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal program of quality control, where procedures, standards, limits and tolerance for variations are clearly established. It is important to note that all measuring systems characteristic analytical variability, which should be monitored by the laboratories themselves. Therefore, it is the use of controls, which allow the assessment of the precision and accuracy of dosages.

PRODUCT PERFORMANCE

CLINICAL SENSITIVITY AND SPECIFICITY

This product has been tested comparing with other methods. A total of 105 samples characterized by canine Leishmaniasis were analyzed; 51 samples were diagnosed as positive and 54 were diagnosed as negative. The results show that the clinical sensitivity of the product is >99% and the clinical specificity is 96.30%.

METHOD		REFERENCE		Total
Leishmaniasis IgG VET FAST	Result	Reactive	Not reactive	
	Reactive	51	2	53
	Not reactive	0	52	52
Total result		51	54	105

Clinical sensitivity:(51/51): > 99% (IC 95% 95% - 100%)
Clinical specificity: (52/54): 96,30% (IC 95% 91,35% - 100%)
Precision: (103/105): 98.09%

REPEATABILITY

Repeatability was calculated from 10 successive determinations, using 2 different samples, obtaining the following results:

Samples	Nº of Repetitions	Expected Result	Found Result
Sample 1	10	Reactive	Reactive
Sample 2	10	Reactive	Reactive

REPRODUCIBILITY

Reproducibility was calculated from 10 successive determinations over 3 consecutive days, using 2 different samples, obtaining the following results:

Samples	Nº of Days	Nº of Repetitions	Expected Result	Found Result
Sample 1	03	10	Reactive	Reactive
Sample 2	03	10	Reactive	Reactive

CLINICAL SIGNIFICANCE

Leishmaniasis is an infectious and zoonotic disease caused by the protozoan *Leishmania sp.*, and transmitted by the bite of the female fo the sand fly *Lutzomyia longipalpis*, also known as the straw mosquito. The incubation period for leishmaniasis ranges from 3 to 7 months to several years, and the main clinical manifestations are generalized lymphadenopathy, exfoliative dermatitis, progressive emaciation, and keratoconjunctivitis. Some dogs may also present ulcerative or nodular mucocutaneous and mucosal lesions, epistaxis, lameness, atrophic myositis of the masticatory muscles. Due to the high frequency of asymptomatic individuals, serological diagnosis is an effective and safe tool. Due to the high sensitivity and specificity, ease of execution, the rapid test is the methodology indicated by the Ministry of Health for screening dogs for *Leishmania* infection.

BIBLIOGRAPHIC REFERENCES

- 1- QUIBASA: Data from the Research and Development Department.
- 2- KILLICK-KENDRICK R, KILLICK-KENDRICK M, PINELLI E, DEL REAL G, MOLINA R, VITUTIA MM, CANAVATE MC, NIETO J. A laboratory model of canine leishmaniasis: the inoculation of dogs with *Leishmania infantum* promastigotes from midguts of ex- perimentally infected phlebotomine sandflies. Parasite, 1994, 7, 311-318.
- 3- RIBEIRO VM. Leishmanioses. In: Associação Nacional de Clínicos Veterinários de Pequenos Animais; 3- DE NARDI AB, ROZA MR, organizadores. PROMOVET Pequenos Animais: Programa de Atualização em Medicina Veterinária: Ciclo 1. Porto Alegre: Artmed Panamericana; 2016.p.107- 50.
- 4- Ministério da Saúde. Manual de Vigilância e Controle da Leishmaniose Visceral. 1ª ed. Brasília- DF, 2006.
- 5- Ministério da Saúde, 2011. Nota Técnica Conjunta Nº1/2011. Esclarecimentos sobre a substituição do protocolo diagnóstico da leishmaniose visceral canina (LVC).

QUALITY GUARANTEE

Before being released for consumption, all **Bioclin** reagents are tested by the Quality Control Department. The quality of the reagents is assured until the validity date mentioned on the packaged, as long as they are stored and transported under the appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda
Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Tel.: +55 31 3439.5454 | E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

CONSUMER SERVICE

Customer Support Service
Tel.: 0800 0315454 | E-mail: sac@bioclin.com.br

Product licensed in the Ministry of Agriculture, Livestock and Supply since 30/12/2019 with the number 10.299/2019.

Technical manager: Dra. Camila Eckstein (CRMV/MG 20611).

Review: February/2024

REF

CATALOG NUMBER

LOT

LOT NUMBER

MANUFACTURING DATE

VALIDITY DATE
(last day of the month)

TEMPERATURE LIMIT
(store)

Σ

CONTENT IS SUFFICIENT
FOR <N> TEST

SEE INSTRUCTIONS
FOR USE

IVD

IN VITRO DIAGNOSTIC
PRODUCT

KEEP AWAY
FROM SUNLIGHT

DO NOT REUSE

CAUTION

MADE BY

CONTROL

CONTROL

CONTROL +

POSITIVE CONTROL

CONTROL -

NEGATIVE CONTROL

BIOLOGICAL RISK

FLAMMABLE

CORROSIVE

TOXIC

DO NOT USE IF
PACKAGE IS
DAMAGED

STERILE

PRODUCT
STERILIZED

DANGER