LEISHMANIOSE VH BIO

REF K218

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FINALIDADE

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- Ritmeijer K. et al. Evaluation of a new recombinant K39 rapid diagnostic test. Am. J. Trop. Med. Hyg., 74(1), pp. 76–80, 2006.
- Costa MM, Penido M, dos Santos MS, Doro D, de Freitas E, et al. Improved Canine and Human Visceral Leishmaniasis Immunodiagnosis Using Combinations of Synthetic Peptides in Enzyme-Linked Immunosorbent Assay. PLoS Negl Trop Dis 6(5): e1622, 2012.
- Om Prakash Singh and Shyam Sundar. Developments in Diagnosis of Visceral Leishmaniasis in the Elimination Era. Journal of Parasitology Research. Volume 2015, Article ID 239469, 2015.
- Srivastava et al. Diagnosis of visceral leishmaniasis. Trans R Soc Trop Med Hyg. January; 105(1): 1–6, 2011.
- Pearson, R. et al. Mem. Inst. Oswaldo Cruz, Rio de Janeiro, Vol. 84 (2): 157-166, 1989.
- QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

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PRINCIPLE OF ACTION

Methodology: Immunochromatography
LEISHMANIASIS VH BIO is an immunochromatographic test intended for the qualitative detection of Anti-*Leishmania spp* IgG Antibodies (*L. donovani*, *L. chagasi* and *L. infantum*) in serum, plasma (EDTA or Heparinized) or whole blood samples. In this kit, the recombinant antigen K39 (rK39) from *Leishmania spp* (*L. donovani*, *L. chagasi* and *L. infantum*) is immobilized in the region of the test line (T). When a sample is added, Anti-*Leishmania spp* Antibodies (*L. donovani*, *L. chagasi* and *L. infantum*) present in the sample react with protein A conjugated with colloidal gold nanoparticles. This complex migrates along the strip and interacts with rK39 in the test line region (T). If the sample contains Anti-*Leishmania spp* Antibodies (*L. donovani*, *L. chagasi* and *L. infantum*), a colored band will appear in the test region (T) indicating a positive result. If the sample does not contain Anti-*Leishmania spp* Antibodies (*L. donovani*, *L. chagasi* and *L. infantum*), this band will not appear indicating a negative result. Protein A is also previously immobilized in the control line (C). When a sample is added, the Antibodies present react with the conjugated protein A forming the Antibody-conjugated Protein A complex, triggering the reaction with the control line reagent (C). At the end of the run, the formation of the colored band in the control region (C) demonstrates the perfect functioning of the test.

REAGENTS

Number 1 (R1) Cassette - Store between 2 and 30°C. **Do not freeze.** Contains:

- A)** 01 plastic cassette with 01 test strip containing:
- Control region (C): Protein A;
 - Test region (T): Recombinant antigen (rK39) from *Leishmania spp* (*L. donovani*, *L. chagasi* and *L. infantum*);
 - Conjugate: Protein A and Colloidal Gold;

B) 01 sachet of silica;

Number 2 (R2) Diluent - Store between 2 and 30°C. **Do not freeze.** **Contains:** Buffer Solution, surfactant and preservative.

Accessories for Sample Collection:

- Disposable plastic pipette
- Alcohol Sachet
- Lancet (Sterile)

PRESENTATION

Reagent	Number of Cassettes per Packaging	Diluent
Cassette + Diluent (Pipette included)*	5	1.0 mL
	10	1.5 mL
	15	2.0 mL
	20	3.0 mL
	25	3.5 mL
	30	4.0 mL
	40	5.0 mL
Cassette + Diluent (Lancet and Pipette included)**	50	6.0 mL
	100	2 x 6.0 mL
	5	1.0 mL
	10	1.5 mL
	15	2.0 mL
	20	3.0 mL
	25	3.5 mL
Cassette + Diluent (Lancet, Pipette and Alcoholic Sachet included)***	30	4.0 mL
	40	5.0 mL
	50	6.0 mL
	100	2 x 6.0 mL
	5	5 x 0.250 mL
	10	10 x 0.250 mL
	15	15 x 0.250 mL
	20	20 x 0.250 mL
	25	25 x 0.250 mL
	30	30 x 0.250 mL
	40	40 x 0.250 mL
	50	50 x 0.250 mL
	100	100 x 0.250 mL

* These presentations are accompanied by disposable plastic pipettes, proportionally to the number of cassettes.
** These presentations are accompanied by disposable plastic lancets and pipettes, proportionally to the number of cassettes.
*** These presentations are accompanied by lancets, disposable plastic pipettes and alcoholic sachets, proportionally to the number of cassettes.

EQUIPMENTS AND OPERATING INPUTS

Pipettes, tips, watch or stopwatch, 70% alcohol (w/w), cotton and garbage for disposal. These items are found in the specialized market for articles for Clinical Analysis Laboratories.

STORAGE AND TRANSPORT CONDITIONS

Transport can be carried out at temperatures up to 45°C for a maximum of 3 weeks. The storage temperature should be 2 to 30°C. Avoid moisture. **Do not freeze.**

SPECIAL CARES

- 1- For *in vitro* diagnostic use only.
- 2- Strictly follow the proposed methodology to obtain exact results.
- 3- The cassettes (test devices) do not contain contaminating substances, but as the samples are biological, the safety routines must be strictly followed and the use of disposable gloves is essential.
- 4- Open the aluminum envelope to remove the cassette, only when using it.
- 5- The kit must not be frozen or exposed to high temperatures (above 45°C), as this will cause deterioration.

- 6- The reagents, as well as the samples, must reach room temperature before starting the test.
- 7- Do not mix reagents from kits with the same batch or with different batch numbers.
- 8- Do not use the kit beyond its expiration date. The expiry date is indicated on the product label.
- 9- The cassette (test device), the plastic pipette and the lancet are disposable and cannot be reused.
- 10- The kit materials that came into contact with the patient's sample must have a special destination. Disposal of lancets must be done in sharps waste. Pipettes, cassettes (test devices) and diluent vials must be disposed of in biowaste.
- 11- Secondary packaging (box), Instructions for Use, aluminum envelope and silica sachet can be disposed of in regular trash.
- 12- We recommend applying the local, state and federal norms for environmental protection so that the disposal of reagents and biological material is carried out in accordance with current legislation.
- 13- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Service Customer Advice) by Quibasa.
- 14- Do not use the product if the packaging is damaged.
- 15- It is imperative that the instruments and equipment used are properly calibrated and submitted to periodic maintenance.

SAMPLES

Serum or Plasma (EDTA or Heparin): Samples can be stored under refrigeration, between 2 and 8°C, for a maximum period of 5 days. If samples cannot be analyzed within 5 days, they can be stored for up to 30 days at -20°C (freezer).

Whole Blood: It can be collected by venipuncture using the appropriate anticoagulant (EDTA or Heparin) or by digital puncture with the aid of a lancet. Whole blood collected by finger prick should be used for immediate testing. Whole blood collected by venipuncture with anticoagulant can be stored at 2 to 8°C for up to 1 day. **Do not freeze.**


PROCESS DESCRIPTION

STABILITY AFTER OPEN


The stability of the LEISHMANIASIS VH BIO cassette (after opening the aluminized sachet) is 2 hours. This stability may vary according to the environmental conditions after opening the aluminized sachet.

TECHNIQUE - VENOPUNCTION

- 1- The sample must be at room temperature before starting the test.
- 2- Remove the cassette from its protective packaging, place it on a clean and level surface and identify it properly.
- 3- For Serum or Plasma: Transfer 10 µL of serum or plasma into the sample well. Hold the Diluent bottle (Reagent 2) vertically and apply 3 drops (80 µL) into the sample well.


- 4- For Whole Blood: Transfer 10 µL of whole blood into the sample well. Hold the Diluent bottle (Reagent 2) vertically and apply 3 drops (80 µL) into the sample well.

- 5- Wait for the lines to form. Interpret results between 10 and 20 minutes. **Do not interpret after 20 minutes.**


TECHNIQUE - DIGITAL PUNCTURE / LANCET

- 1- Remove the cassette from its protective packaging, place it on a clean and level surface and identify it properly.
- 2- Wash and dry the patient's hands. Apply 70% alcohol (w/w) on the tip of the finger to antiseptics the area used.
- 3- Press the finger so that there is an accumulation of blood on the tip of the finger that will be perforated by the lancet.
- 4- Remove the protective cap from the lancet.
- 5- Position and press the lancet firmly over the area to be punctured. Then the blood will come out through the pierced area.
- 6- Collect the blood using the disposable plastic pipette that comes with the kit. Squeeze the plastic pipette above the marked line and place the open well of the plastic pipette against the blood drop. Release the pressure on the plastic pipette so that the blood is aspirated. Aspirate blood to the marked line on the disposable plastic pipette (equivalent to 10 µL of blood).
- 7- In the sample well, dispense the blood by pressing the pipette again.
- 8- Hold the Diluent bottle (Reagent N° 2) vertically and apply 3 drops (80 µL) in the sample well.
- 9- Wait for the lines to form. Interpret results between 10 and 20 minutes. **Do not interpret after 20 minutes.**

INTERPRETATION OF THE RESULT


Reagent Test: Formation of two red lines, one in the control region (C) and one in the test region (T), after 10 minutes. **Do not interpret after 20 minutes.**

Test Non-Reactive: Formation of a red line in the control region (C) and complete absence of line in the test region (T) after 10 minutes. **Do not interpret after 20 minutes.**

Invalid Result: The absence of line formation in the control region (C), indicates an error in the procedure or deterioration of the cassette. In this case, repeat the test using a new cassette.

Notes:

- 1- The line in the control region will appear before the 10 minutes of incubation of the test, and this does not mean that the results can be interpreted before the expected time.
- 2- It is important to use the correct sample volume, as lower or higher volumes may determine erroneous results.
- 3- The results provided by this kit must be interpreted by the responsible professional, not being the only criterion for determining the diagnosis and/or treatment of the patient.
- 4- Isolated test results cannot be used for definitive diagnosis.
- 5- A non-reactive result in the Leishmaniasis VH Bio kit does not exclude infection by *Leishmania spp* (*L. donovani*, *L. chagasi* and *L. infantum*). If the suspicion of infection persists, a new sample must be collected after 30 days for a repeat test.

INTERFERING

No significant interferences were found up to 10 g/L of Hemoglobin, up to 30 g/L of lipids, 0.2 g/L of Bilirubin, 60 g/L of protein and up to 280 IU/mL of Rheumatoid Factor. No cross-reaction with Chagas disease, malaria and toxoplasmosis was verified.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control program, where procedures, standards, limits and tolerance for variations are clearly established. It is important to point out that all measurement systems have a characteristic analytical variability, which must be monitored by the laboratories themselves. For that, it is recommended the use of controls, which allow evaluating the precision and accuracy of the dosages.

PRODUCT PERFORMANCE

CLINICAL SENSITIVITY AND SPECIFICITY

The LEISHMANIASIS VH BIO kit has been tested using clinically known samples. For this test, 210 samples were analyzed, being 110 reagents for Leishmaniasis and 100 non-reagents for Leishmaniasis. The following results were found:

	Expected Result	LEISHMANIASIS VH BIO
Reagent Sample	110	109
Non-Reactive Sample	100	100
Total Samples	210	

Clinical sensitivity: 99.09% (109/110) (IC 95%: 97.33 - 100%)
Clinical specificity: >99.99% (100/100) (IC 95%: 91 - 100%)

PRECISION

Repeatability

The repeatability was calculated from 10 successive determinations, using 2 different samples, obtaining the following results:

Sample	Nº of Repetitions	Expected Result	Found Result
01	10	Reagent	Reagent
02	10	Reagent	Reagent
03	10	Reagent	Reagent

Reproducibility

The reproducibility was calculated from 10 successive determinations during 3 consecutive days, using 2 different samples, obtaining the following results:

Sample	Nº of Days	Nº of Repetitions per Day	Expected Result	Found Result
01	03	10	Reagent	Reagent
02	03	10	Reagent	Reagent
03	03	10	Reagent	Reagent

DIAGNOSTIC SIGNIFICANCE

Visceral Leishmaniasis, also known as kala-azar, is a non-contagious systemic parasitic infectious disease caused by the protozoan *Leishmania spp* (*L. donovani*, *L. chagasi* and *L. infantum*). It is transmitted through the bite of sandflies, with *Lutzomyia longipalpis* being the main vector species in Brazil. The disease is characterized by a progressive evolution that can lead to death if not treated. Its main symptoms are: weight loss, irregular fever, anemia, weakness (reduced muscle strength), abdominal swelling, enlarged liver and spleen. A patient is considered clinically suspicious when he has a prolonged history of fever (more than two weeks) with splenomegaly and weight loss.

According to the World Health Organization, only 50 to 60% of patients with the aforementioned clinical findings have Visceral Leishmaniasis. Therefore, a laboratory differential diagnosis is essential. The main diseases with symptoms similar to kala-azar are: malaria, schistosomiasis, brucellosis, typhoid fever, splenic abscess, myeloproliferative diseases, leukemias, lymphomas and chronic hemolytic anemia. The main laboratory diagnostic methods for Visceral Leishmaniasis are serological, including direct agglutination test, ELISA and rapid test (immunochromatography), and parasitological tests such as microscopic analysis of lymph node aspirate, splenic aspirate and bone marrow aspirate.

BIBLIOGRAPHIC REFERENCES

1. Ritmeijer K. et al. Evaluation of a new recombinant K39 rapid diagnostic test. Am. J. Trop. Med. Hyg., 74(1), pp. 76–80, 2006.
2. Costa MM, Penido M, dos Santos MS, Doro D, de Freitas E, et al. Improved Canine and Human Visceral Leishmaniasis Immunodiagnosis Using Combinations of Synthetic Peptides in Enzyme-Linked Immunosorbent Assay. PLoS Negl Trop Dis 6(5): e1622, 2012.
3. Om Prakash Singh and Shyam Sundar. Developments in Diagnosis of Visceral Leishmaniasis in the Elimination Era. Journal of Parasitology Research. Volume 2015, Article ID 239469, 2015.
4. Srivastava et al. Diagnosis of visceral leishmaniasis. Trans R Soc Trop Med Hyg. January; 105(1): 1–6, 2011.
5. Pearson, R. et al. Mem. Inst. Oswaldo Cruz, Rio de Janeiro, Vol. 84 (2): 157-166, 1989.
6. QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

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