MAGNÉSIO MONORREAGENTE

REF K314

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FINALIDADE

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REF K314

INSTRUCTIONS FOR USE

TO OBTAIN THE INSTRUCTIONS FOR USE IN PRINTED FORMAT, AT NO ADDITIONAL COST,
CONTACT CUSTOMER ADVISORY SERVICE:

SAC: +55 (31) 3439 5454 / 0800 031 5454 / sac@bioclin.com.br

FUNCTION

Method for determination of Magnesium in biological samples of serum, plasma (Heparin), liquor and urine. Colorimetric test, for *in vitro* diagnostics use only.

PRINCIPLE OF ACTION

Methodology: Mann and Yoe/ Magon Sulfonate/ Xilidil Blue Mann and Yoe dye in alkaline pH and in the presence of Magnesium develop a red color. The intensity of the red color of the complex is proportional to the concentration of Magnesium. The method does not require deproteinization, being more sensitive than the Titan Yellow method, allowing therefore the use of only 10 µL of the sample to be analyzed, what turns the technique excellent for pediatric use. Besides its simplicity, the method does not suffer interference from Calcium Gluconate. The presence of surfactants rules out the interference of lipemic serum.

Reagents

Number 1 - Color Reagent - Store between 2 and 8°C. Contains: Sulphonated Magon < 1 mmol/L, Anhydrous Potassium Carbonate < 300 mmol/L, stabilizers, surfactant and preservative.

Number 2 - Standard* - Store between 2 and 8°C. Contains: Magnesium 2 mg/dL and preservative.

PRESENTATION

Reagent N° 1	Reagent N° 2*
1 x 100 mL	1 x 3 mL
2 x 100 mL	1 x 3 mL
4 x 100 mL	1 x 3 mL
5 x 20 mL	1 x 3 mL
5 x 40 mL	1 x 3 mL
10 x 40 mL	1 x 3 mL

* The R2 Standard is not present in all presentations.

EQUIPMENTS AND OPERATIONAL INPUTS

Spectrophotometer or colorimeter, watch or stopwatch, pipettes, test tubes, Biocal Kit, Biocontrol N and Biocontrol P Bioclin. Can be found in markets specialized in Laboratories of Clinical Analysis.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be between 2 to 8°C. The transport at temperatures up to 30°C should not exceed 5 days. Protect from light and avoid moisture. **Do not freeze.**

SPECIAL CARE

- 1- For professional *in vitro* diagnostic use only.
- 2- Strictly follow the methodology proposed to obtain exact results.
- 3- Water used in the cleansing of the material should be recent and free of contaminants.
- 4- Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers that may alter the results significantly.

5- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

6- To obtain information related to biosafety or in case of accidents with the product, consult the SDS (Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.

7- Do not use the product in case of damaged packaging.

8- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Serum obtained free of hemolysis, plasma collected with Heparin, liquor and 24h urine.

The serum and plasma sample is stable for 7 days at 2 to 8°C, and 1 year at -20°C.⁶

The urine sample is stable for 3 days at 2 to 8°C and 1 year at -20°C.⁶

PROCESS DESCRIPTION

The calibration stability of the Magnesium Single Reagent kit installed on refrigerated equipment is at least 03 days. This stability may vary depending on the conditions of the test, equipment and environment. Therefore, it is suggested to follow the product performance using control serum.

TECHNIQUE

Bioclin recommends, use as calibrator the Biocal Bioclin Kit and control serum, Biocontrol N and P Bioclin Kits.

Mark 3 test tubes: B (Blank), A (Sample), S (Standard) and proceed as follows:

	Blank	Standard	Sample
Reagent N° 1	1.0 mL	1.0 mL	1.0 mL
Sample	----	----	10 µL
Reagent N° 2	----	10 µL	----

Homogenize well and read the absorbance from Standard and from the Sample in 500 nm, correcting the zero with the Blank. The immediate color is stable for 30 minutes.

24 hour urine – Sample preparation:

Homogenize urine and measure its volume. Add 1 drop of concentrated HCl in 20 mL of this urine. Dilute 1.0 mL of acidified urine with 4.0 mL of distilled water (1:5 dilution). Homogenize it well by agitating and proceed the determination of Magnesium the same way that was proposed to the serum. Multiply the results obtained by 5 (dilution factor).

CALCULATIONS

Serum

$$\text{Magnesium (mg/dL)} = \frac{\text{Sample Absorbance} \times 2}{\text{Standard Absorbance}}$$

Since the reaction follows the Lambert-Beer law, the Calibration Factor may be used.

$$\text{Calibration Factor} = \frac{\text{Standard Concentration (2 mg/dL)}}{\text{Standard Absorbance}}$$

$$\text{Magnesium (mg/dL)} = \text{Sample absorbance} \times \text{Calibration Factor}$$

Urine

$$\text{Magnesium (mg/24h)} = \frac{\text{mg/dL} \times \text{Urinary Volume (mL)}}{100}$$

Results are expressed in mg/dL for serum and mg/24h for urine.

INTERFERENCES

No interference was observed by Bilirubin until 8 mg/dL and Hemoglobin until 250 mg/dL.

Usage of hemolyzed samples may result in falsely elevated results, since the erythrocytes contain 3 times more Magnesium than serum or plasma.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

TRACEABILITY

The kit's standard is traceable to the reference material NIST (National Institute of Standards and Technology) SRM 929A.

REFERENCE VALUES

Reference values expressed in mg/dL, for the present method, were obtained through the determination of Magnesium in healthy populations of male and female.

Serum or plasma..... 1.6 to 2.4 mg/dL (All ages)
Liquor..... 2.4 to 3.4 mg/dL
Urine..... 32 to 150 mg/24h
(Varies with according to eating habits)

In order to convert values of mg/dL into mmol/L (SI) multiply it by 0.41.

These values should be used as an orientation, being that each laboratory should create its own range of reference values, according to the population served.

The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

PRODUCT PERFORMANCE

ACCURACY

Comparison of Methods

The Magnesium Single Reagent was compared with other commercially available methods for measurement of Magnesium. 42 analyzes were performed and the results were evaluated. The linear equation obtained was $Y = 1.0091X - 0.0683$ and the correlation coefficient 0.9923. With these results, we can conclude that the kit shows good methodological specificity.

PRECISION

Repeatability

The repeatability was calculated from 40 successive determinations, using 3 samples with different concentrations, obtaining the following results:

Repeatability	Sample 1	Sample 2	Sample 3
Average Concentration (mg/dL)	0.96	2.01	5.13
Standard Deviation (mg/dL)	0.05	0.03	0.10
Coefficient of Variation (%)	5.17	1.66	1.97

Reproducibility

The reproducibility was calculated from 40 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

Reproducibility	Sample 1	Sample 2	Sample 3
Average Concentration (mg/dL)	0.96	2.00	5.04
Standard Deviation (mg/dL)	0.03	0.01	0.09
Coefficient of Variation (%)	3.41	0.66	1.79

SENSITIVITY

The sensitivity was calculated from 40 determinations of a sample free the Magnesium. The average found, that indicates the method detection limit, was 0.01 mg/dL and the standard deviation of 0.01 mg/dL. The sensitivity, which indicates the method quantification limit, corresponds the average plus 3 times the standard deviation and is equal to 0.04 mg/dL.

LINEARITY

Reaction is linear up to the concentration of 5.0 mg/dL. For higher values, dilute the sample with NaCl 0.85%, repeat the determination and multiply the results by the dilution factor.

DIAGNOSTIC SIGNIFICANCE

Magnesium is an electrolyte found primarily in intracellular liquid and bone. Participates as a cofactor in several enzymatic systems, in the metabolism of carbohydrates, muscle contraction, blood clotting and is essential in preserving the molecular structure of DNA, RNA and ribosomes. Magnesium levels decreased in plasma are associated with tetany, weakness, disorientation and drowsiness, which effect the ionized magnesium deficiency. Clinical manifestations seizures associated with hypocalcemia and hypomagnesemia, attributed to a selective defect in intestinal absorption of Magnesium, stressing the importance of dosage this ion in pediatrics. The most frequent causes of low concentrations of Magnesium are: chronic diarrhea, acute pancreatitis, alcoholism, chronic hepatitis, diabetes mellitus, hypoparathyroidism, hyperthyroidism, hyperaldosteronism. Hypermagnesemia, although rare, can be registered in cases of renal failure, severe dehydration, intensive treatment with Magnesium salts.

BIBLIOGRAPHIC REFERENCES

- 1 - MANN, C. K. & Yoe, J. H.: Ann. Chem. 28:202, 1956.
- 2 - MANN, C. K. & Yoe, J. H.: Anal. Chim. Acta. 16:155, 1957.
- 3 - RICE, E. W. & LAPARA, C. Z.: Clin. Chim. Acta. 10:260, 1964.
- 4 - WEISSMANN, N. & PILEGGI, V. J. (1974) in Clinical Chemistry Principles and Technics 2nd. Ed. Henry R., Cannon, D. C. e Winkelman, J. W. p.678. Haper and Row Publishers.
- 5 - TONKS, D. B., Clin. Chem. 9:217, 1963.
- 6 - WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.
- 7 - QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

 **QUIBASA QUÍMICA BÁSICA Ltda**
Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone: + 55 31 3439.5454 | E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service
Phone.: 0800 0315454 | E-mail: sac@bioclin.com.br

ANVISA resgistration for Magnesium Single Reagent Family kit: 10269360373

Review: May/2025

UNIVERSAL SYMBOLOGY

 CATALOG NUMBER

 LOT NUMBER

 MANUFACTURING DATE

 VALIDITY DATE
(last day of the month)

 TEMPERATURE LIMIT
(store)

 CONTENT IS SUFFICIENT
FOR <N> TEST

 SEE INSTRUCTIONS
FOR USE

 IN VITRO DIAGNOSTIC
PRODUCT

 KEEP AWAY
FROM SUNLIGHT

 DO NOT REUSE

 CAUTION

 MADE BY

 CONTROL

 POSITIVE CONTROL

 NEGATIVE CONTROL

 BIOLOGICAL RISK

 FLAMMABLE

 CORROSIVE

 TOXIC

 DO NOT USE IF
PACKAGE IS
DAMAGED

 PRODUCT
STERILIZED

 DANGER