FINALIDADE

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1 - KILLINGSWORTH, L. M; SAVORY, J.: Automated Immunochemical Procedures for Measurement of Immunoglobulins IgG, IgA and IgM in human serum, Clin. Chem., 1971; 17: 936.

2 - PESCE, A. J.; KAPLAN, L. A. : Methods in Clinical Chemistry C. V. Mosby Company, 1987.

3 - HELLSING, K., Profiles in the Biological Fluids, 1973; 23: 579.

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USAGE INSTRUCTIONS

FUNCTION

Reagent used for calibration of turbidimetrics tests, for *in vitro* diagnostic use only.

REAGENTS

Number 1 - Multicalibrator - Store between 2 and 8°C. Contains: Serum Proteins and preservative. **Potentially infectious.**

Warning: The concentration of the analytes varies according to the batch.

See table that accompanies the product.

PRESENTATION

Reagent N° 1	1 mL
--------------	------

EQUIPMENTS AND OPERATIONAL INPUTS

Automated or semi-automated biochemical devices, watch or stopwatch, pipettes, test tubes, kits and water bath at 37°C. They can be found at markets specialized on Laboratories of Clinical Analysis.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be between 2 to 8°C. The transport, at temperatures up to 30°C, should not exceed 5 days. Protect from light and avoid moisture.

Do not freeze.

SPECIAL CARE

1- For *in vitro* diagnostic use only.

2- Strictly follow the methodology proposed to obtain exact results.

3- Water used in material cleaning must be recent and free of contaminants.

4- Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers that may alter the results significantly.

5- Reagent 1 and 2 should be handled cautiously since they are subject to biological contamination.

6- Handle with care, as the reagent contains Sodium Azide, which is irritant to skin and mucous.

7- The reagent has been tested for antibodies anti-HIV and HBs Antigen using state of the art methods and found negative results. The risk of infection cannot be excluded and the reagent should be handled as cautiously as observed for the patient's serum. **Potentially infectious.**

8- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

9- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.

10- Do not use the product in case of damaged packaging.

11- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

PROCESS DESCRIPTION

The utilization of the calibrator must comply with specific technique for the equipment to be used in calibrations.

The concentration of analytes varies each lot. See table that accompanies the product.

BIBLIOGRAPHIC REFERENCES

1 - KILLINGSWORTH, L. M; SAVORY, J. : Automated Immunochemical Procedures for Measurement of Immunoglobulins IgG, IgA and IgM in human serum, Clin. Chem., 1971; 17: 936.

2 - PESCE, A. J.; KAPLAN, L. A. : Methods in Clinical Chemistry C. V. Mosby Company, 1987.

3 - HELLSING, K., Profiles in the Biological Fluids, 1973; 23: 579.

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QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca

CEP 31565-130 - Belo Horizonte - MG - Brasil

Phone.: +55 31 3439.5454| E-mail: bioclin@bioclin.com.br

CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service

Phone.: 0800 0315454 | E-mail: sac@bioclin.com.br

ANVISA registration for the Multical kit: 10269360077

Review: August/2022

UNIVERSAL SYMBOLOGY

CATALOG NUMBER

LOT NUMBER

MANUFACTURING DATE

VALIDITY DATE
(last day of the month)

TEMPERATURE LIMIT
(store)

CONTENT IS SUFFICIENT
FOR <N> TEST

SEE INSTRUCTIONS
FOR USE

IN VITRO DIAGNOSTIC
PRODUCT

KEEP AWAY
FROM SUNLIGHT

DO NOT REUSE

CAUTION

MADE BY

CONTROL

POSITIVE CONTROL

NEGATIVE CONTROL

BIOLOGICAL RISK

FLAMMABLE

CORROSIVE

TOXIC

DO NOT USE IF
PACKAGE IS
DAMAGED

PRODUCT
STERILIZED

DA