MULTICONTROL 1

REF **K246-1**

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INSTRUCTIONS FOR USE

FUNCTION

Reagent used as a control for turbidimetric tests, for *in vitro* diagnostic use only.

REAGENTS

Number 1 (R1) - Multicontrol 1 - Store between 2 and 8 C. Contains: Serum Proteins and preservative. **Potentially infectious.**

Warning: The concentration of analytes varies each lot. See table that accompanies the product.

PRESENTATION

Reagent Nº 1	1 mL
Reagent Nº 1	2 mL
Reagent Nº 1	3 mL
Reagent Nº 1	4 mL
Reagent Nº 1	5 mL
Reagent Nº 1	10 mL
Reagent Nº 1	15 mL
Reagent Nº 1	20 mL

EQUIPMENTS AND OPERATIONAL INPUTS

Automated biochemical devices and kits. They can be found at markets specialized in laboratories of clinical analysis.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be between 2 to 8°C.
The transport at temperatures up to 30°C should not exceed 5 days.
Protect from light and avoid moisture. **Do not freeze.**

SPECIAL CARE

- 1- For *in vitro* diagnostic use only.
- 2- Strictly follow the methodology proposed to obtain exact results.
- 3- Water used in material cleaning must be recent and free of contaminants.
- 4- Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers that may alter the results significantly.
- 5- Handle with care, as the reagent contains Sodium Azide, which is irritant to skin and mucous.
- 6- The reagent has been tested for anti-HIV antibodies, and HBs antigen using state of the art methods and found negative results. The risk of infection cannot be excluded and the reagent should be handled as cautiously as observed for the patient's serum. **Potentially Infectious.**
- 7- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.
- 8- To obtain information related to biosafety or in case of accidents with the product, consult the SDS (Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.
- 9- Do not use the product in case of damaged packaging.
- 10- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

PROCESS DESCRIPTION

The use of the control must obey the specific technique for the equipment to be used.

The concentration of the analytes varies with each batch. See table accompanying the product.

TO OBTAIN THE INSTRUCTIONS FOR USE IN PRINTED FORMAT, AT NO ADDITIONAL COST, CONTACT CUSTOMER ADVISORY SERVICE:

SAC: +55 (31) 3439 5454 / 0800 031 5454 / sac@bioclin.com.br

BIBLIOGRAPHIC REFERENCES

- 1 - KILLINGSWORTH, L. M.; SAVORY, J.: Automated Immunochemical Procedures for Measurement of Immunoglobulins IgG, IgA and IgM in human serum, Clin. Chem., 1971; 17, 936.
- 2 - PESCE, A. J.; KAPLAN, L. A.: Methods in Clinical Chemistry C. V. Mosby Company, 1987.
- 3 - HELLSING, K., Profiles in the Biological Fluids, 1973 23: 579.
- 4 - BLOM, M. and HJOME, H., Clinical Chemistry, 1976; 22: 657.
- 5 - HILLS, L. P. and TIFFANY, T. I., Clinical Chemistry, 1980; 26: 1459.
- 6 - QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

QUALITY ASSURANCE

Before being released for consumption, all Bioclin reagents are tested by the Quality Control Department. The quality of the reagents is ensured until the expiration date mentioned on the presentation packaging, provided they are stored and transported under the appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte/MG - Brasil
Tel.: +55 31 3439 5454 | E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Indústria Brasileira

CUSTOMER SERVICE

Customer Advisory Service
Phone: +55 31 3439 5454 | E-mail: customerservice@bioclin.com.br

ANVISA registration for the Multicontrol 1 kit: 10269360340

Review: October/2025

UNIVERSAL SYMBOLOGY

	CATALOG NUMBER		MADE BY
	LOT NUMBER		CONTROL
	MANUFACTURING DATE		POSITIVE CONTROL
	VALIDITY DATE (last day of the month)		NEGATIVE CONTROL
	TEMPERATURE LIMIT (store)		BIOLOGICAL RISK
	CONTENT IS SUFFICIENT FOR <N> TEST		FLAMMABLE
	SEE INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC PRODUCT		TOXIC
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED
	DO NOT REUSE		PRODUCT STERILIZED
	CAUTION		DANGER