PADRÃO HEMOGLOBINA

REF K029

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FINALIDADE

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USAGE INSTRUCTIONS

FUNCTION

Reagent for hemoglobin dosage standardization, for *in vitro* diagnostic use only.

REAGENTS

Number 1 - Standard - Store between 2 and 8°C. Contains: Hemoglobin, surfactant and stabilizer. **Potentially infectious.**

Warning: Hemoglobin concentration may vary according to the lot - See flask label.

PRESENTATION

Reagent	Volume
Number 1	1.0 mL

EQUIPMENTS AND OPERATIONAL INPUTS

Spectrophotometer or colorimeter, stopwatch, pipettes, test tubes and Hemoglobin Bioclin kit. These items are found at markets specialized in Laboratories of Clinical Analysis.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be between 2 to 8°C. The transport at temperatures up to 30°C should not exceed 5 days. Protect from light and avoid moisture.

SPECIAL CARE

- For *in vitro* diagnostic use only.**
- Strictly follow the methodology proposed to obtain exact results.
- Water used in material cleaning must be recent and free of contaminants.
- Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers that may alter the results significantly.
- Standard was analyzed for the detection of antibodies anti HIV, anti HCV and HBs antigens, with negative results. However for a greater security consider handling it as **potentially infectious**.
- The tips from automatic pipettes must be clean, dry and in good conditions of usage.
- Working reagents and the reaction residues must never be deployed in the presence of acid reagents, for that takes to the formation of toxic gas.
- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.
- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.
- Do not use the product in case of damaged packaging.
- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

PROCESS DESCRIPTION

TECHNIQUE

Homogenize before each use.

Use Working Reagent* from **Bioclin Hemoglobin** kit. To obtain Calibration Factor, dose the Standard in triplicate.

To do this, mix the reagents according to the information in the table below.

Working Reagent*	2.5 mL
Standard	10 µ

Homogenize and wait for 5 minutes in temperatures between 15 and 30°C, reading it at 540 nm (520 - 550 nm) settling the zero with the Working Reagent*. Colored reaction is stable for 60 minutes.

CALCULATIONS

Calibration Factor = $\frac{\text{Standard Concentration}}{\text{Standard Absorbance}}$

Hemoglobin (g/dL) = $\frac{\text{Sample Absorbance}}{\text{Absorbance}} \times \text{Calibration Factor}$

Example:

Sample Absorbance = 0.430

Standard Average Absorbance = 0.330

Standard Value = 11.6 g/dL

Calibration Factor = $\frac{11.6}{0.330} = 35.15$

Hemoglobin (g/dL) = $0.430 \times 35.15 = 15.11$

PROCEDURE LIMITATIONS

The absorbances of the Standard should not vary more than 5% so that the average of these can be used to calculate the calibration factor.

TRACEABILITY

The Hemoglobin Standard is traceable to the reference material NIST (National Institute of Standards and Technology) SRM 931.

PRODUCT PERFORMANCE

ACCURACY

Recovery

A recovery analysis was performed with 05 determinations of standard. Accuracies were calculated and found in good agreement with the reference values, obtaining a recovery between 98 and 102%.

PRECISION

Repeatability

The repeatability was calculated from 20 successive determinations, using 01 sample, obtaining the following results:

Repeatability	Standard
Average Concentration (g/dL)	11.48
Standard Deviation (g/dL)	0.25
Coefficient of variation (%)	2.22

Reproducibility

The reproducibility was calculated from 20 successive determinations for 3 consecutive days, using 01 samples, obtaining the following results:

Reproducibility	Standard
Average Concentration (g/dL)	11.50
Standard Deviation (g/dL)	0.03
Coefficient of variation (%)	0.30

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of control serum, allowing the precision and accuracy of the dosages.

BIBLIOGRAPHIC REFERENCES

- 1 - HENRY, R. J., Clinical Chemistry – Principles and Technics – Hoeber Division, Haper and Row, New York.
- 2 - TONKS, D. B., Clin. Chem., 9:217, 1983.
- 3 - DRABKIN, D. L., Austin, J. H., J. Biol. Chem. , 112:51, 1935.
- 4 - QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone: +55 31 3439.5454 | E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service
Phone: 0800 0315454 | E-mail: sac@bioclin.com.br

ANVISA registration for Hemoglobin Standard kit: 10269360087

Review: October/2022

UNIVERSAL SYMBOLOGY

	CATALOG NUMBER		MADE BY
	LOT NUMBER		CONTROL
	MANUFACTURING DATE		POSITIVE CONTROL
	VALIDITY DATE (last day of the month)		NEGATIVE CONTROL
	TEMPERATURE LIMIT (store)		BIOLOGICAL RISK
	CONTENT IS SUFFICIENT FOR <N> TEST		FLAMMABLE
	SEE INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC PRODUCT		TOXIC
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED
	DO NOT REUSE		PRODUCT STERILIZED
	CAUTION		DANGER