Bioclin

PROTEÍNAS TOTAIS MONOREAGENTE

REF

K031

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USAGE INSTRUCTIONS

FUNCTION

Method for determination of Total in biological serum, plasma (EDTA or Heparin), pleural, synovial and ascitic fluids samples. Colorimetric test. For *in vitro* diagnostic only.

PRINCIPLE OF ACTION

Methodology: Biuret

The peptidic bindings of the protein (-CONH-) reacts with the cupric ions, in an alkaline environment, forming a violet colored complex proportional to the protein content in the environment. The Sodium Tartrate and Potasyum stablishes the reagent and the appropriate concentration of Potasyum lodate prevines its self-reduction.

REAGENTS

Number 1 - Biuret - Store between 2 and 8°C. Contains: Sodium Hydroxide < 500 mmol/L, Sodium and Potasyum Tartrate < 100 mmol/L, Potasyum lodate < 20 mmol/L and Cooper Sulphate < 200 mmol/L.

Number 2 - Standard - Store between 2 and 8°C. Shake before use. Contains: Albumin 4g/dL, stabilizer and preservative.

PRESENTATION

Presentation	Reagent Nº 1	Reagent Nº 2
1	250 mL	2 mL
2	500 mL	2 mL
3	10 x 20 mL	2 mL
4	10 x 40 mL	2 mL
5	5 x 40 mL	2 mL

EQUIPMENTS AND OPERATIONAL INPUTS

Spectrophotometers or colorimeter, watch and stopwatch, pipettes and test tubes, Biocontrol N and Biocontrol P Bioclin. They can be found at markets specialized in clinical analysis laboratories.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature shall be between 2 and 8°C. The transport at temperatures up to 30°C should not exceed 5 days. Protect from light and avoid moisture. **Do not Freeze.**

SPECIAL CARE

- 1- For *in vitro* diagnostic use only.**
- 2-** Strictly follow the methodology proposed to obtain exact results.
- 3-** Water used in material cleaning must to be recent and free of contaminants.
- 4-** Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers that may alter the results significantly.
- 5-** Reagent Nº 2 is caustic, handle carefully. Reagent Nº 1 contains sodium azide, irritating to skin and mucous membrane.
- 6-** Do not use plasma.
- 7-** Determine the factor periodically and in each product batch.
- 8-** We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.
- 9-** To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Advisory Service Customer) of Quibasa.

- 10-** Do not use the product in case of damaged packaging.
- 11-** It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Serum and plasma (EDTA or Heparin) obtained free of hemolysis and fluids (pleural, synovial and ascitic). The analyte is stable for 7 days between 2 and 8°C and 60 days at negative 10°Cº, 9.

PROCESS DESCRIPTION

The calibration stability of the Total Protein Monoreagent kit installed on refrigerated equipment is up to 15 days. This stability may vary depending on the conditions of the test, equipment and environment. Therefore, it is suggested to follow the product performance using control serum.

REAGENT PREPARATION

The reagents are ready for use.

TECHNIQUE

Bioclin recommends, as control serum, Biocontrol N and P Bioclin Kits.

Mark 3 test tubes: B (Blank), P (standard), A (Sample) and proceed as follows:

	Blank	Standard	Sample
Sample	---	---	50 µL
Standard	---	50 µL	---
Biuret	2.5 mL	2.5 mL	2.5 mL

Homogenize well and let it rest for 10 minutes. Read the sample and standard absorbances at 545 nm (510 - 550 nm), matching the zero with the blank. The color is stable for 30 minutes.

CALCULATIONS

Total Protein = $\frac{\text{Sample absorbance}}{\text{Standard Absorbance}} \times 4$ (g/dL)

As the reaction follows the Lambert-Beer law, the Calibration Factor can be used.

Calibration factor = $\frac{\text{Standard Concentration}}{\text{Standard Absorbance}}$ (4 g/dL)

Total Protein = Sample Absorbance x Calibration Factor (g/dL)

The results are express in g/dL

PROCEDURE LIMITATIONS

The described method is not applied to the Protein dosage in urine and liquor. For this samples use the Bioclin Bioprot U/LCR kit.

INTERFERENCES

The hemolized sample, jaundice, lipemic or medicine such as: digitalytics, corticosteroids, anabolic steroids interfere in the dosages. The estrogen, laxatives and oral contraceptives interfere, producing reduced values of the proteins. For samples with bilirubin above 40 mg/dL and Hemoglobin above 180 mg/dL make the Blank of the Sample as follows: mixture 2.5 mL of NaCl 0.85% with 0.05 mL of sample. Determine the absorbance at 545 nm matching the zero with distilled or deionized water. Reduce the absorbance obtained from the absorbance test and calculate the concentration.

No interference was observed by Ascorbic Acid until 60 mg/dL and Rheumatoid Factor until 300 UI/mL.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

TRACEABILITY

The kit's standard is traceable to the reference material NIST (National Institute of Standards and Technology) SRM 927D.

REFERENCE VALUES

The reference values in g/dL, for this method were obtained through the determination of the Total Protein in healthy populations of male and female.

Children and Teenagers	
Premature	3.6 - 6.0 g/dL
Newborn	4.6 - 7.0 g/dL
7 days - 1 year	4.4 - 7.6 g/dL
1 - 2 years	5.6 - 7.5 g/dL
Over 3 years	6.0 - 8.0 g/dL
Adults	
6.0 - 8.0 g/dL	

These values should be used as guidelines, each laboratory should establish its range of reference values, according to the population served.

The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

PRODUCT PERFORMANCE

ACCURACY

Comparison of Methods

The Total Protein Monoreagent kit was compared wit other method for total protein dosage, commercially available. 42 analyses were performed and the results were evaluated. The linear equation obtained was Y = 0.984X + 0.0315 and the correlation coefficient 0.9978. With these results we can conclude that the kit shows good methodological specificity.

PRECISION

Repeatability

The repeatability was calculated from 40 successive determinations, using 3 samples with different concentrations, obtaining the following results:

Repeatability	Sample 1	Sample 2	Sample 3
Average Concentration (g/dL)	4.99	8.87	6.94
Standard Deviation (g/dL)	0.04	0.03	0.02
Coefficient of Variation (%)	0.75	0.39	0.24

Reproducibility

The reproducibility was calculated from 40 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

Reproducibility	Sample 1	Sample 2	Sample 3
Average Concentration (g/dL)	4.95	8.84	6.74
Standard Deviation (g/dL)	0.04	0.06	0.35
Coefficient of Variation (%)	0.77	0.73	5.22

SENSITIVITY

The sensitivity was calculated from 40 determinations of a sample free of Total Protein. The average found was 0.036 g/dL with a standard deviation of 0.002 g/dL. The sensitivity, that indicates the method detection limit, corresponds to the average plus 3 times the standard deviation, and is equal to 0.043 g/dL.

LINEARITY

The reaction is linear up to 12 g/dL. For higher values dilute the serum with saline 0.85%, repeat the dosage and multiply the result obtained per the dilution factor.

DIAGNOSTIC SIGNIFICANCE

The plasmatic levels of the Total Proteins suffer variation according to the changes of its various fractions. Its concentration is influenced by the nutritional state, hepatic, kidney and metabolic errors.

Increased values of the Total Protein can occur in dehydration, multiple myeloma, collagen diseases (erythematous lupus, rheumatoid arthritis), bacterial sub-acute endocarditis.

The hypoproteinemia can happen in the following clinical conditions: nephrotic syndrome, hepatic insufficiency (cirrhosis, chronicle hepatitis, neoplasia), severe innutrition, severe anemia, extended fever state, fever infections, massive bleeding.

The Proteins in serum liquids (pleural and ascitic) is important to difference between transudates (Protein values below 50% of serum concentration) and exudates (Protein values superior to 50% of serum concentration).

The Protein levels in sensorial liquids are inferior to the blood varying from 1.2 a 2.5 g/dL.

BIBLIOGRAPHIC REFERENCES

1 - GORNALL, A. G.; BARDAWILL, C. J.; DAVID, M. M. J., Biol. Chem., 1977, 751.
2 - WEICHSELBAUM, T. E.; AMER,J., Clin. Pathol., 1946, 16,40.
3 - SLATER, L.; CARTER, P. M.; HOBBS, J. R., Ann. Clin. Biochem., 1975, 12,333.
4 - BATSAKIS, J. G.; AROUSOHN, R. S.; WALKER, W. A.; BARNES, B.; AMER, J., Clin. Pathol., 1.976, 66,238.
5 - HOEL. P. G., em Estatística Elementar, Ed. Fundo de Cultura S/A, 1969.
6 - TONKS, D. B., Clin. Chem., 1983, 9,217.
7 - CARL, A. B. and EDWARD, R. A., Tietz Textbook of Clinical Chem. 2nd ed, 1994, 695-697.
8 - WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.
9 - QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUIMICA BASICA Ltda.

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone: +55 31 3439.5454 | E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service
Phone.: 0800 0315454 | E-mail: sac@bioclin.com.br

ANVISA registration for Total Protein Monoreagent kit:
10269360116

Review: August/2022

UNIVERSAL SYMBOLOGY

CATALOG NUMBER

LOT NUMBER

MANUFACTURING DATE

VALIDITY DATE
(last day of the month)

TEMPERATURE LIMIT
(store)

CONTENT IS SUFFICIENT
FOR <N> TEST

SEE INSTRUCTIONS
FOR USE

IN VITRO DIAGNOSTIC
PRODUCT

KEEP AWAY
FROM SUNLIGHT

DO NOT REUSE

CAUTION

MADE BY

CONTROL

POSITIVE CONTROL

NEGATIVE CONTROL

BIOLOGICAL RISK

FLAMMABLE

CORROSIVE

TOXIC

DO NOT USE IF
PACKAGE IS
DAMAGED

PRODUCT
STERILIZED

DANGER