PROTEÍNA C REATIVA

REF K059

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FINALIDADE

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INSTRUCTIONS FOR USE

FUNCTION

Method for quantitative determination of C Reactive Protein (CRP) in biological samples of serum or plasma (Heparin or EDTA). Immunoturbidimetric test, for *in vitro* diagnostic only.

PRINCIPLE OF ACTION

Methodology: Immunoturbidimetry

The reaction permits to quantify by a turbidimetric method, the CRP concentration in the sample. The polystyrene particles coated with anti-CRP mix with the sample forming aggregates in the presence of CRP. The agglutination process causes an increase of particle size and consequently an increase in absorbance, which is measured by comparison with the calibrator of known concentration.

REAGENTS

Number 1 - Buffer - Store between 2 and 8°C. Contains: Buffer < 200 mmol/L, stabilizers, surfactant and preservative.

Number 2 - CRP Latex - Store between 2 and 8°C. Contains: Suspension of Latex particles sensitized with anti-CRP antibodies < 10 µg/dL, stabilizers and preservative.

Number 3 - Calibrator - Store between 2 and 8°C. Contains: C Reactive Protein, stabilizers and preservative. **Potentially infectious.**

Warning: Concentrations of C Reactive Protein vary according to the lot - See bottle label.

PRESENTATION

Presentations	Reagent		
	Nº 1	Nº 2	Nº 3
1	18 mL	2 mL	1,0 mL
2	45 mL	5 mL	1,0 mL
3	2 x 45 mL	2 x 5 mL	2 x 1,0 mL
4	20 mL	5 mL	1,0 mL
5	2 x 20 mL	2 X 5 mL	1,0 mL
6	3 x 20 mL	3 x 5 mL	2 x 1,0 mL
7	4 x 20 mL	4 x 5 mL	2 x 1,0 mL
8	40 mL	10 mL	1,0 mL
9	2 x 40 mL	2 x 10 mL	2 x 1,0 mL
10	3 x 40 mL	3 x 10 mL	3 x 1,0 mL
11	4 x 40 mL	4 x 10 mL	4 x 1,0 mL
12	60 mL	15 mL	2 x 1,0 mL
13	2 x 60 mL	2 x 15 mL	3 x 1,0 mL

EQUIPMENTS AND OPERATIONAL INPUTS

Biochemistry automated equipment, Multical and Multicontrol Bioclin. These itens are found at markets specialized on Laboratories of Clinical Analysis.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be between 2 to 8°C. The transport at temperatures up to 30°C should not exceed 5 days. Protect from light and avoid moisture. **Do not freeze.**

SPECIAL CARE

- 1- For *in vitro* diagnostic use only.**
- Strictly follow the methodology proposed to obtain exact results.
- Water used in material cleaning must be recent and free of contaminants.
- Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers that may alter the results significantly.
- Reagent Nº 3 should be handled cautiously, since it is subject to biological contamination.
- Handle with care, as the reagent contains Sodium Azide, which is irritant to skin and mucous.
- The biological materials were tested for HIV and HBsAg using state of the art methods and found negative results. The risk of infection cannot be excluded and the reagent should be handled as cautiously as observed for the patient's serum. **Potentially Infectious.**
- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.
- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.
- Do not use the product in case of damaged packaging.
- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Serum or plasma (Heparin or EDTA) obtained free of hemolysis and intense lipemia. Samples are stable between 2 and 8°C for up to 8 days or 3 months at -20°C. Samples with presence of fibrin should be centrifuged before testing.

PROCESS DESCRIPTION

The calibration stability of the C Reactive Protein kit installed on refrigerated equipment is 22 days. This stability may vary depending on the conditions of the test, equipment and environment. Therefore, it is suggested to follow the product performance using control serum.

TECHNIQUE

The reagents are ready to use. Homogenize well the Reagent Nº 2 before use.

For calibration of the reaction, utilize preferentially the kit's internal calibrator or the Bioclin Multical kit. For control use the Bioclin Multicontrol Kit.

This kit is for use in automated biochemical analyzers only. **Check the application sheet for this equipment on the website www.bioclin.com.br or through SAC 0800 031 5454.**

PROCEDURE LIMITATIONS

Hemolysis, icterus and lipemia low or moderate do not interfere in test performance. Gross lipemia and turbid samples should be processed with a greater dilution. Values below 2 mg/L lead to few reproductive outcomes.

INTERFERENT

No interference was observed for Bilirubin up to 60 mg/dL, Lipids up to 900 mg/dL Hemoglobin up to 500 mg/dL and Rheumatoid Factor up to 300 IU/mL.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

TRACEABILITY

The kit's calibrator is traceable to the reference material NIBSC 85/506 (Human C Reactive Protein 1st International Standard - WHO International Standard).

REFERENCE VALUES

The reference values for this method were obtained through the determination of C-Reactive Protein in healthy populations of male and female. Less than 8 mg/L. These values should be used as guidelines, each laboratory should establish its range of reference values, according to the population served. The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

PRODUCT PERFORMANCE ACCURACY

Recovery

The recovery analysis was performed with 05 determinations of samples. Accuracies were calculated and were found in good agreement with the reference values, obtaining a recovery between 97 and 106%.

Comparison of Methods

The C Reactive Protein kit was compared with other commercially available method for dosage of C Reactive Protein. 41 analyzes were performed and the results were evaluated. The linear equation obtained was Y = 0.9998X - 1.2981 and correlation coefficient 0.9976. With these results we can conclude that the kit shows good methodological specificity.

PRECISION

Repeatability

The repeatability was calculated from 40 successive determinations, using 3 samples with different concentrations, obtaining the following results:

Repeatability	Sample 1	Sample 2	Sample 3
Concentration (mg/L)	9,47	17,22	77,82
Standard Deviation (mg/L)	0,08	0,10	0,51
Coefficient of Variation (%)	0,81	0,60	0,65

Reproducibility

The reproducibility was calculated from 40 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

Reproducibility	Sample 1	Sample 2	Sample 3
Concentration (mg/L)	9.47	16.73	77.15
Standard Deviation (mg/L)	0.09	0.40	0.70
Coefficient of Variation (%)	0.96	2.39	0.90

SENSITIVITY

Sensitivity was calculated from 40 determinations of a sample free of C Reactive Protein. The average found was 0.06 mg/L, with a standard deviation of 0.05 mg/L. The sensitivity, which indicates the Method Detection Limit, corresponds the average plus 3 times the standard deviation, and is equal to 0.21 mg/L.

LINEARITY

The reaction is linear to the highest point on the curve. For samples with higher values, it is recommended to dilute the sample with Sodium Chloride 0.85%, repeat the dosage and multiply the result by the dilution factor.

HIGH DOSE HOOK EFFECT

Hook effect wasn't observed with a high concentration of PCR up to 400 mg/L.

DIAGNOSTIC SIGNIFICANCE

The C-reactive protein is a useful indicator of ongoing inflammatory process, either of infectious origin (pneumonia, tuberculosis) or noninfectious (active rheumatic fever, rheumatoid arthritis, lupus erythematosus). Is also present in several other pathological conditions, such as acute myocardial infarction, neoplasm, intense trauma, viruses, burns. The determination of its plasma concentration is an effective test to monitor the therapy used and the prognosis of inflammation.

BIBLIOGRAPHIC REFERENCES

1 - PESCE, A. J.; KAPLAN, L. A.: Methods in Clinical Chemistry, C. V. Mosby Company, 1987.
2 - HELLSING, K., Profides in the Biological Fluids, 1973; 23, 579.
3 - BLOM, M. and HJOME, H., Clinical Chemistry, 1976; 22, 657.
4 - HILLS, L. P. and TIFFANY, T. I., Clinical Chemistry, 1980; 26, 1459.
5 - GALVIN, J. P. Et al Clin. Lab Assays 73, 4th, 1983.
6 - PASSING, H.; BABLOCK, W., J. Clin. Chem. Clin. Biochem 21, 709; 1983.
7 - NOKAMA, Y.; NAKAMURA, R. M., J. Clin. Lab. Anal., 1:15; 1987.
8 - OKOMURA, J. M.; MIYAGI, J. M.; TERADA, K. J., Clin. Lab. Anal. 4, 231; 1990.
9 - QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaig, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca

CEP 31565-130 - Belo Horizonte - MG - Brasil

Phone.: +55 31 3439.5454 | E-mail: bioclin@bioclin.com.br

CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service

Phone.: 0800 0315454 | E-mail: sac@bioclin.com.br

ANVISA registration number for C Reactive Protein kit:

10269360070

Review: May/2023

UNIVERSAL SYMBOLOGY

REF

CATALOG NUMBER

LOT

LOT NUMBER

MANUFACTURING DATE

VALIDITY DATE
(last day of the month)

TEMPERATURE LIMIT
(store)

CONTENT IS SUFFICIENT
FOR <N> TEST

SEE INSTRUCTIONS
FOR USE

IVD

IN VITRO DIAGNOSTIC
PRODUCT

KEEP AWAY
FROM SUNLIGHT

DO NOT REUSE

CAUTION

MADE BY

CONTROL

CONTROL

CONTROL +

POSITIVE CONTROL

CONTROL -

NEGATIVE CONTROL

BIOLOGICAL RISK

FLAMMABLE

CORROSIVE

TOXIC

DO NOT USE IF
PACKAGE IS
DAMAGED

STERILE R

PRODUCT
STERILIZED

DANGER