SÍFILIS BIO

REF K181

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FINALIDADE

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INSTRUCTIONS FOR USE

FUNCTION

Immunochromatographic method for rapid test for qualitative determination of total anti-*Treponema pallidum* antibodies (IgG, IgM and IgA) in serum, plasma or whole blood samples.

PRINCIPLE OF ACTION

Methodology: Immunochromatography.
 Syphilis Bio is an immunochromatographic assay (rapid test) for the qualitative detection of total antibodies (IgG, IgM and IgA) anti-*Treponema pallidum* in serum, plasma or whole blood. In this kit, recombinant *T. pallidum* antigens are immobilized in the region of the test line. When a sample is added, it reacts with colored particles conjugated to antigens of *T. pallidum*. This complex migrates along the test strip and interacts with the immobilized antigens. If the sample shows anti-*Treponema pallidum* antibodies, a colored line will appear in the test region indicating a reagent result. If the sample does not present anti-*Treponema pallidum* antibodies, this line will not appear indicating a non-reactive result.

REAGENTS:

Number 1 (R1) Cassette - Store between 2 and 30°C. **Do not freeze.** Contains:

A) 01 plastic cassette with 01 test strip containing:
 - Control region (C): Anti-*T. pallidum* polyclonal antibody;
 - Test region (T): Recombinant *T. pallidum* antigen;
 - Conjugate: Recombinant antigen conjugate of *T. pallidum* and colloidal gold;

B) 01 sacht of silica;

Number 2 (R2) Diluent - Store between 2 and 30°C. **Do not freeze.** Contains: Buffer solution and preservative.

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PRESENTATION

Reagent	Number of Cassettes for Package	Diluent
Cassette + Diluent (Pipette included) *	5	1.0 mL
	10	1.5 mL
	15	2.0 mL
	20	3.0 mL
	25	3.5 mL
	30	4.0 mL
	40	5.0 mL
	50	6.0 mL
Casette + Diluent (Lancet and pipette included) **	100	2 x 6.0 mL
	5	1.0 mL
	10	1.5 mL
	15	2.0 mL
	20	3.0 mL
	25	3.5 mL
	30	4.0 mL
	40	5.0 mL
	50	6.0 mL
	100	2 x 6.0 mL

* These presentations are accompanied by plastic pipettes disposables proportionally to the number of cassettes.

** These presentations are accompanied by lancets and plastic pipettes disposables, in proportion to the number of cassettes.

EQUIPMENTS AND OPERATIONAL INPUTS

Pipettes, tips, watches and stopwatches, alcohol 70% (p/p), cotton and waste for disposal. They are in the specialized market of articles for Clinical Analysis Laboratories. The kit allows reading without the aid of equipment.

TRANSPORTATION AND STORAGE CONDITIONS

Transport can be carried out at temperatures up to 45°C by maximum, 3 weeks. The storage temperature should be from 2 to 30°C. Avoid moisture. **Do not freeze.**

SPECIAL CARE

- 1- For *in vitro* diagnostic use only.
- 2- Strictly follow the methodology proposed for obtaining results.
- 3- The cassettes (test devices) do not contain contaminating substances, but since the samples are biological, the safety routine must be followed strictly and the use of disposable gloves is essential.
- 4- Open the aluminum envelop, to remove the cassette (test devices), only at the time of use.
- 5- The kit should not be frozen or exposed to high temperatures (above 45°C), as this will cause deterioration of the kit.
- 6- The reagents, as well as the samples, shold reach room temperature before starting the test.
- 7- Do not mix reagents with the same lot or with different lot numbers.
- 8- Do not use kit beyond its expiration date. The expiry date is stated on the product label.
- 9- The cassette (test device), the disposable plastic pipette and the disposable lancet and can not be reused.

10- The kit materials that had contact with the sample of the patient must have a special destination. Disposal of lancets should be made into sharps. The pipettes, the cassettes (test devices) and diluent vials should be disposed of in organic waste.

11- Secondary packaging (box), Instruction for use, envelope aluminum and silica sachet can be disposed of in trash common.
 12- We recommend to apply the local, state and federal environmental protection so that the disposal of reagents and biological material is done in accordance with the current legislation.

13- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or through request by the Customer Service Department (SAC) of the Quibasa.

14- Do not use the product in case of damage to the packaging.
 15- It is essential that instruments and equipment are properly calibrated and periodic maintenance.

SAMPLES

Serum or Plasma: Samples may be stored under refrigeration, between 2 and 8°C, for a maximum period of 3 days. If samples can not be analyzed within 3 days, can be stored for up to 30 days at -20°C (freezer). To obtain the plasma, collect data using such as EDTA, Heparin or Citrate.

Whole Blood: May be collected by venipuncture using adequate anticoagulant (EDTA, Heparin or Citrate) or digital puncture with the aid of a lancet. Whole blood should be used for immediate testing. Whole blood collected by venipuncture with anticoagulant may be stored between 2 and 8°C for up to 7 days. **Do not freeze.**

DESCRIPTION OF PROCESS
 TECHNIQUE – VENOPUNCTION

- 1- The sample should be at room temperature before starting the test.
- 2- Remove the cassette (test device) from the protective packaging, place it on a clean and level surface and identify it properly.
- 3- For Serum, Plasma or Whole Blood: Transfer 10 µL of serum, plasma or whole blood in sample well, smaller well (A). Hold the Diluent (Reagent N° 2) vertically and apply 2 drops (70 µL) in diluent well, bigger well (B).
- 4- Wait for the lines to form. Interpret the results between 15 and 30 minutes. **Do not interpret after 30 minutes.**

TECHNIQUE - DIGITAL PUNCTURE / LANCET


- 1- Remove the cassette (test device) from the protective packaging, place it on a clean and level surface and identify it properly.
- 2- Wash and dry the patient's hands. Pass alcohol 70% (p/p) at the tip of the finger for antisepsis of the area used.
- 3- Pressing the finger so that there is accumulation of blood in the tip of the finger that will be perforated by the lancet.
- 4- Remove the protective cap from the lancet.
- 5- Position and press the lancet firmly on the area to be thrown. Then the blood will come out through the perforated area.
- 6- Collect the blood with a disposable plastic pipette that accompanies the kit. Tighten the plastic pipette above the marked trace and touch the open cavity of the plastic pipette in the drop of blood.

Relieve the pressure in the plastic pipette so the blood is aspirated. Aspirate the blood until the trace marked on the disposable plastic pipette (equivalent to 10 µL of blood).

7- In sample well, smaller well (A) dispense the blood by pressing the pipette again.


8- Hold the Diluent bottle (Reagent N° 2) vertically and apply 2 drops (70 µL) in diluent well, bigger well (B).

9- Wait for the lines to form. Interpret the results between 15 and 30 minutes. **Do not interpret after 30 minutes.**

INTERPRETATION OF RESULTS


Reagent Test: Formation of two red lines, one in control region (C) and one in the test region (T) after 15 minutes. **Do not interpret after 30 minutes.**

Non-Reagent Test: Formation of a red line in the control region (C) and complete absence of line in the test region (T) after 15 minutes. **Do not interpret after 30 minutes.**

Invalid Testing: Absence of line formation in the region (C) indicates an error in the procedure or deterioration of the cassette. In this case, repeat the test using a new cassette.

Notes:

- 1- The line in the control region will appear before 15 minutes test incubation, and this does not mean that the results can be interpreted ahead of schedule.
- 2- It is important to use the correct sample volume because lower or higher volumes can determine results mistakes.
- 3- The results provided by this kit should be interpreted responsible professional, and is not the sole criterion for the determination of the diagnosis and/or treatment of the patient.
- 4- The isolated test results can not be used to definitive diagnosis. For conclusion of diagnosis follow, the guidelines of Administrative Rule N°. 2.012 of October 19 2016, approving the Technical Manual for the Syphilis Diagnosis.
- 5- The Syphilis Bio kit detects antibodies against *Treponema pallidum* at all stages of the disease, without being restricted to a stage specific.
- 6- In primary syphilis, antibodies begin to appear in the the bloodstream about 7 to 10 days after the of hard cancer. However, it is important to note that time for seroconversion may vary from individual to individual. Therefore, a Non-Reagent Result in the Syphilis Bio Kit Does Not Preclude Infection by *T. pallidum*. If a suspicion of infection persists, a new sample should be collected after 30 days for repeat testing.
- 7- Serological tests for syphilis are classified as non - treponemics and treponemics. The kit Sifilis Bio - Bioclin is a treponemic test.

INTERFERENT

Hemolysis with hemoglobin values up to 10 g/L interfere with the dosage of syphilis.
No significant interferences were found up to 30 g/L of lipids, 0.2 g/L of bilirubin and 60 g/L of proteins.
No cross-reactivity with HBsAg, HCV and HIV.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal program of quality control, where procedures, standards, limits and tolerance for variations are clearly established. It is important to note that all measuring systems characteristic analytical variability, which should be monitored by the laboratories themselves. Therefore, it is the use of controls, which allow the assessment of the precision and accuracy of dosages.

PRODUCT PERFORMANCE

ACCURACY

Repeatability

The repeatability was calculated from 10 successive determinations, using 2 differents samples, obtaining the following results:

Samples	Nº of Repetitions	Expected Results	Found Results
01	10	Positive	Positive
02	10	Positive	Positive

Reproducibility

The reproducibility was calculated from 10 successive determinations for 3 consecutive days, obtaining the following results:

Samples	Nº of days	Nº of Repetitions	Expected Results	Found Results
01	03	10	Positive	Positive
02	03	10	Positive	Positive

SENSITIVITY AND CLINICAL SPECIFICITY

The Syphilis Bio kit was tested using clinically known samples. For this test, 990 samples were analyzed, of which 396 reagents for Syphilis and 594 non-reagents for Syphilis. We found the following results:

	Total Expected	Syphilis Bio
Sample Positive	396	396
Sample Negative	594	593
Total samples tested	990	

Clinical Sensitivity: > 99.9% (396/396)
Specificity: 99.8% (594/593)

DIAGNOSTIC SIGNIFICANCE

Syphilis is an infectious disease caused by a spirochete called *Treponema pallidum* that evolves slowly in three stages, characterized by lesions of the skin and mucous membranes. You can be transmitted through sexual contact, thus configuring sexually transmitted disease and, more rarely, for fetal-placental contamination (congenital syphilis).
The first manifestations occur after a period of average incubation time of 21 days, and may vary from 3 up to 90 days.

The disease presents three distinct phases, with manifestations characteristics in each and a period of latency (without symptoms) between the second and third phases.
In the primary phase, the classic lesion occurs in the genitals, the entrance of *Treponema*, called hard chancre. The chancre leads in average 3 to 6 weeks to heal, and may not leave a mark some.
The secondary phase occurs after 4 to 8 weeks of onset of chancre, and this lesion may still be present.
At this stage, the greatest amount of treponemes occurs in the body.
The tertiary phase is the phase of slow progressive (chronic) inflammation, with symptoms related to organs predominantly committed, it is destructive and incapacitating.

BIBLIOGRAPHIC REFERENCES

- 1- GRADWOHL'S Clinical Laboratory Methods and Diagnosis - Sonnen Wriith and Jerret, Mosby Cd., 1980.
- 2- Manual de Reacciones Para el Diagnóstico de la Sífilis - nº 311 - Organización Mundial de Saúde.
- 3- Ferreira, A.W.; Diagnóstico Laboratorial das Principais Doenças Infecciosas e auto-imunes, 2 ed., Editora Guanabara Koogan SA, 2001.
- 4- Braverman PK: Herpes, Syphilis and other Ulcergenital Conditions, Adolesc Med. 1996; 7 (1): 93-118.
- 5- George R, Pope V, Fears M, Morrill B, Larsen S: An analysis of the value of some antigen-antibody interactions used as diagnostic indicators in a treponemal Western Blot (TWB) test for syphilis. J Clin Lab Immunol. 1998; 50 (1):27-44.
- 6- Schachter J: Office and Laboratory Testing for STDs: Present and future. Adolesc Med. 1990 Oct;1(3):583-596.
- 7- Portaria Nº 2.012, de 19 de Outubro de 2016 – Aprova o Manual Técnico para Diagnóstico da Sífilis e dá outros providências.
- 8- QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda
Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone.: +55 (31) 3439.5454 | E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service
Phone.: 0800 0315454 | E-mail: sac@bioclin.com.br

ANVISA registration for kit Syphilis Bio: 10269360308

Review: March/2023

UNIVERSAL SYMBOLOGY

CATALOG NUMBER

LOT NUMBER

MANUFACTURING DATE

VALIDITY DATE
(last day of the month)

TEMPERATURE LIMIT
(store)

CONTENT IS SUFFICIENT
FOR <N> TEST

SEE INSTRUCTIONS
FOR USE

IN VITRO DIAGNOSTIC
PRODUCT

KEEP AWAY
FROM SUNLIGHT

DO NOT REUSE

CAUTION

MADE BY

CONTROL

POSITIVE CONTROL

NEGATIVE CONTROL

BIOLOGICAL RISK

FLAMMABLE

CORROSIVE

TOXIC

DO NOT USE IF
PACKAGE IS
DAMAGED

PRODUCT
STERILIZED

DANGER