SILICATO

REF F007

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FINALIDADE

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- 1- ASTM - Standard Specification for Reagents Water. ASTM document D 1193-91(1991).
- 2- NCCLS-National Committee for Clinical Laboratory Standards: Preparation and Testing of Reagent Water in the Clinical Laboratory - Third Edition - 1997, vol. 17 nr.18.
- 3- Kaplan, L.A.; Pesce, A.J. - Methods in Clinical Chemistry- Second Edition. St. Louis, Missouri, C.V. Mosby Co., 1989.
- 4- Silicates in Water and Wastewater, Industrial Method No. 105-71W, Technicon Instruments Corporation, Tarrytown, NY (1973).
- 5- MANNO, B.R.; ABUKHALAF, I.K.; MANNO, J.E. A Sensitive Spectrophotometric Assay for the Measurement of Soluble Silica in Water - Journal of Analytical Toxicology, Volume 21, Number 6, October 1997, pp. 503–505.
- 6- QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

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- 1- ASTM - Standard Specification for Reagents Water. ASTM document D 1193-91(1991).
- 2- NCCLS-National Committee for Clinical Laboratory Standards: Preparation and Testing of Reagent Water in the Clinical Laboratory - Third Edition – 1997, vol. 17 nr.18.
- 3- Kaplan, L.A.; Pesce, A.J. - Methods in Clinical Chemistry- Second Edition. St. Louis, Missouri, C.V. Mosby Co., 1989.
- 4- Silicates in Water and Wastewater, Industrial Method Nº. 105-71W, Technicon Instruments Corporation, Tarrytown, NY (1973).
- 5- MANNO, B.R.; ABUKHALAF, I.K.; MANNO, J.E. A Sensitive Spectrophotometric Assay for the Measurement of Soluble Silica in Water - Journal of Analytical Toxicology, Volume 21, Number 6, October 1997, pp. 503–505.
- 6- QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

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INSTRUCTIONS FOR USE

FUNCTION

Method for qualitative determination of silicate dissolved in water samples, through colorimetric reactions. Manual appliance.

PRINCIPLE OF ACTION

Methodology: Colorimetric

The soluble silica reacts with molybdate in acid forming silica-molybdate complex of yellow color, whose intensity is proportional to the concentration of silica in the sample. Product absorbance is measured at wavelengths between 405 and 410 nm.

REAGENTS

Number 1 - Acidifier - Store between 15 and 30°C. Contains Hydrochloric Acid $\leq 6,25N$. **Corrosive.**

Number 2 - Catalyst - Store between 15 and 30°C. Contains Molybdate ≤ 81 mmol/L (pH 6,0).

PRESENTATION

Reagent	Volume
Nº 1	3.5 mL
Nº 2	3.5 mL

EQUIPMENTS AND OPERATIONAL INPUTS

Semi-automated biochemical devices, watches, stopwatches, test tubes and kits. Can all be found at stores specialized in Laboratories of Clinical Analysis.

TRANSPORTATION AND STORAGE CONDITIONS

The storage and transportation temperatures should be between 15 to 30°C. Protect from light and avoid moisture.

SPECIAL CARE

1- Reagent Nº 1 is corrosive, therefore handle with care. Reagent Nº 2 may occur alterations (color shift, formation of precipitates) not indicating contamination or improper for use.

2- Strictly follow the methodology proposed to obtain exact results.

3- Water used in material cleansing must be recent and free of contaminants.

4- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

5- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS

(Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.

SAMPLES

Deionized water sample must be collected at the moment of determination.

PROCESS DESCRIPTION

TECHNIQUE

Select wavelength between 405 nm to 410 nm. Take 1 (one) test tube and proceed accordingly:

Deionized water	5.0 mL
Reagent Nº 1	2 drops
Reagent Nº 2	2 drops

Homogenize and wait for 2 minutes.

Read the sample absorbance against the deionized water blank. The color is stable for 20 minutes.

QUALITY CONTROL OF WATER

Should be a routine practice of Clinical Laboratory the quality control of deionized water used in the laboratory. One of the first contaminants emitted from the ion exchange resin, when it approaches saturation, is Silica (Silicate). Then highly ionized substances, which increases the conductivity and, consequently, reduces the resistivity are released.

It is suggested that the water used in the laboratory must have conductivity < 1 microsiemens/cm, resistivity > 1 megaohm.cm and silicate concentration $< 0,10$ mg/L. A high concentration of Silicate in deionized water may indicate that the ion exchange resin deionizer already reached the saturation point and, therefore, it is necessary to regenerate its resin beds.

REFERENCE VALUES

The absorbance of the reaction, read between 405 and 410 nm, should not exceed 0.010, which indicates Silicate concentrations above 0.10 mg/L. If the absorbance is exceeded, it is suggested the regeneration of resins or exchange of the deionizer columns. If water is obtained by another process, it must undergo an evaluation to detect the reason for the high level of silica. The formation of visible yellow color of the reaction indicates higher silicate concentration than acceptable levels.

PRODUCT PERFORMANCE
SENSITIVITY

Was calculated from 20 determinations of a sample with 0 (zero) concentration of silicate. The average found of 0.0016 mg/L with standard deviation of 0.0007 mg/L. Sensitivity, that indicates the limits of detection from the method, corresponds the average plus 3 times the standard deviation and is equal to 0.0037 mg/L.

BIBLIOGRAPHIC REFERENCES

- 1- ASTM - Standard Specificacion for Reagents Water. ASTM document D 1193-91(1991).
- 2- NCCLS-National Committee for Clinical Laboratory Standards: Preparation and Testing of Reagent Water in the Clinical Laboratory - Third Edition – 1997, vol. 17 nr.18.
- 3- Kaplan, L.A.; Pesce, A.J. - Methods in Clinical Chemistry- Second Edition. St. Louis, Missouri, C.V. Mosby Co., 1989.
- 4- Silicates in Water and Wastewater, Industrial Method No. 105-71W, Technicon Instruments Corporation, Tarrytown, NY (1973).
- 5- MANNO, B.R.; ABUKHALAF, I.K.; MANNO, J.E. A Sensitive Spectrophotometric Assay for the Measurement of Soluble Silica in Water - Journal of Analytical Toxicology, Volume 21, Number 6, October 1997, pp. 503–505.
- 6- QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda
Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone: +55 31 3439.5454 | E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service
Phone.: 0800 0315454 | E-mail: sac@bioclin.com.br

Product not eligible for registration with ANVISA.

Review: October/2023

UNIVERSAL SYMBOLOGY

CATALOG NUMBER

LOT NUMBER

MANUFACTURING DATE

VALIDITY DATE
(last day of the month)

TEMPERATURE LIMIT
(store)

CONTENT IS SUFFICIENT
FOR <N> TEST

SEE INSTRUCTIONS
FOR USE

IN VITRO DIAGNOSTIC
PRODUCT

KEEP AWAY
FROM SUNLIGHT

DO NOT REUSE

CAUTION

MADE BY

CONTROL

POSITIVE CONTROL

NEGATIVE CONTROL

BIOLOGICAL RISK

FLAMMABLE

CORROSIVE

TOXIC

DO NOT USE IF
PACKAGE IS
DAMAGED

PRODUCT
STERILIZED

DANGER