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## WORKSTATION

REF L009

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## CLEANSING SOLUTION WORKSTATION

TO OBTAIN THE INSTRUCTIONS FOR USE IN PRINTED FORMAT,  
AT NO ADDITIONAL COST, CONTACT CUSTOMER ADVISORY SERVICE:

SAC: +55 (31) 3439 5454 / 0800 031 5454 / sac@bioclin.com.br

REF **L009**

### INSTRUCTIONS FOR USE

#### FUNCTION

Solution for cleaning automatic ELISA equipment.

#### Suitable for:

\*Family of Automated ELISA Analyzers (BIOLISA WORKSTATION)

#### APPLICATION

The cleaning solution acts as a degreasing and sanitizing agent, removing reagent residues, salts, and proteins that may interfere with equipment performance and the reproducibility of ELISA assay results.

#### REAGENTES

**1 – Cleaning Solution Workstation** – Store between 15 °C and 30 °C. Contains: Isopropyl Alcohol.

#### PRESENTATION

| Presentation | Volume     |
|--------------|------------|
| 1            | 1 x 200 mL |
| 2            | 2 x 200 mL |

#### STORAGE AND TRANSPORT CONDITIONS

The storage and transport temperature should be between 15 °C and 30°C. Keep away from light and avoid moisture.

#### SPECIAL PRECAUTIONS

- 1- For *in vitro* diagnostic use only.
- 2- Strictly follow the proposed methodology to achieve adequate cleaning.
- 3- Avoid contact with eyes and skin. In case of contact, wash immediately with plenty of water.
- 4- Use laboratory gloves when handling the product.
- 5- The product may be harmful if swallowed or inhaled; if swallowed, seek medical attention.
- 6- The Cleaning Solution contains Isopropyl Alcohol. Therefore, handle it with due care.
- 7- For information related to biosafety or in case of accidents with the product, consult the SDS (Safety Data Sheet) available on the website [www.bioclin.com.br](http://www.bioclin.com.br) or by requesting it from Quibasa's Customer Service (SAC).
- 8- Do not use the product if the packaging is damaged.

#### PREPARATION OF THE WORKING REAGENT CLEANSING SOLUTION WORKSTATION

Cleaning Solution Workstation

Dilute the contents of the Workstation Cleaning Solution bottle in 800 mL of distilled or deionized water; for example: 200 mL Workstation Cleaning Solution in 800 mL of distilled or deionized water. After preparation, the solution can be stored at room temperature 15–30°C for up to 30 days.

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#### QUALITY ASSURANCE

Before being released for consumption, all Bioclin reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

#### QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca  
CEP 31565-130 - Belo Horizonte - MG - Brasil  
Phone.: +55 31 3439.5454 | E-mail: [bioclin@bioclin.com.br](mailto:bioclin@bioclin.com.br)  
CNPJ: 19.400.787/0001-07 - Made in Brazil

#### CUSTOMER SERVICE

Customer Advisory Service  
Phone.: 0800 0315454 | E-mail: [sac@bioclin.com.br](mailto:sac@bioclin.com.br)

Product exempt from ANVISA registration, in accordance with current legislation.

**Review:** February/2026

## UNIVERSAL SYMBOLOGY

|   |  |   |  |
|---|--|---|--|
|    | CATALOG NUMBER                           |    | MADE BY                                |
|    | LOT NUMBER                               |    | CONTROL                                |
|    | MANUFACTURING DATE                       |    | POSITIVE CONTROL                       |
|   | VALIDITY DATE<br>(last day of the month) |   | NEGATIVE CONTROL                       |
|  | TEMPERATURE LIMIT<br>(store)             |  | BIOLOGICAL RISK                        |
|  | CONTENT IS SUFFICIENT<br>FOR <N> TEST    |  | FLAMMABLE                              |
|  | SEE INSTRUCTIONS<br>FOR USE              |  | CORROSIVE                              |
|  | IN VITRO DIAGNOSTIC<br>PRODUCT           |  | TOXIC                                  |
|  | KEEP AWAY<br>FROM SUNLIGHT               |  | DO NOT USE IF<br>PACKAGE IS<br>DAMAGED |
|  | DO NOT REUSE                             |  | PRODUCT<br>STERILIZED                  |
|  | CAUTION                                  |  | DANGER                                 |