TIPAGEM SANGUÍNEA FELINA VET FAST

REF VET085

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FINALIDADE

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- 1- FELDMAN, B. F. Blood Transfusion Guidelines. In: WSAVA 2003 Congress Proceedings. Blacksburg, VA: Virginia-Maryland Regional College of Veterinary Medicine, 2003.
- 2- KNOTTENBELT, C. M. The feline AB blood group system and its importance in transfusion medicine. Journal of Feline Medicine and Surgery, v. 4, p. 69-76, 2002.
- 3- SETH, M.; JACKSON, K. V.; GIGER, U. Comparison of five blood-typing methods for the feline AB blood group system. American Journal of Veterinary Research, v. 72, n. 2, p. 203-209, fev. 2011.
- 4- SPADA, E.; PEREGO, R.; BAGGIANI, L. et al. Prevalence of Blood Types

and Alloantibodies of the AB Blood Group System in Non-Pedigree Cats from Northern (Lombardy) and Southern (Sicily) Italy. Animals, v. 10, art. 1129, 03 jul. 2020-1-

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- 1- FELDMAN, B. F. Blood Transfusion Guidelines. In: WSAVA 2003 Congress Proceedings. Blacksburg, VA: Virginia-Maryland Regional College of Veterinary Medicine, 2003.
- 2- KNOTTENBELT, C. M. The feline AB blood group system and its importance in transfusion medicine. Journal of Feline Medicine and Surgery, v. 4, p. 69–76, 2002.
- 3- SETH, M.; JACKSON, K. V.; GIGER, U. Comparison of five blood-typing methods for the feline AB blood group system. American Journal of Veterinary Research, v. 72, n. 2, p. 203–209, fev. 2011.
- 4- SPADA, E.; PEREGO, R.; BAGGIANI, L. et al. Prevalence of Blood Types

and Alloantibodies of the AB Blood Group System in Non-Pedigree Cats from Northern (Lombardy) and Southern (Sicily) Italy. Animals, v. 10, art. 1129, 03 jul. 2020.1-

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INSTRUCTIONS FOR USE

FUNCTION

Rapid immunochromatographic test for the qualitative determination of feline erythrocyte antigens of types A and B in feline whole blood samples.

PRINCIPLE OF ACTION

Methodology:

Immunochromatography
Feline Blood Typing VET FAST is an *in vitro* diagnostic test for the rapid and qualitative detection of feline erythrocyte antigens of types A and B in feline whole blood samples.

The test consists of two independent cassettes, joined by a support (combo cassette): the Feline-A cassette detects Antigens A and the Feline-B cassette detects Antigens B. Each cassette consists of a nitrocellulose membrane, impregnated on the test line with anti-antigen A antibodies (Feline-A cassette) and anti-antigen B antibodies (Feline-B cassette). When the sample is added to the cassettes, it migrates by capillarity through the test strip until it reaches the test line. In this region, the following occurs:

In the Feline-A cassette: red blood cells with antigen A, possibly present in the sample, interact with the anti-antigen A antibody, forming the complex: Antigen A – Immobilized anti-antigen A antibody. As a result, a red line is formed on the test line, indicating a reactive result.

In the Feline-B cassette: red blood cells with antigen B, possibly present in the sample, interact with the anti-antigen B antibody, forming the complex: Antigen B – Immobilized anti-antigen B antibody. As a result, a red line is formed on the test line, indicating a reactive result.

In both cassettes, the excess sample continues to migrate through the strip until it reaches the control line, and a second line is formed (control line), indicating that the test presented a valid result.

REAGENTS

Number 1 (R1) Cassette - Store between 2 and 30 °C. **Do not freeze.** Contains:

- A) Cassette Feline-A
- B) Cassette Feline-B
- C) 01 sachet of silica

Number 2 (R2) Diluent - Store between 2 and 30 °C. **Do not freeze.** Contains: Buffered solution and preservative.

Accessories:

Disposable plastic pipette and Collection Tube with EDTA.
Note: The kit comes with a veterinary report and an identification sticker. It is suggested that you stick the identification sticker on the cat's vaccination card or other document that you deem relevant.

PRESENTATION

Reagent	Presentation	Number of Cassettes per Pack	Number of Diluent per Pack
Cassette	VET085-1	10	10 x 1.0 mL
	VET085-2	25	25 x 1.0 mL
Diluent	VET085-3*	10	10 x 1.0 mL
	VET085-4*	25	25 x 1.0 mL
	VET085-5	5	5 x 1.0 mL

*Individualized presentation.

EQUIPMENTS AND OPERATIONAL INPUTS

Clock or stopwatch, PPE (gloves, mask, protective glasses) and trash for disposal. All listed inputs can be found in the specialized market for articles for Clinical Analysis Laboratories.

STORAGE AND TRANSPORT CONDITIONS

The storage and transport temperature of the product should be between 2 and 30 °C. Keep protected from direct sunlight. **Do not freeze.**

The cassette is sensitive to air humidity. Therefore, it must always be kept inside its primary packaging (aluminum pouch) until the moment of use.

SPECIAL CARES

- 1- Product for *in vitro* diagnosis, for veterinary use only.
- 2- Strictly follow the proposed methodology to ensure accurate results.
- 3- The cassettes (test devices) do not contain contaminating substances,

but since the samples are biological, strict safety protocols must be followed. The use of disposable gloves and other personal protective equipment is mandatory.

- 4- Use suitable containers for storing biological samples.
- 5- Open the aluminum pouch to remove the cassette (test device) only at the time of use.
- 6- The kit must not be frozen or exposed to high temperatures, as this will cause deterioration of its components.
- 7- Reagents and samples must reach room temperature before starting the test.
- 8- Do not use reagents from different batches or kits.
- 9- Do not use the kit past its expiration date. The expiration date is indicated on the product label.
- 10- All components of the kit that come into contact with the sample are disposable and must not be reused.
- 11- Components that have been in contact with the patient's sample and the diluent included in the kit must be discarded as biological waste.
- 12- Secondary packaging (box), Instructions for Use, aluminum pouch, and silica gel packet can be disposed of in regular waste.
- 13- Follow local, state, and federal environmental protection regulations to dispose of reagents and biological material in compliance with current legislation.
- 14- For biosafety information or in case of accidents with the product, consult the SDS (Safety Data Sheet) available on the website www.bioclin.com.br or request it via Quibasa's Customer Support Service.
- 15- Do not use the product if the packaging is damaged.







SAMPLES

Whole Blood

Collect the blood sample into a tube containing EDTA and Sodium Citrate following standard laboratory procedures. Whole blood can be stored between 2 and 8 °C and tested within 24 hours after collection. **Do not freeze.** Samples with clots, lipemia, hemolysis or other abnormal situations should not be used.


PROCESS DESCRIPTION


- 1- The sample must be at a temperature between 15 and 30 °C before starting the test.
- 2- Remove the cassette (test device) from its protective packaging, place it on a clean, level surface, and label it appropriately.
- 3- Make sure the diluent (R2) is at the bottom of the vial and remove the dropper stopper.
- 4- Collect the sample using the disposable plastic pipette provided in the kit.
- 5- Press the bulb of the pipette and touch it to the sample. Release the pressure on the bulb so that the sample is aspirated.
- 6- To dispense the sample, press the bulb again and apply 1 drop of the sample (approximately 50 µL) into the vial containing the diluent.
Note: Alternatively, an automatic pipette can be used to add the sample.
Note: The diluted sample must be tested within 20 minutes of dilution. Otherwise, the sample must be diluted again.
- 7- Insert the dropper stopper into the diluent bottle.
- 8- Homogenize the sample by inverting the bottle for 5 seconds.
- 9- Remove the bottle cap and slowly add three drops of the mixture (diluent + sample) (75 µL) to each sample well in both cassettes (Feline-A and Feline-B).
- 10- Wait for the lines to form. Interpret the results in 10 minutes. **Do not interpret after 15 minutes.**

INTERPRETATION OF THE RESULT

<p>Reactive Test: Formation of a red line in the test region (T) and another line in the control region (C). Do not interpret after 15 minutes. Attention: Presence of white line is a non-reactive result.</p>	
<p>Non-Reactive Test: Formation of a red line in the control region (C) and complete absence of a red line in the test region (T). Do not interpret after 15 minutes. Attention: Presence of white line is a non-reactive result.</p>	
<p>Invalid Result: Complete absence of a line in the control region (C), with or without a red line in the test region (T). Test the sample again.</p>	

B) Comparison with a commercial kit for Feline Blood Typing

A total of 83 samples characterized by another commercial kit for Feline Blood Typing were analyzed based on the different blood types, namely type A, type B and type AB. For this test, 49 samples of type A, 22 samples of type B and 12 samples of type AB were analyzed. The results show that the clinical sensitivity of the product is > 99% and the clinical specificity is > 99%.

METHOD	REFERENCE - Commercial kit				Total
	Result	Type A	Type B	Type AB	
Feline Blood Typing VET FAST	Type A	49	0	0	49
	Type B	0	22	0	22
	Type AB	0	0	12	12
	Result Total	49	22	12	83

Clinical Sensitivity: > 99% - 95% CI 95% - 100%

Clinical Specificity: > 99% - 95% CI 95% - 100%

Accuracy: 100%

PRECISION

Repeatability

Repeatability was calculated from 10 successive determinations, using 3 different positive samples, 1 type A sample, 1 type B sample and 1 type AB sample. The following results were obtained:

Sample	N° of Repetitions	Expected result	Result found
Type A sample	10	Reactive to A	Reactive to A
Type B sample	10	Reactive to B	Reactive to B
Type AB sample	10	Reactive to AB	Reactive to AB

Reproducibility

Reproducibility was calculated from 10 successive determinations over 3 consecutive days, using 3 different positive samples, 1 type A sample, 1 type B sample and 1 type AB sample. The following results were obtained:

Sample	N° of Days	N° of Repetitions	Expected result	Result found
Type A sample	03	10	Reactive to A	Reactive to A
Type B sample	03	10	Reactive to B	Reactive to B
Type AB sample	03	10	Reactive to AB	Reactive to AB

CLINICAL SIGNIFICANCE

The feline AB blood system includes three groups: A, B, and AB. Type B cats have anti-A alloantibodies, while most type A cats have anti-B alloantibodies, making transfusions between incompatible groups high-risk.

Transfusion compatibility:

- AB cats can receive either A or AB blood
- B cats should receive only B blood
- A cats should preferably receive A blood

Incompatible transfusions can cause potentially fatal acute hemolysis or reduce red blood cell survival to less than 48 hours. Therefore, it is recommended to type the donor and recipient before the first transfusion and perform crossmatching in subsequent transfusions.

In emergency situations, rapid typing increases therapeutic success and avoids wasting blood components. In breeding, prepartum typing prevents neonatal isoerythrocyte lysis in litters of B females mated with A or AB males. Routine use of this test significantly reduces morbidity and mortality associated with feline transfusions and neonatal isoerythrolysis.

BIBLIOGRAPHIC REFERENCE

1- FELDMAN, B. F. Blood Transfusion Guidelines. In: WSAVA 2003 Congress Proceedings. Blacksburg, VA: Virginia-Maryland Regional College of Veterinary Medicine, 2003.

2- KNOTTENBELT, C. M. The feline AB blood group system and its importance in transfusion medicine. Journal of Feline Medicine and Surgery, v. 4, p. 69-76, 2002.

3- SETH, M., JACKSON, K. V.; GIGER, U. Comparison of five blood-typing methods for the feline AB blood group system. American Journal of Veterinary Research, v. 72, n. 2, p. 203-209, fev. 2011.

4- SPADA, E.; PEREGO, R.; BAGGIANI, L. et al. Prevalence of Blood Types and Alloantibodies of the AB Blood Group System in Non-Pedigree Cats from Northern (Lombardy) and Southern (Sicily) Italy. Animals, v. 10, art. 1129, 03 jul. 2020-1-

5- QUIBASA: Data from the Research and Development Department.

QUALITY GUARANTEE

Before being released for consumption, all Bioclin reagents are tested by the Quality Control Department. The quality of the reagents is assured until the validity date mentioned on the packaged, as long as they are stored and transported under the appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda
Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Tel.: (31) 3439.5454 | E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

CONSUMER SERVICE

Customer Support Service
Tel.: 0800 0315454 | E-mail: sac@bioclin.com.br

Product Exempt from Registration at the Ministry of Agriculture, Livestock and Supply.

Technical Responsible: Dr. Camila Eckstein (CRMV/MG 20.611).

Review: January/2026

UNIVERSAL SYMBOLOGY

	CATALOG NUMBER		MADE BY
	LOT NUMBER		CONTROL
	MANUFACTURING DATE		POSITIVE CONTROL
	VALIDITY DATE (last day of the month)		NEGATIVE CONTROL
	TEMPERATURE LIMIT (store)		BIOLOGICAL RISK
	CONTENT IS SUFFICIENT FOR <N> TEST		FLAMMABLE
	SEE INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC PRODUCT		TOXIC
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED
	DO NOT REUSE		PRODUCT STERILIZED
	CAUTION		DANGER

Notes:

1- The line in the control region may appear before the 10 minutes of incubation; however, this does not mean results can be interpreted earlier than the stipulated time.

2- Results must not be interpreted after 15 minutes to avoid incorrect interpretations.

3- Results provided by this kit must be interpreted by a licensed veterinary professional and should not be the sole criterion for diagnosing or treating the patient.

4- The reddish coloration of the reading window occurs because this kit uses a whole blood sample, and the migration of red blood cells through the reading window is essential for the detection of antigens A and B through the binding of feline red blood cells to the test and control lines. Therefore, the presence of red blood cells gives the window a reddish tint without affecting the interpretation of the results.

TEST LIMITATIONS

1- Using the correct sample volume is crucial, as smaller or larger volumes may produce erroneous results.

2- Results may be compromised in samples with very low hematocrit (<15%), intense hemolysis, jaundice or lipemia.

3- Intense autoagglutination (possible in cases of immune-mediated hemolytic anemia) may block migration and invalidate the test.

4- Reading should only be done within the indicated interval; late or subjective interpretations may produce false results.

INTERNAL QUALITY CONTROL

The clinical laboratory should implement an internal quality control program with clearly established procedures, standards, limits, and tolerances for variations. It is essential to recognize that all measurement systems have inherent analytical variability, which should be monitored by the laboratories themselves. For this purpose, the use of controls is recommended to assess precision and accuracy.

PRODUCT PERFORMANCE

ACCURACY

CLINICAL SENSITIVITY AND SPECIFICITY

This product was tested in comparison with other methods, being Real-Time PCR and another commercial kit, the data can be seen below:

A) Comparison with Real-Time PCR

A total of 115 samples characterized by Real-Time PCR were analyzed based on the different genes for blood type A, type B and type AB, being 90 samples of type A, 14 samples of type B and 11 samples of type AB.

The results show that the clinical sensitivity of the product is > 99% and the clinical specificity is > 99%.

METHOD	REFERENCE - Real-Time PCR				Total
	Result	Type A	Type B	Type AB	
Feline Blood Typing VET FAST	Type A	90	0	0	90
	Type B	0	14	0	14
	Type AB	0	0	11	11
	Result Total	90	14	11	115

Clinical Sensitivity: > 99% - 95% CI 95% - 100%

Clinical Specificity: > 99% - 95% CI 95% - 100%

Accuracy: 100%