TP BIOCLIN

 **K089**

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FINALIDADE

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- 1 - Hirsh J , Deykin D, Poller L: Therapeutic Range For Oral Anticoagulant Therapy. Chest 89/12 (Suppl) 115-155, 1986.
- 2 - International Committee For Standardization In Haematology And International Committee On Trombosis And Haemostasis. Amer. J. Clin. Path . 88:779 1985.
- 3 - National Committee For Clinical Laboratory Standards: Collecion, Transport And Processing Of Blood Specimens For Coagulation Testing And Performance Of Coagulation Assays, 1991. Nccls.
- 4 - Douketis Jd, Lane A , Milne J , Ginsberg Js. Trombosis Research , 1998 ; 92 : 11-17.
- 5 - QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

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- 1 - Hirsh J , Deykin D, Poller L: Therapeutic Range For Oral Anticoagulant Therapy. Chest 89/12 (Suppl) 115-155, 1986.
- 2 - International Committee For Standardization In Haematology And International Committee On Trombosis And Haemostasis. Amer. J. Clin. Path . 88:779 1985.
- 3 - National Committee For Clinical Laboratory Standards: Coleccion, Transport And Processing Of Blood Specimens For Coagulation Testing And Performance Of Coagulation Assays, 1991. Nccls.
- 4 - Douketis Jd, Lane A, Milne J, Ginsberg Js. Trombosis Research, 1998; 92: 11-17.
- 5 - QUIBASA: Datos do Departamento de Pesquisa e Desenvolvimento.

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INSTRUCTIONS FOR USE

FUNCTION

Method for determination of manual or automated of Prothrombin Time in citrated plasma. For *in vitro* diagnostic use only.

PRINCIPLE OF ACTION

**Methodology:** Coagulometry  
Consists of measuring the Clotting Time of plasma after addition a source of Thromboplastin (factor III, tissue factor) and Calcium. The factor VII is then activated to form a complex (complex FT-FVIIa) that will activate factors IX and factor Xa XO with Phospholipids factor tissue, Calcium and factor Va form the Prothrombin activator complex which transforms Prothrombin (factor II) into Thrombin (factor IIa). Thrombin acts on Fibrinogen, transforming it into Fibrin. This test assesses the factors of Prothrombin Complex (factors II, V, VII and X).

REAGENT

**Number 1 - Clot Former** - Store between 2 and 8°C. Contains: Rabbit brain extract 2.6%, Calcium Chloride, buffer and preservative.  
**Presentation ready for use.**

PRESENTATION

| Presentation | Reagent 1 |
|--------------|-----------|
| 1            | 10 x 5 mL |
| 2            | 20 x 2 mL |
| 3            | 10 x 2 mL |
| 4            | 20 x 5 mL |

EQUIPMENTS AND OPERATIONAL INPUTS

Log-log paper, pipette and tips, tubes, coagulometer or stopwatch, Biocontrol N and P Coagulation Bioclin. These itens are found at markets specialized on Laboratories of Clinical Analysis.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature is between 2 and 8°C and the transport temperature shall be up to 30°C, not exceeding 5 days. Protect from light and avoid moisture. **Do not freeze.**

SPECIAL CARE

- 1- For professional *in vitro* diagnostic use only.
- 2- Follow strictly the proposed methodology to obtain results accurate.
- 3- Handle with care reagents because they contain Sodium Azide, irritating to skin and mucous membranes.
- 4- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.
- 5- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Safety Data Sheet Safety of Chemicals) available on the website [www.bioclin.com.br](http://www.bioclin.com.br) or upon request by the SAC (Advisory Service Customer) of Quibasa.
- 6- Do not use the product in case of damaged packaging.
- 7- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

**Material:** Plasma  
**Anticoagulant:** Use the kit as anticoagulant Citrate Bioclin A001 as written in usage instructions of the same, or 3.8% Sodium Citrate (0.130 M) at a ratio of 9 parts blood to 1 part anticoagulant (Example: 4.5 mL blood + 0.5 mL of anticoagulant).  
**Collection:** Blood should be obtained by venipuncture, avoiding hemolysis, prolonged tourniquet, blistering and aspiration tissue fluid.  
**Preparation:** Blood should be mixed with the anticoagulant soon after collection. Centrifuge immediately at 3000 rpm for 15 minutes. Remove the plasma without cellular pipette red or yellow layer. The samples should be tested in less than three hours. If the test can not be done in this period the plasma should be frozen for at maximum 1 week at -20°C.

PROCESS DESCRIPTION

The calibration stability of the PT Bioclin kit installed on refrigerated equipment is at least 14 days. This stability may vary depending on the conditions of the test, equipment and environment. Therefore, it is suggested to follow the product performance using control serum.

TECHNIQUE

Bioclin recommends, as control serum, Biocontrol Coagulation N and P Bioclin kits.

- 1- Homogenize Reagent N°1 before use.
- 2- Preheat Reagent N°1 to 37°C for 4 to 5 minutes. This procedure should be performed only with the volume that will be used to perform the tests. The reagent must remain warm for the entire procedure.
- 3- Separate the reaction tubes that will be used to perform the tests.
- 4- Pipette 100 µL of the Sample/Control into the respective reaction tube and incubate at 37°C for 2 to 3 minutes in a water bath or thermal block.
- 5- Pipette 200 µL of Reagent N°1 into the Sample/Control tube (preheated to 37°C) in the Sample/Control tube, simultaneously starting the timer.
- 6- Leave the tube in the water bath at 37°C, shaking gently. Keep in the bath for about 09 seconds.
- 7- Remove the tube from the bath, gently tilt once per second and stop the timer at the time of the clot formation. Note the clotting time.

**Note:** After opened, reagent N° 1 is stable for 14 days at 2 to 8°C.

CALCULATIONS AND RESULTS

The results of Prothrombin Time (PT) may be expressed as:  
**Quick Time:** In seconds.  
**Percentage of activity:** Related to a pool of normal plasmas (considered 100% activity), and for this you must draw a curve activity, using a normal plasma pool and freshly harvested and diluted with saline (0.9%).

Curve

| Dilution                  | 1:1  | 1:2 | 1:3 | 1:4 | 1:8   | 1:10 |
|---------------------------|------|-----|-----|-----|-------|------|
| Pool (mL)                 | 0.2  | 0.2 | 0.2 | 0.2 | 0.2   | 0.2  |
| Physiologic Solution (mL) | ---  | 0.2 | 0.4 | 0.6 | 1.4   | 1.8  |
| Activity                  | 100% | 50% | 33% | 25% | 12.5% | 10%  |

In a log-log paper, plot the graph by putting the seconds value on the “y” axis and the percentage value on the “x” axis. This procedure should be repeated on every batch.

INR (International Normalized Ratio)

The INR is calculated using the ratio between the patient's PT and the PT of the normal pool, raised to a power, called the International Sensitivity Index (ISI). The ISI value is stated in each kit.

INR = (patient PT / normal pool PT)<sup>ISI</sup>

Example:  
Patient PT: 14.5 s  
Normal pool PT: 14.0 s  
ISI: 1.17  
INR = (14.5 / 14.0)<sup>1.17</sup> = 1.04

INTERFERENCES

Excessive hemolysis and lipemia. Use of plasma with Heparin and EDTA. The PT can be increased by substances such as Anticoagulants Heparin and Warfarin, Oral Contraceptives, Ethanol, Tetracycline, Corticosteroids and EDTA. The reduction of PT is observed in patients taking: Vitamin K, Anti-Histamines, Caffeine and Phenobarbital.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

TRACEABILITY

The PT Bioclin kit is traceable to the reference material NIBSC RBT/05 (International Standard for Thromboplastin - WHO International Standard).

REFERENCE VALUES

Normal time values vary between laboratories, depending on the reagent, equipment and validated methodology. Therefore, it is recommended that each laboratory establish its reference values, according to its population and the technique used.  
- **Indicative range of prothrombin time (quick)** - 10 to 14 seconds. It should be noted that the reference value varies depending on the measurement system.  
- **Activity (%)** - 70 to 100%. Values greater than 100% have no clinical significance.  
- **INR** values for monitoring oral anticoagulants vary according to pathology and therapeutic indication and each laboratory must establish its own reference values.  
The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

PRODUCT PERFORMANCE

ACCURACY

**Comparison of Methods and Specificity Methodology**  
PT Bioclin kit was compared with another method for evaluating Prothrombin Time commercially available. Were performed 07 analysis and the results were evaluated. The obtained linear equation was Y = 1.0804X - 1.2695 and correlation coefficient 0.9954. With these results, we can conclude that the kit shows good methodological specificity.

SPECIFICITY

The specificity of the kit is provided by the addition of an exclusive source thromboplastin, specific activator of factor VII, which eliminates the possibility of direct activation of the cascade of other factors coagulation.

PRECISION

Repeatability

The repeatability was calculated from 20 successive determinations, using 3 samples with different concentrations, obtaining the following results:

|                              | Sample 1 | Sample 2 | Sample 3 |
|------------------------------|----------|----------|----------|
| Time (s)                     | 12.65    | 24.85    | 17.45    |
| Standard Deviation (%)       | 0.49     | 0.75     | 0.51     |
| Coefficient of Variation (%) | 3.87     | 3.00     | 2.93     |

Reproducibility

The reproducibility was calculated from 20 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

|                              | Sample 1 | Sample 2 | Sample 3 |
|------------------------------|----------|----------|----------|
| Time (s)                     | 12.57    | 24.77    | 17.40    |
| Standard Deviation (%)       | 0.19     | 0.08     | 0.05     |
| Coefficient of Variation (%) | 1.51     | 0.31     | 0.29     |

SENSITIVITY

The analytical sensitivity of the kit is given by the International Sensitivity Index (ISI) value.

DIAGNOSTIC SIGNIFICANCE

The test evaluates changes related to congenital and acquired deficiencies of factors of the extrinsic and common coagulation pathways, and may be altered in isolated congenital deficiencies of factors VII, X, V, II or I, diseases that affect the liver, vitamin deficiency K, disseminated intravascular coagulation, hemorrhagic diseases and use of oral anticoagulants.

BIBLIOGRAPHIC REFERENCES

- 1 - Hirsh J, Deykin D, Poller L: Therapeutic Range For Oral Anticoagulant Therapy. 89/12 Chest (Suppl) 115-155, 1986.
- 2 - International Committee For Standardization In Haematology International Committee On Thrombosis And Haemostasis. Amer. J. Clin. Path. 88:779 1985.
- 3 - National Committee For Clinical Laboratory Standards: Coleccion, Transport And Processing Of Blood Coagulation Specimens For Testing And Performance Of Coagulation Assays, 1991. NCCLS.
- 4 - Douketis JD, Lane A, Milne J, Ginsberg JS. Thrombosis Research, 1998, 92: 11-17.
- 5 - Quibasa: Data from the Department of Research and Development.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca  
CEP 31565-130 - Belo Horizonte - MG - Brasil  
Phone: +55 31 3439.5454 | E-mail: bioclin@bioclin.com.br  
CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service  
Phone: 0800 0315454 | E-mail: sac@bioclin.com.br

ANVISA resgistration for PT Bioclin kit: 10269360151

Review: December/2023

UNIVERSAL SYMBOLOGY

CATALOG NUMBER

LOT NUMBER

MANUFACTURING DATE

VALIDITY DATE  
(last day of the month)

TEMPERATURE LIMIT  
(store)

CONTENT IS SUFFICIENT  
FOR <N> TEST

SEE INSTRUCTIONS  
FOR USE

IN VITRO DIAGNOSTIC  
PRODUCT

KEEP AWAY  
FROM SUNLIGHT

DO NOT REUSE

CAUTION

MADE BY

CONTROL

POSITIVE CONTROL

NEGATIVE CONTROL

BIOLOGICAL RISK

FLAMMABLE

CORROSIVE

TOXIC

DO NOT USE IF  
PACKAGE IS  
DAMAGED

PRODUCT  
STERILIZED

DANGER