Bioclin

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REF

K207

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- 1 - Hirsh J , Deykin D, Poller L: Therapeutic Range For Oral Anticoagulant Therapy. Chest 89/12 (Suppl) 115-155, 1986.
- 2 - International Committee For Standardization In Haematology And International Committee On Trombosis And Haemostasis. Amer. J. Clin. Path . 88:779 1985.
- 3 - National Committee For Clinical Laboratory Standards: Colleccion, Transport And Processing Of Blood Specimens For Coagulation Testing And Performance Of Coagulation Assays, 1991. Nccls.
- 4 - Douketis Jd, Lane A , Milne J , Ginsberg Js. Trombosis Research , 1998 ; 92 : 11-17.
- 5 - QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

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REFERÊNCIAS BIBLIOGRÁFICAS
1 - Hirsh J , Deykin D, Poller L: Therapeutic Range For Oral Anticoagulant Therapy. Chest 89/12 (Suppl) 115-155, 1986.
2 - International Committee For Standardization In Haematology And International Committee On Trombosis And Haemostasis. Amer. J. Clin. Path . 88:779 1985.
3 - National Committee For Clinical Laboratory Standards: Colleccion, Transport And Processing Of Blood Specimens For Coagulation Testing And Performance Of Coagulation Assays, 1991. Nccls.
4 - Douketis Jd, Lane A , Milne J , Ginsberg Js. Trombosis Research , 1998 ; 92 : 11-17.
5 - QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

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INSTRUCTIONS FOR USE

FUNCTION

Test to determine the Prothrombin Time (PT) in biological samples (citrated plasma), using the coagulometric test. Only for *in vitro* diagnostic use.

PRINCIPLE OF ACTION

Methodology: Coagulometry

It consists of measuring the plasma clotting time after the addition of a source of thromboplastin (factor III, tissue factor), phospholipids and calcium. Factor VII is then activated, forming the FT-FVIIa complex that will activate factors IX and X. Factor Xa together with factor Va and calcium transform prothrombin (factor II) into thrombin (factor IIa). Thrombin acts on fibrinogen, converting it into fibrin. This test evaluates the factors of the extrinsic and common coagulation pathways (factors II, V, VII, IX and X).

REAGENT

Number 1 - Clot Former - Store between 2 and 8°C. Contains: Thromboplastin < 5% obtained from rabbit brain extract, buffer solution, activator and preservative.
Ready-to-use presentation.

PRESENTATION

Presentation	Reagent 1
1	10 x 5 mL
2	20 x 2 mL
3	10 x 2 mL
4	20 x 5 mL
5	6 x 5 mL
6	12 x 10 mL
7	12 x 15 mL
8	12 x 20 mL

OPERATIONAL EQUIPMENT AND CONSUMPTIONS

Log-log paper, tips and pipettes, tubes, coagulometer or stopwatch. They are found in the specialized market for articles for Clinical Analysis Laboratories.

STORAGE AND TRANSPORT CONDITIONS

The storage temperature should be 2 to 8°C. Transport, at temperatures up to 30°C, should not exceed 5 days. Keep away from light and avoid moisture. **Do not freeze.**

SPECIAL CARES

- 1- For *in vitro* diagnostic use only.**
- Strictly follow the proposed methodology to obtain exact results.
- Handle the reagents with care as they contain Sodium Azide, which irritates the skin and mucous membranes.
- We recommend applying local, state and federal environmental protection standards so that reagents and biological material are disposed of in accordance with current legislation.
- To obtain information related to biosafety or in the event of accidents with the product, consult the SDS (Safety Data Sheet) available on the website www.bioclin.com.br or through a request through the SAC (Service Customer Advisory) from Quibasa.
- Do not use the product if the packaging is damaged.

7- It is essential that the instruments and equipment used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Material: Plasma

Anticoagulant: Use 3.8% Sodium Citrate (0.130 M) as an anticoagulant, in a proportion of 9 parts of blood to 1 part of anticoagulant (for example: 4.5 mL of blood + 0.5 mL of anticoagulant).

Collection: Blood must be obtained by venipuncture, avoiding hemolysis, prolonged tourniquet, blister formation and aspiration of tissue fluid.

Preparation: The blood must be mixed with the anticoagulant immediately after collection. Centrifuge immediately at 3000 rpm for 15 minutes. Remove plasma without pipetting red cells or the yellow layer. Samples must be tested in less than 3 hours. If the test cannot be performed within this period, the plasma must be frozen for a maximum of 1 week at -20°C.

Interferers:

Excessive hemolysis and lipemia. Use of plasmas with Heparin and EDTA. PT can be increased by substances such as: Anti-Coagulants Heparin and Warfarin, Oral Contraceptives, Ethanol, Tetracycline, Corticosteroids and EDTA. The reduction in PT is observed in patients using: Vitamin K, Antihistamines, Caffeine and Phenobarbital.

PROCESS DESCRIPTION

The stability of the PT ULTRASENSITIVE kit installed in refrigerated equipment is up to 14 days. This stability may vary according to test, equipment and environmental conditions. Therefore, it is suggested to monitor the product's performance using control plasmas.

TECHNIQUE

Bioclin recommends, to use the kit, use the Biocontrol Coagulation N and P Bioclin's kits as control plasma.

- 1- Preheat reagent No. 1 (Clot Former) to 37°C for 4 to 5 minutes (heat only as much as necessary to carry out the test).
- 2- Homogenize the reagent before use.
- 3- In a clean tube, pipette 100 µL of the plasma to be measured (controls or test samples) and incubate at 37°C for 2 to 3 minutes (in a water bath or thermal blocks).
- 4- Add 200 µL of reagent No. 1 (Clot Former), previously heated, to the tube containing the plasma and start the timer.
- 5- Leave the tube in the heat, shaking it slowly.
- 6- Remove the tube before the expected clotting time (approximately 9 seconds), observe the clot formation. Immediately stop the timer when you observe clot formation and record the time.
- 7- Calculate the average of the results obtained by the duplicate.

Note: Once opened, reagent No. 1 is stable for 14 days at 2 to 8°C.

CALCULATIONS AND RESULTS

Prothrombin Time (PT) results can be expressed as:

Quick time: In seconds.

Percentage of activity: Related to a pool of normal plasma (considered 100% activity), and for this an activity curve must be drawn, using a pool of normal plasma recently collected and diluted with physiological saline (0.9%).

Curve

Dilutions	1:1	1:2	1:3	1:4	1:8	1:10
Pool (mL)	0.2	0.2	0.2	0.2	0.2	0.2
Physiological solution (mL)	---	0.2	0.4	0.6	1.4	1.8
Activity	100%	50%	33%	25%	12.5%	10%

Determine the Prothrombin Time for each dilution in duplicate, calculating the average value of the results.
On log-log paper, plot the graph by placing the seconds value on the "y" axis and the percentage value on the "x" axis. This procedure must be repeated for each batch.

INR (International Normalized Ratio)

The INR is calculated using the ratio between the patient's PT and the PT of the normal pool, raised to a power, called the International Sensitivity Index (ISI). The ISI value is stated in each kit.

$$INR = (\text{patient PT} / \text{normal PT})^{ISI}$$

Example:
Patient PT: 14.5 s
Normal PT (reference): 14 s
ISI: 1.00

$$INR = (14.5 / 14)^{1.00} = 1.0$$

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control program, where procedures, standards, limits and tolerance for variations are clearly established. It is important to highlight that all measurement systems present a characteristic analytical variability, which must be monitored by the laboratories themselves. To this end, it is recommended to use controls, which allow evaluating the precision and dosage accuracy.

TRACEABILITY

The PT ULTRASENSITIVE kit is traceable to the NIBSC RBT/05 reference material (International Standard for Thromboplastin - WHO International Standard).

REFERENCE VALUES

Normal time values vary between laboratories, depending on the reagent, equipment and validated methodology. Therefore, it is recommended that each laboratory establish its reference values, according to its population and the technique used.

- **Indicative range of prothrombin time (quick)** - 12 to 15 seconds. It should be noted that the reference value varies depending on the measurement system.
 - **Activity (%)** - 70 to 100%. Values greater than 100% have no clinical significance.
 - **INR** values for monitoring oral anticoagulants vary according to pathology and therapeutic indication and each laboratory must establish its own reference values.
- The results provided by this kit must be interpreted by the responsible medical professional, and are not the only criteria for determining the patient's diagnosis and/or treatment.

PRODUCT PERFORMANCE

ACCURACY

Comparison of Methods and Methodological Specificity

The PT ULTRASENSITIVE kit was compared with another commercially available method for evaluating Prothombine Time. 20 analyzes were performed and the results were evaluated. The linear equation obtained was $Y = 0.9591X + 0.5508$ and the correlation coefficient was 0.9811. With these results, it can be concluded that the kit presents good methodological specificity.

Interferences (Analytical Specificity)

The specificity of the kit is guaranteed by the addition of an exclusive source of thromboplastin, a specific activator of factor VII, which eliminates the possibility of direct activation of other factors in the coagulation cascade.

According to the scientific literature, some substances are known as possible interferers for the coagulometry methodology for measuring thromboplastin time, which are causes of falsely elevated or decreased results.

Hemolysis and lipemia are interfering factors that, in excess, can alter the results in an unpredictable way.

Furthermore, PT can be increased by substances such as: anticoagulants (Heparin and EDTA), warfarin, oral contraceptives, ethanol, tetracycline and corticosteroids.

The reduction in PT can be observed in patients using: Vitamin K, antihistamines, caffeine and phenobarbital.

PRECISION

Repeatability

The repeatability was calculated from 10 successive determinations, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Time (Sec)	12.50	38.35	14.59
Standard Deviation	0.28	0.30	0.27
Coefficient of variation (%)	2.26	0.77	1.87

Reproducibility

Reproducibility was calculated from 10 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2	Sample 3
Time (Sec)	12.48	38.33	14.50
Standard Deviation	0.02	0.14	0.09
Coefficient of variation (%)	0.17	0.36	0.59

SENSITIVITY

The analytical sensitivity of the kit is given by the International Sensitivity Index (ISI) value.

DIAGNOSTIC MEANING

The test evaluates changes related to congenital and acquired deficiencies of factors of the extrinsic and common coagulation pathways, and may be altered in isolated congenital deficiencies of factors VII, X, V, II or I, diseases that affect the liver, vitamin deficiency K, disseminated intravascular coagulation, hemorrhagic diseases and use of oral anticoagulants.

BIBLIOGRAPHIC REFERENCES

- 1 - Hirsh J , Deykin D, Poller L: Therapeutic Range For Oral Anticoagulant Therapy. Chest 89/12 (Suppl) 115-155, 1986.
- 2 - International Committee For Standardization In Haematology And International Committee On Trombosis And Haemostasis. Amer. J. Clin. Path . 88:779 1985.
- 3 - National Committee For Clinical Laboratory Standards: Collecion, Transport And Processing Of Blood Specimens For Coagulation Testing And Performance Of Coagulation Assays, 1991. Nccls.
- 4 - Douketis Jd, Lane A , Milne J , Ginsberg Js. Trombosis Research , 1998 ; 92 : 11-17.
- 5 - QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Quality Control Department. The quality of the reagents is guaranteed until the expiration date mentioned on the presentation packaging, as long as they are stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda
Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Tel.: (31) 3439.5454 | E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service
Tel.: 0800 0315454 | E-mail: sac@bioclin.com.br

Registration number of the PT ULTRASENSITIVE kit at ANVISA: 10269360313

Revision: April/2024

UNIVERSAL SYMBOLOGY

	CATALOG NUMBER		MADE BY
	LOT NUMBER		CONTROL
	MANUFACTURING DATE		POSITIVE CONTROL
	VALIDITY DATE (last day of the month)		NEGATIVE CONTROL
	TEMPERATURE LIMIT (store)		BIOLOGICAL RISK
	CONTENT IS SUFFICIENT FOR <N> TEST		FLAMMABLE
	SEE INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC PRODUCT		TOXIC
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED
	DO NOT REUSE		PRODUCT STERILIZED
	CAUTION		DANGER