# Bioclin

## TRANSAMINASE AST (TGO) CINÉTICA

REF K048

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#### FINALIDADE

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- 1 - The Committee on Enzymes of the Scandinavian Society for Clinical Chemistry and Clinical Physiology Scand J. Clin Lab Invest., 1974, 33, 291 - 306.
- 2 - BERGMAYER, Bowers and cols., Clin. Chim., Acta., 1.976, 70, 19 - 42, 1977, 21 - 22.
- 3 - BERGMAYER, HV.; SCHEIBE, P.; WAHLEFELD, A.W., Clin., Chem., 1.978, 24, 58 - 73.
- 4 - Expert Panel on Enzymes of the International Federation of Clinical Chemistry: Part 3. Revised IFCC Method for Aspartate Aminotransferase, Clin.,Chem., 1.978, 24, 720.
- 5 - Scandinavian Committee on Enzymes Scand J. Clin Lab Invest, 1.981, 41, 107 - 116.
- 6 - BURTIS; CARL, A.; ASHWOOD; EDWARD, R., Clin. Chem., Tietz TEXT BOOK of; 2a ed., 1.986, 788 - 797.
- 7 - PESCE, A., J.; KAPLAN, L., A., Methods in Clin. Chem., 1.987.
- 8 - SOLDIN, S.J., BRUGNARA, C., WONG, E.C.: Pediatric Reference Intervals, 5.ed. Washington: AACC Press, 2005.p.3-4.
- 9 - WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31
- 10 - ELLIS S. Benson M.D. & D. B. Tonks PH.D. (1963) Laboratory medicine, Postgraduate Medicine, 34:4, A-58-A-70, DOI: 10.1080/00325481.1963.11694890
- 11 - QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

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- 1 - The Committee on Enzymes of the Scandinavian Society for Clinical Chemistry and Clinical Physiology Scand J. Clin Lab Invest., 1974, 33, 291 - 306.
- 2 - BERGMAYER, Bowers and cols., Clin. Chim., Acta., 1.976, 70, 19 - 42, 1977, 21 - 22.
- 3 - BERGMAYER, HV.; SCHEIBE, P.; WAHLEFELD, A.W., Ciin., Chem., 1.978, 24, 58 - 73.
- 4 - Expert Panel on Enzymes of the International Federation of Clinical Chemistry: Part 3. Revised IFCC Method for Aspartate Aminotransferase, Clin.,Chem., 1.978, 24, 720.
- 5 - Scandinavian Committee on Enzymes Scand J. Clin Lab Invest, 1.981, 41, 107 - 116.
- 6 - BURTIS; CARL, A.; ASHWOOD; EDWARD, R., Clin. Chem., Tietz TEXT BOOK of; 2a ed., 1.986, 788 - 797.
- 7 - PESCE, A., J.; KAPLAN, L., A., Methods in Clin. Chem., 1.987.
- 8 - SOLDIN, S.J., BRUGNARA, C., WONG, E.C.: Pediatric Reference Intervals, 5.ed. Washington: AACC Press, 2005.p.3-4.
- 9 - WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31
- 10 - ELLIS S. Benson M.D. & D. B. Tonks PH.D. (1963) Laboratory medicine, Postgraduate Medicine, 34:4, A-58-A-70, DOI: 10.1080/00325481.1963.11694890
- 11 - QUIBASA: Datos del Departamento de Investigación y Desarrollo.

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### INSTRUCTIONS FOR USE

#### FUNCTION

Test for determination of Aspartate Amino Transferase (AST or TGO) in biological sample (serum or plasma), through kinetic test. Only for *in vitro* diagnostic use.

#### PRINCIPLE OF ACTION

**Methodology:** Kinetic (UV)

AST catalyzes the transfer of Amine groups from Aspartate to Alfa - Ketoglutarate, leading to the formation of Glutamate and Oxalacetate. Oxalacetate in the presence of MDH reacts with NADH, reducing to Malate and NADH oxidizes to NAD<sup>+</sup>. The oxidation rate is proportional to the activity of AST in the sample.

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#### REAGENTS

**Number 1 - Substrate** - Store between 2 e 8 °C. Contain Buffer < 200 mmol/L, Sodium Aspartate < 450 mmol/L, D-Lactate Dehydrogenase < 5 KU, Malate Dehydrogenase < 5 KU, chelator, surfactant and preservative.

**Number 2 - Coenzyme** - Store between 2 e 8 °C. Contém: Tris Buffer < 300 mmol/L, Alpha Ketoglutaric Acid < 180 mmol/L, NADH < 5 mmol/L, surfactant and preservative.

#### PRESENTATION

Presentation	Reagent N° 1	Reagent N° 2
1	2 x 40 mL	2 x 10 mL
2	4 x 40 mL	4 x10 mL

#### EQUIPMENTS AND OPERATIONAL INPUTS

Thermostatic spectrophotometer, pipettes and tips, test tubes, clock or stopwatch, Biocal, Biocontrol N e Biocontrol P Bioclin. These items are found in the specialized market for articles for Clinical Analysis Laboratories.

#### TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be between 2 to 8 °C. The transport, at temperatures up to 30 °C, should not exceed 5 days. Protect from light and avoid moisture. **Do not freeze.**

#### SPECIAL CARE

**1- For *in vitro* diagnostic use only.**

**2-** Strictly follow the methodology proposed to obtain exact results.

**3-** Water used in material cleaning must be recent and free of contaminants.

**4-** Saturated deionizing columns release alkaline water, various ions and oxidizing and reducing agents, which can significantly alter the results.

**5-** Visible hemolysis can cause slight variations in results.

**6-** We recommend applying local, state and federal environmental protection standards so that the disposal of reagents and biological material is carried out in accordance with current legislation.

**7-** To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Safety Information Sheet for Chemical Products) available on the website [www.bioclin.com.br](http://www.bioclin.com.br) or upon request by the Customer Service Customer) of Quibasa

**8-** Do not use the product in case of damage to the packaging.

**9-** It is essential that the instruments and equipment used are properly calibrated and subjected to periodic maintenance.

**10-** Reagent 2 contains Sodium Azide and must be handled with care.

**11-** It is important, for the good performance of the test, a strict control of time, temperature and pH.

**12-** Lipemic and icteric samples increase the absorbance by 340 nm. In this case, dilute the sample 1: 2 with physiological saline multiply the result obtained by 2.

**13-** Samples stored without centrifugation for long periods may have their pyruvate concentration increased, which may alter the results. See more information in Interferents.

**14-** Absorbance readings below 0.800 of the Working Reagent indicate loss of the same. In this case, do not use the reagent.

#### SAMPLES

**Serum or Plasma:** collected with EDTA or Heparin, obtained free from hemolysis. The enzyme is stable for 3 days at 2 to 8 °C and 3 months at -20 °C.

#### PROCESS DESCRIPTION

##### PREPARATION OF WORKING REAGENT

Mix 4 parts of Reagent N° 1 with 1 part of Reagent N° 2. Working Reagent is stable during 72 hours between 15 and 30 °C and 14 days between 2 and 8 °C.

The stability of the Transaminase AST (TGO) product Kinetics after opening (In use / On board - Product installed in equipment) is 30 days. This stability can vary depending on the condition of the problem. Therefore, it is suggested to monitor the performance of the product using control numbers.

#### REACTION CONDITIONS

It is indispensable to use a thermostated cuvette at 37 °C, 1 cm optical path and reading at 340 nm (334 - 365 nm).

#### TECHNIQUE

Bioclin recommends using the Biocal kit as a calibrator and the Biocontrol N and P Bioclin kits as the control serum.

Add 100 µL of Sample to 1.0 mL of the Working Reagent, mix and transfer to a thermostated cell at 37 °C and wait 1 minute. Take the initial reading, simultaneously triggering the stopwatch. Repeat the readings after 1, 2 and 3 minutes. Calculate the average of the absorbance differences per minute ( $\Delta\text{Abs.} / \text{Min.}$ ) And use it to calculate the result.

#### DESCRIPTION OF THE CALCULATIONS

AST - Results expressed in U/L

340 nm =  $\Delta\text{Abs.}/\text{min.} \times 1746$

334 nm =  $\Delta\text{Abs.}/\text{min.} \times 1780$

365 nm =  $\Delta\text{Abs}/\text{min.} \times 3235$

#### PROCESS LIMITATIONS

The specifications below refer to semi-automatic equipment:

The kinetic method is based on molar absorptivity, for this reason, the readings must be performed on a spectrophotometer that fulfills the following conditions:

- Wavelength 340 nm
- Semi-trajectory of the 10 nm passband
- Spurious light less than 0.5%
- Thermostatized 1cm bucket

#### INTERFERENT

No interference was observed by Bilirubin until 19 mg/dL, Hemoglobin until 45 mg/dL, Triglycerides until 650 mg/dL and Pyruvate until 0.2 mmol/L.

#### INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

#### TRACEABILITY

The kit calibration can be made using the theoretical calibration factor based on the molar absorptivity of the NADH, or through the BIOCAL calibrator. Bioclin recommends the usage of the BIOCAL calibrator which is traceable to the reference material ERM-AD457.

#### REFERENCE VALUES

The reference values in U/L for the present method were obtained by determining AST in healthy male and female populations.

Age	Male (U/L)	Female (U/L)
1 - 7 days	26 - 98	20 - 93
8 - 30 days	16 - 67	20 - 69
1 - 6 months	16 - 62	16 - 61
7 - 12 months	16 - 52	16 - 60
1 - 3 years	16 - 57	16 - 57
4 - 6 years	10 - 47	10 - 47
7 - 15 years	10 - 41	5 - 36
Adults	< 38	< 31

These values should be used as a guide, and each laboratory should create its range of reference values, according to the population served. The results provided by this kit must be interpreted by the responsible medical professional, not being the only criterion for determining diagnosis and / or treatment of the patient.

#### PRODUCT PERFORMANCE

##### ACCURACY

##### Comparison of Methods

The Transaminase AST (TGO) Kinetic kit was compared with another commercially available method for measuring aspartate aminotransferase. 40 analyzes were carried out and the results were evaluated. The linear equation obtained was  $Y = 0.9765X + 1.0283$ , with a linear correlation coefficient equal to 0.9999. With these results, it can be concluded that the kit has good methodological specificity.

#### PRECISION

##### Repeatability

The repeatability was calculated from 10 successive determinations, using 3 samples with different concentrations, obtaining the following results:

Repeatability	Sample 1	Sample 2	Sample 3
Average Concentration(U/L)	35.28	41.92	58.86
Standard Deviation (U/L)	0.52	0.06	0.11
Coefficient of Variation (%)	1.46	0.13	0.19

##### Reproducibility

Reproducibility or was calculated from the analysis of 3 samples of different levels during 15 days, not necessarily consecutive, totaling 60 determinations for each sample, obtaining the following results:

Reproducibility	Sample 1	Sample 2	Sample 3
Average Concentration(U/L)	40.34	45.71	58.53
Standard Deviation (U/L)	0.4	0.5	0.6
Coefficient of Variation (%)	1.01	1.01	1.02

### SENSITIVITY

Sensitivity was calculated from 40 determinations of a sample free from Aspartate aminotransferase. The average found was 2.10 U/L with a standard deviation of 0.585 U/L. The sensitivity, which indicates the detection limit of the method, corresponds to 3 times the standard deviation, and is equal to 1.756.

### LINEARITY

The reaction is linear up to a concentration of 400 U/L. For samples with values above 400 U/L, it is recommended to dilute the sample with 0.85% Sodium Chloride, repeat the dosage and multiply the result obtained by the dilution factor. For example, if the sample is diluted 1:10 and the result obtained is 80 U/L, multiply by 10 (80 x 10) U/L and consider the activity value of 800 U/L as the final result.

### DIAGNOSTIC SIGNIFICANCE

The increase in the activity of the enzymes Aspartate AminoTransferase - AST (located in the cytomitochondrial region) reflects changes in several tissues. This enzyme is found in high concentration in the heart, liver, skeletal muscle, kidneys and pancreas. Its plasma activity increases after 6 to 8 hours after myocardial infarction, reaching a peak in 24 to 48 hours after the onset. Considerable increases occur in viral, toxic hepatitis, liver necrotic diseases - 3 to 50 or 100 times the reference values (VR), mononucleoses (20 times the VR), intrahepatic cholestasis (20 times the VR) and muscular dystrophies (8 times the VR). In chronic liver diseases associated with cell necrosis, due to the increase in AST, inversion of the ALT (TGP) / AST ratio may occur.

### BIBLIOGRAPHIC REFERENCES

- 1 - The Committee on Enzymes of the Scandinavian Society for Clinical Chemistry and Clinical Physiology Scand J. Clin Lab Invest., 1974, 33, 291 - 306.
- 2 - BERGMAYER, Bowers and cols., Clin. Chim., Acta., 1.976, 70, 19 - 42, 1977, 21 - 22.
- 3 - BERGMAYER, HV.; SCHEIBE, P.; WAHLEFELD, A.W., Clin., Chem., 1.978, 24, 58 - 73.
- 4 - Expert Panel on Enzymes of the International Federation of Clinical Chemistry: Part 3. Revised IFCC Method for Aspartate Aminotransferase, Clin.,Chem., 1.978, 24, 720.
- 5 - Scandinavian Committee on Enzymes Scand J. Clin Lab Invest, 1.981, 41, 107 - 116.
- 6 - BURTIS; CARL, A.; ASHWOOD; EDWARD, R., Clin. Chem., Tietz TEXT BOOK of; 2a ed., 1.986, 788 - 797.
- 7 - PESCE, A., J.; KAPLAN, L., A., Methods in Clin. Chem., 1.987.
- 8 - SOLDIN, S.J., BRUGNARA, C., WONG, E.C.: Pediatric Reference Intervals, 5.ed. Washington: AACC Press, 2005.p.3-4.
- 9 - WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2,2002:31
- 10 - ELLIS S. Benson M.D. & D. B. Tonks PH.D. (1963) Laboratory medicine, Postgraduate Medicine, 34:4, A-58-A-70, DOI: 10.1080/00325481.1963.11694890
- 11- QUIBASA: Data from the Research and Development Department.

### QUALITY ASSURANCE

Before being released for consumption, all Bioclin reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

**QUIBASA QUÍMICA BÁSICA Ltda**

Rua Teles de Menezes, 92 - Santa Branca

CEP 31565-130 - Belo Horizonte - MG - Brasil

Phone.: +55 31 3439.5454 | E-mail: bioclin@bioclin.com.br

CNPJ: 19.400.787/0001-07 - Made in Brazil

### CUSTOMER SERVICE

Customer Advisory Service

Phone: 0800 0315454 | E-mail: customerservice@bioclin.com.br

ANVISA registration for Transaminase AST Kinetic kit: 10269360119

Review: July/2023

## UNIVERSAL SYMBOLOGY

CATALOG NUMBER

MADE BY

LOT NUMBER

CONTROL

MANUFACTURING DATE

POSITIVE CONTROL

VALIDITY DATE  
(last day of the month)

NEGATIVE CONTROL

TEMPERATURE LIMIT  
(store)

BIOLOGICAL RISK

CONTENT IS SUFFICIENT  
FOR <N> TEST

FLAMMABLE

SEE INSTRUCTIONS  
FOR USE

CORROSIVE

IN VITRO DIAGNOSTIC  
PRODUCT

TOXIC

KEEP AWAY  
FROM SUNLIGHT

DO NOT USE IF  
PACKAGE IS  
DAMAGED

DO NOT REUSE

PRODUCT  
STERILIZED

CAUTION

DANGER