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INSTRUCTIONS FOR USE

FUNCTION

Method for determination of transferrin in biological samples of serum or plasma. Immunoturbidimetric test, for *in vitro* use only.

PRINCIPLES OF ACTION

Methodology: Immunoturbidimetry

Turbidimetric method in which the specific antibodies present in reagent form insoluble complexes with serum transferrin, generating turbidity, with intensity proportional to the quantity of transferrin in sample.

REAGENTS

Number 1 - Buffer - Store between 2 and 8°C. Contains: PEG 6.000 < 15%, buffer and preservative.

Number 2 - Antiserum - Store between 2 and 8°C. Contains: Goat antiserum for Transferrin < 3,8 g/L, buffer and preservative

PRESENTATION

Presentation	Reagent	
	Nº 1	Nº 2
1	25 mL	2 mL
2	2 x 25 mL	2 x 2 mL
5	1 x 40 mL	1 x 10 mL
7	20 mL	5 mL
8	2 x 20 mL	2 x 5 mL
9	3 x 20 mL	3 x 5 mL
10	4 x 20 mL	4 x 5 mL

EQUIPMENTS AND OPERATIONAL INPUTS

To perform this technique it is required the biochemistry automated equipment, Multical kit and Multicontrol Bioclin.

TRANSPORTATION AND STORAGE CONDITIONS

Store temperatures must be between 2 and 8°C. Transportation in temperatures up to 30°C should not exceed 5 days. Keep away from light and avoid humidity. **Do not freeze.**

SPECIAL CARE

- For professional *in vitro* diagnostic use only.
- Water used in material cleansing must be recent and free of contaminants.
- Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers that may alter the results significantly.
- Reagents Nº 1 and Nº 2 contain Sodium Azide, irritant to skin and mucous. Handle with care.
- Do not freeze the reagents, the freezing of antibody or buffer can affect their performance.
- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.
- To obtain information related to biosafety or in case of accidents with the product, consult the SDS (Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.
- Do not use the product in case of damaged packaging.
- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Fresh serum or plasma obtained with Heparin or EDTA. The analyte is stable for 2 days at 2 to 8°C and 6 months at -20°C ².

TO OBTAIN THE INSTRUCTIONS FOR USE IN PRINTED FORMAT, AT NO ADDITIONAL COST, CONTACT CUSTOMER ADVISORY SERVICE:

SAC: +55 (31) 3439 5454 / 0800 031 5454 / sac@bioclin.com.br

PROCESS DESCRIPTION

The calibration stability of the Transferrin kit installed on refrigerated equipment is at least 15 days. This stability may vary depending on the conditions of the test, equipment and environment. Therefore, it is suggested to follow the product performance using control serum.

REAGENTS PREPARATION

Reagents are ready to be used.

TECHNIQUE

For calibration and control of the reaction, use the Multical and Multicontrol Kits Bioclin.

This kit is indicated to use in automated biochemical analyzers only. **Verify equipment set up through the website www.bioclin.com.br or through SAC.**

INTERFERENCE

No interferences to Hemoglobin up to 1000 mg/dL, Sodium Citrate up to 1000 mg/dL, Heparin up to 50 mg/dL, Bilirubin up to 20 mg/dL, Triglycerides up to 2500 mg/dL and Rheumatoid Factor up to 1700 UI/mL.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

REFERENCE VALUES

Reference values were obtained through the determination of transferrin in healthy populations of male and female.

Reference Values: 170 - 340 mg/dL (IFCC)

These values must be used only as a reference. Each laboratory must create his own range of reference, according to the population served. The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

PRODUCT PERFORMANCE

ACCURACY

Comparison of Methods

The kit Transferrin was compared with other methods commercially available for measurement of transferrin. For this test, 41 clinical samples were analyzed. The linear equation $Y = 0.976X + 3.228$ with correlation coefficient equal to 0.992. With these results, we can conclude that the kit has good methodological specificity.

PRECISION

Repeatability

The repeatability was calculated from 10 successive determinations, using 2 samples with different concentrations, obtaining the following results:

Repeatability	Samples	
	1	2
Average Concentration (mg/dL)	217.57	129.33
Standard Deviation (mg/dL)	2.03	0.93
Coefficient of Variation (%)	0.93	0.72

Reproducibility

The reproducibility was calculated from 10 successive determinations for 3 consecutive days, using 2 samples with different concentrations, obtaining the following results:

Reproducibility	Samples	
	1	2
Average Concentration (mg/dL)	217.82	129.33
Standard Deviation (mg/dL)	0.34	1.13
Coefficient of Variation (%)	0.16	0.87

SENSITIVITY

The sensitivity was calculated from 10 determinations of a sample free of Transferrin. The average found, which indicates the method detection limit, was 1.50 mg/dL, with standard deviation of 0.08 mg/dL. The sensitivity, which indicates the method quantification limit, corresponds to the average plus 3 times the standard deviation and is equal to 1.74 mg/dL.

LINEARITY

The reaction is linear up to the level of the highest calibrator concentration.

HIGH DOSE HOOK EFFECT

Hook effect wasn't observed with a high concentration of Transferrin up to 2000 mg/L.

DIAGNOSTIC SIGNIFICANCE

Transferrin is a glycoprotein synthesized in the liver and is main plasma protein for iron transport. Its concentration is correlated with the total capacity of serum iron binding. Dosage levels of transferrin in plasma are useful for differential diagnosis of anemia and to monitor their treatment. Increased levels are found in iron deficiency in pregnancy, administration of estrogen and hereditary hemochromatosis. Decreased levels can be found in hereditary deficiencies, testosterone administration, infection, acute inflammation, some forms of necrosis, tumors, acute malaria and in cases of malnutrition.

BIBLIOGRAPHIC REFERENCES

1. Dati, F. et al., Lab. Med. 13, 87 (1989)
2. WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31
3. QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone.: +55 31 3439.5454 | E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service
Phone.: 0800 0315454 | E-mail: sac@bioclin.com.br

ANVISA registration for Transferrin kit: 10269360186

Review: March/2026

UNIVERSAL SYMBOLOGY

	CATALOG NUMBER		MADE BY
	LOT NUMBER		CONTROL
	MANUFACTURING DATE		POSITIVE CONTROL
	VALIDITY DATE (last day of the month)		NEGATIVE CONTROL
	TEMPERATURE LIMIT (store)		BIOLOGICAL RISK
	CONTENT IS SUFFICIENT FOR <N> TEST		FLAMMABLE
	SEE INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC PRODUCT		TOXIC
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED
	DO NOT REUSE		PRODUCT STERILIZED
	CAUTION		DANGER