# Bioclin

## TRIGLICÉRIDES MONOREAGENTE

REF K117

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#### FINALIDADE

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

### USAGE INSTRUCTIONS

#### FUNCTION

Method for determination of Triglycerides in biological samples of serum or plasma (Heparin or EDTA). Enzymatic colorimetric test, only for *in vitro* diagnostic use.

#### PRINCIPLE OF ACTION

**Methodology:** Enzymatic Colorimetric (Trinder)

O<sub>2</sub>, 4 Aminoantipyrine and p-Chlorophenol in the presence of Peroxidase result in a cherry colored compound whose color intensity is proportional to the concentration of Triglycerides.

#### REAGENTS

**Number 1 - Enzyme Reagent** - Store between 2 and 8°C. Contains: 4-Chlorine Phenol < 5 mmol/L, Lipoprotein Lipase < 5000 U/L, Glycerol Kinase < 3000 U/L, Peroxidase < 5000 U/L, Glycerol-3-Phosphate Oxidase < 5000 U/L, 4-Amineantipyrine < 1 mmol/L, ATP < 5 mmol/L, buffers, activators, stabilizer, surfactant and preservative.

**Number 2 - Standard** - Store between 2 and 8°C. Contains: Balstilled Glycerin 0,098 g/L (equivalent to 100,0 mg/dL Triglycerides) and preservative.

#### PRESENTATION

Presentation	Reagent N° 1	Reagent N° 2
1	1 x 100 mL	1 x 3 mL
2	2 x 100 mL	1 x 3 mL
3	4 x 100 mL	1 x 3 mL
4	5 x 20 mL	1 x 3 mL
5	5 x 40 mL	1 x 3 mL
6	10 x 40 mL	1 x 3 mL
7	2 x 60 mL	1 x 3 mL
8	4 x 60 mL	1 x 3 mL
9	6 x 60 mL	1 x 3 mL
10	8 x 60 mL	1 x 3 mL
11	10 x 60 mL	1 x 3 mL

#### EQUIPMENTS AND OPERATIONAL INPUTS

Spectrophotometer or colorimeter, water bath at 37°C, watch or stopwatch, pipettes, test tubes, Biocontrol N and Biocontrol P Bioclin. They can be found at markets specialized on Laboratories of Clinical Analysis.

#### TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be between 2 to 8°C. The transport at temperatures up to 30°C should not exceed 5 days. Protect from light and avoid moisture. **Do not freeze.**

#### SPECIAL CARE

**1- For *in vitro* diagnostic use only.**

**2-** Strictly follow the methodology proposed to obtain exact results.

**3-** Water used in material cleaning must be recent and free of contaminants.

**4-** Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers that may alter the results significantly.

**5-** Reagents N° 1 and N° 2 contains Sodium Azide, and should be handled with care.

**6-** The Sodium Hypochlorite is a contaminating agent that can significantly alter the results, so the materials used to perform the tests must be properly washed and free from this type of residue.

**7-** Water level in water bath should be above the level of reagents in test tubes.

**8-** We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

**9-** To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website [www.bioclin.com.br](http://www.bioclin.com.br) or upon request by the SAC (Customer Advisory Service) of Quibasa.

**10-** Do not use the product in case of damaged packaging.

**11-** It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

#### SAMPLES

Serum obtained free of hemolysis (to avoid falsely elevated results) or plasma collected with EDTA or Heparin. The compound is stable for up to 3 days between 2 to 8°C, and 30 days at -20°C. The blood can be collected without fasting or with a fasting period of 12 to 14 hours. Lipemic samples must first be diluted with Sodium Chloride 0,85%, in proportion 1:2.

#### PROCESS DESCRIPTION

The calibration stability of the Triglycerides Monoreagent kit installed on refrigerated equipment is at least 30 days. This stability may vary depending on the conditions of the test, equipment and environment. Therefore, it is suggested to follow the product performance using control serum.

#### TECHNIQUE

The Bioclin recommends, as control serum, Biocontrol N and P Bioclin kits.

Mark 3 test tubes: B (Blank), A (Sample), P (Standard), and proceed as follows:

	Blank	Standard	Sample
Sample	---	---	10 µL
Reagent N°1	1.0 mL	1.0 mL	1.0 mL
Reagent N°2	---	10 µL	---

Homogenize well and put in a water bath at 37°C for 10 minutes. Read the absorbance of Sample and Standard at 500 nm (490 - 540 nm), hitting the zero with the Blank. The color is stable for 30 minutes.

#### CALCULATIONS

$$\text{Triglycerides (mg/dL)} = \frac{\text{Sample Absorbance} \times 100}{\text{Standard Absorbance}}$$

As the reaction follows the Beer-Lambert Law, the Calibration Factor can be used.

$$\text{Calibration Factor} = \frac{\text{Standard Concentration (100 mg/dL)}}{\text{Standard Absorbance}}$$

$$\text{Triglycerides (mg/dL)} = \text{Sample Absorbance} \times \text{Calibration Factor}$$

Results are expressed in mg/dL.

#### INTERFERENCES

Reducing substances such as Ascorbic Acid, even at low concentrations, and icteric samples with Bilirubin above 5 mg/dL, interfere in the methodology, leading to falsely decreased results. Some substances such as Alcohol, Oral Contraceptives and Estrogens raise Triglyceride values.

#### INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

#### TRACEABILITY

The kit's standard is traceable to the reference material NIST (National Institute of Standards and Technology) SRM 1951.

#### REFERENCE VALUES <sup>6</sup>

The reference values in mg/dL, for this method, were obtained through the determination of Triglycerides in healthy populations of male and female.

Reference values below are for individuals with 12 hours fasting and follow Adult Treatment Panel III (ATP III).

Adults (≥ 20 years)	Desirable	< 150 mg/dL
	Limit	150 – 199 mg/dL
	Elevated	200 – 499 mg/dL
	Very Elevated	≥ 500 mg/dL

Childrens and Adolescents (2 to 19 years) <sup>7</sup>	Desirable	<100 mg/dL
	Limit	100 - 129 mg/dL
	Elevated	≥ 130 mg/dL

To convert the values from mg/dL to mmol/L (SI), multiply by 0.0113.

These values should be used as a guideline only, and each laboratory should establish its range of reference values, according to the population served.

The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

## PRODUCT PERFORMANCE

### ACCURACY

#### Comparison of Methods

Triglycerides Monoreagent kit was compared with another method for measurement of Triglycerides commercially available. 42 tests were performed and the results were evaluated. The linear equation obtained was  $Y = 0.9694X + 3.5955$  and coefficient of correlation 0.9978. With these results we can conclude that the kit shows good methodological specificity.

### PRECISION

#### Repeatability

The repeatability was calculated from 40 successive determinations, using 3 samples with different concentrations, obtaining the following results:

Repeatability	Sample 1	Sample 2	Sample 3
Average (mg/dL)	102.47	173.31	139.40
SD (mg/dL)	0.63	1.08	0.75
CV (%)	0.61	0.62	0.53

#### Reproducibility

The reproducibility was calculated from 40 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

Reproducibility	Sample 1	Sample 2	Sample 3
Average (mg/dL)	103.00	176.67	139.87
SD (mg/dL)	0.48	1.36	0.54
CV (%)	0.46	0.78	0.38

### SENSITIVITY

The sensitivity was calculated from 40 determinations of a sample free of Triglycerides. The average was 1.34 mg/dL, with standard deviation of 0.41 mg/dL. The sensitivity that indicates the method detection limit, corresponds the average plus 3 times the standard deviation, and is equal to 2.58 mg/dL.

### LINEARITY

The reaction is linear up to the concentration of 900 mg/dL. For samples with values above 900 mg/dL dilute the sample with Sodium Chloride 0.85%, repeat the dosage and multiply the result by dilution factor.

### DIAGNOSTIC SIGNIFICANCE

Triglycerides, constituents of many lipoproteins, are found in different concentrations and are of great importance classification and phenotyping of hyperlipoproteinaemia. In several pathologies that occur in hyperlipidemia, triglycerides only are not elevated in type IIa. Values are checked increased in several diseases such as diabetes, diseases cardiovascular disease, pancreatitis, nephrotic syndrome, uremia, hypothyroidism, chronic alcoholism. In screening procedures, the determination of triglycerides, cholesterol, evidence of cooling serum and other parameters can provide consistent data to classification and phenotyping of hyperlipoproteinaemia.

### BIBLIOGRAPHIC REFERENCES

- 1 - BUCOLO, G.; DAVID, H., Clin. Chem., 1973, 476.
- 2 - TONKS, D. B., Quality Control in Clinical Laboratories, 1983.
- 3 - MC GOWAN, M. W.; ARTISS, J. D.; STRANDBERGH, D. R.; ZAK, B., Clin. Chem., 1983, 29-538.
- 4 - TRINDER, P.; Ann. Clin. Biochem, 1969, 6-24.
- 5- Reunião Conjunta – Laudos Laboratoriais. SBC/DA, SBAC, SBPC/ML, SBBM; 15 de Outubro de 2013.
- 6- Executive Summary of the Third report of the National Cholesterol Education Program (NCEP) Expert Panel on Detection, Evaluation and Treatment of High Blood Cholesterol in Adults (Adult Treatment Panel III). JAMA 2001: 285:2486-97.
- 7 - Expert Panel on Integrated Guidelines for cardiovascular H, Risk Reducon in C, Adolescents, Naonal Heart L, Blood I. Expert panel on integrated guidelines for cardiovascular health and risk reducon in children and adolescents: summary report. Pediatrics. 2011; 128 Suppl 5: S213&56.
- 8- QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

### QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

#### QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca  
CEP 31565-130 - Belo Horizonte - MG - Brasil  
Phone: +55 31 3439.5454 | E-mail: bioclin@bioclin.com.br  
CNPJ: 19.400.787/0001-07 - Made in Brazil

### CUSTOMER SERVICE

Customer Advisory Service  
Phone: 0800 0315454 | E-mail: sac@bioclin.com.br

ANVISA registration for Triglycerides Monoreagent kit:  
10269360193

Review: February/2023

## UNIVERSAL SYMBOLOGY

CATALOG NUMBER

MADE BY

LOT NUMBER

CONTROL

MANUFACTURING DATE

POSITIVE CONTROL

VALIDITY DATE  
(last day of the month)

NEGATIVE CONTROL

TEMPERATURE LIMIT  
(store)

BIOLOGICAL RISK

CONTENT IS SUFFICIENT  
FOR <N> TEST

FLAMMABLE

SEE INSTRUCTIONS  
FOR USE

CORROSIVE

IN VITRO DIAGNOSTIC  
PRODUCT

TOXIC

KEEP AWAY  
FROM SUNLIGHT

DO NOT USE IF  
PACKAGE IS  
DAMAGED

DO NOT REUSE

PRODUCT  
STERILIZED

CAUTION

DANGER