Bioclin

TTPA BIOCLIN

REF **K090**

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FINALIDADE

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INSTRUCTIONS FOR USE

FUNCTION

Method for measuring manual or automatic time of Activated Partial Thromboplastin Time in biological samples of citrated plasma. For *in vitro* diagnostics use only.

PRINCIPLE OF ACTION

Methodology: Coagulometry

Consist on the determination of time of citrated plasma after adding reagents that contain a Plasmatic Activator (Ellagic Acid) and Phospholipids, that will act as substitutes to platelets. Initially plasma is incubated with the Contact Activator, at 37°C for 3 minutes. Following, add the Clot Forming (Calcium Chlorite) for recalcification, triggering the intrinsic coagulation mechanism. This step should be timed until the formation of the clot.

REAGENTS

Number 1 - Contact Activator (Activated Cephalin) - Store between 2 and 8°C. Contains: Rabbit brain extract, Ellagic Acid and preservative.

Number 2 - Clot Former (Calcium Chlorite Solution) - Store between 2 and 8°C. Contains: Calcium Chloride and preservative.

Presentation ready for use.

PRESENTATION

Presentation	Reagent 1	Reagent 2
1	6 x 2.5 mL	4 x 4 mL
2	6 x 5 mL	8 x 4 mL
3	12 x 2.5 mL	8 x 4 mL
4	12 x 5 mL	16 x 4 mL

EQUIPMENTS AND OPERATIONAL INPUTS

Coagulometer, water bath at 37°C, watch or stopwatch, pipettes, test tubes, Biocontrol N Coagulation and Biocontrol P Coagulation. These itens are found at markets specialized on Laboratories of Clinical Analysis.

TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be 2 to 8°C. Transport, at temperatures up to 30°C, should not exceed 5 days. Keep away from light and avoid moisture. **Do not freeze.**

SPECIAL CARE

1- For *in vitro* diagnostic use only.

2- Strictly follow the methodology proposed to obtain exact results.

3- The level of water in water bath must be superior to the level of reagents in test tubes.

4- Reagent Nº 1 contains Sodium Azide, irritant to skin and mucous, handle with care.

5- Avoid the warming of reagent that will not be used. Already heated reagent should not be returned to the original bottle.

6- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.

7- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.

8- Do not use the product in case of damaged packaging.

9- It is essential that the instruments and equipmants used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Material: Plasma

Anticoagulant: Use as coagulant the Bioclin Citrate kit, according to instructions of usage, or Sodium Citrate 3.8% (0.130 M) or 3.2% (0.109 M), proportionally in 9 parts of Blood to 1 of Anticoagulant (for example: 4.5 mL of Blood + 0.5 mL of Anticoagulant).

Collection: Blood must be obtained through vein puncture, avoiding hemolysis, prolonged tourniquet, bubble formation and tissue liquid aspiration. When evaluating only the coagulation test, always collect 2 tubes of citrate and discard the first one

Preparation: Blood should be mixed with anticoagulant right after collection. Centrifuge immediately at 3000 rpm, during 15 minutes.

Remove plasma without pipetting red cells and the yellow layer.

Samples must be tested in less than 3 hours. If tests can not be made during this period plasma should be frozen for at maximum 1 week at -20°C.

TECHNIQUE

Bioclin recommends, for use of the kit, to use the Biocontrol N Coagulation and Biocontrol P Coagulation kits as the control serum.

1- Preheat the Reagent Nº 2 to 37°C. This procedure should be performed only with the volume that will be used to perform the tests. The reagent must remain warm for the entire procedure.

2- Separate the reaction tubes that will be used to perform the tests.

3- Pipette 100 µL of the Sample/Control into the respective reaction tube and preheat to 37°C for 2 minutes in a water bath.

4- Pipette 100 µL of Reagent Nº 1 into the Sample/Control tube and incubate at 37°C, for 2 to 3 minutes.

5- Pipette 100 µL of Reagent Nº 2 (preheated to 37°C) in the Sample/Control tube, simultaneously triggering the timer.

6- Leave the tube in the water bath at 37°C, shaking gently. Keep in the bath for about 20 seconds.

7- Remove the tube from the bath, gently tilt it once a second and stop the timer at the time of the clot formation. Note the clotting time.

CALCULATIONS

Calculate the average between coagulation time of each samples/controls. Release all results in seconds.

CURVE OF SENSIBILITY TO HEPARIN

This procedure is useful to monitor therapies with Heparin.

A curve of sensibility to Heparin must be built adding the quantity of Heparin known to a mixture of normal plasmas (pool).

As a suggestion the following values of Heparin can be used:

0.0 / 0.1 / 0.2 / 0.4 / 0.6

Determine the dilution time dosages, plotting in a graphic the times and concentrations.

Each laboratory must establish its own calibration curve.

PROCEDURE LIMITATIONS

The concentration of citrate may interfere with coagulation time because it directly affects the concentration of calcium. Thus, APTT values are higher in tubes of 3.8% than in tubes of 3.2%.⁴

Citrate at 3.8% concentration will release results higher than 3.2%. Therefore, it is necessary that the laboratory standardize so that the monitoring of the patient's diagnosis is always performed with the same anticoagulant concentration.

The concentration of sodium citrate should be adjusted in patients with hematocrit values above 55%. In samples with high hematocrit, the blood/anticoagulant ratio (9: 1) is not maintained and may cause excess citrate to the volume of blood present in the tube. To calculate the citrate volume, use the formula described in item 4 of the Laboratory Diagnostic Manual for Hereditary Coagulopathies and Plaquetopathies of the Ministry of Health.⁵

INTERFERENCES

Icteric, lipemic or hemolyzed samples.

Usage of plasma containing Heparin, EDTA or Sodium Oxalate.

Conditions and substances that may alter TTPA: Oral contraceptives, Estrogen, pregnancy, drugs Coumarin, Heparin, Asparaginase, Acetylsalicylic Acid, Metronidazole, Digitalis and Tetracycline.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

REFERENCE VALUES

Interval should be used only as guidance. It is recommended that each laboratory establish, in population served, its own reference range, using a pool of normal plasmas or reference plasmas.

Age	Interval (seconds)
2 months	26 – 47
5 months	26 – 46
Children older than 6 months and Adults	27 - 38

The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

PRODUCT PERFORMANCE

ACCURACY

Comparison of Methods and Methodological Specificity

APTT Bioclin kit was compared with other method commercially available for dosage of Activated Partial Thromboplastin Time. 07 analyzes were performed and the results were evaluated. The linear equation obtained was Y = 1.0094X - 1.5788 and the correlation coefficient 0.9974. With these results, we can conclude that the kit shows good methodological specificity.

PRECISION

Repeatability

The repeatability was calculated from 20 successive determinations, using 3 samples with different concentrations, obtaining the following results:

Repeatability	Sample 1	Sample 2	Sample 3
Time (s)	32.15	64.90	44.55
Standard Deviation (s)	1.14	0.97	0.94
Coefficient of Variation (%)	3.54	1.49	2.12

Reproducibility

The reproducibility was calculated from 20 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

Reproducibility	Sample 1	Sample 2	Sample 3
Time (s)	32.67	64.32	44.65
Standard Deviation (s)	0.48	0.93	0.22
Coefficient of Variation (%)	1.47	1.44	0.49

SENSITIVITY

The kit has been evaluated for its sensitivity to the factors of deficiency in coagulation and obtained the following results:

Factor	VIII		IX		XI		XII		Pre-kallikr.
% Activ.	< 1	20	< 1	20	< 1	20	< 1	20	< 1
Result. in sec.	82.3	44.7	83.7	41.0	134.3	47.6	< 200	36.1	69.4

SPECIFICITY

The specificity of the kit is guaranteed by the addition of one single activation source of coagulant factors, Ellagic Acid.

DIAGNOSTIC SIGNIFICANCE

The Activated Partial Thromboplastin Time (APTT) is a test that monitors both the intrinsic pathway and a common pathway of coagulation and is useful in monitoring of Heparin therapy and to verify the deficiency of all coagulation factors except factor VII. Are widely used in deficiencies identified in stage 1 of the coagulation factors VIII, IX, XI, XII, Fletcher factor.

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QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone.: +55 31 3439.5454 | E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service
Phone.: 0800 0315454 | E-mail: sac@bioclin.com.br

ANVISA registration number for APTT Bioclin kit: 10269360152

Review: August/2023

UNIVERSAL SYMBOLOGY

CATALOG NUMBER

LOT NUMBER

MANUFACTURING DATE

VALIDITY DATE
(last day of the month)

TEMPERATURE LIMIT
(store)

CONTENT IS SUFFICIENT
FOR <N> TEST

SEE INSTRUCTIONS
FOR USE

IN VITRO DIAGNOSTIC
PRODUCT

KEEP AWAY
FROM SUNLIGHT

DO NOT REUSE

CAUTION

MADE BY

CONTROL

POSITIVE CONTROL

NEGATIVE CONTROL

BIOLOGICAL RISK

FLAMMABLE

CORROSIVE

TOXIC

DO NOT USE IF
PACKAGE IS
DAMAGED

PRODUCT
STERILIZED

DANGER