## URÉIA UV

REF **K056**

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## UREA UV

TO OBTAIN THE INSTRUCTIONS FOR USE IN PRINTED FORMAT, AT NO ADDITIONAL COST,  
CONTACT CUSTOMER ADVISORY SERVICE:

REF **K056**

SAC: +55 (31) 3439 5454 / 0800 031 5454 / sac@bioclin.com.br

### INSTRUCTIONS FOR USE

#### FUNCTION

Method for determination of the Urea in biological samples of serum, plasma (EDTA, Heparin or Citrate) and urine. Kinetic test, only for *in vitro* diagnostic use.


#### PRINCIPLE OF ACTION

**Methodology:** Fixed Time Kinetic

Urea is hydrolyzed to NH<sub>3</sub> and CO<sub>2</sub> by Urease.

The Glutamate Dehydrogenase (GLDH) in the presence of NH<sub>3</sub> and α-Ketoglutarate, oxidizes NADH to NAD<sup>+</sup>.

The oxidation of NADH to NAD<sup>+</sup>, as measured by the decrease in absorbance is proportional to the concentration of urea in the sample.

#### REAGENTS

**Number 1 - Buffer** - Store between 2 and 8°C. Contains: Tris Buffer < 100 mmol/L, NADH < 3 mmol/L and preservative.

**Number 2 - Enzymatic Reagent** - Store between 2 and 8°C. Contains: Tris Buffer < 500 mmol/L, ADP < 5 mmol/L, α-Ketoglutarate < 100 mmol/L, Urease < 50 KU/L, Glutamate Dehydrogenase < 5 KU/L, stabilizers and preservative.

**Number 3 - Standard** - Store between 2 and 8°C. Contains: Urea 70,0 mg/dL and stabilizer.

#### PRESENTATION

Presentation	Reagent N° 1	Reagent N° 2	Reagent N° 3
1	2 x 80 mL	1 x 40 mL	1 x 3 mL
2	1 x 40 mL	1 x 10 mL	1 x 3 mL
3	2 x 40 mL	2 x 10 mL	1 x 3 mL
4	4 x 40 mL	4 x 10 mL	1 x 3 mL
5	2 x 40 mL	1 x 20 mL	1 x 3 mL
6	4 x 40 mL	2 x 20 mL	1 x 3 mL
7	4 x 60 mL	4 x 15 mL	1 x 3 mL
8	5 x 60 mL	5 x 15 mL	1 x 3 mL
9	6 x 60 mL	6 x 15 mL	1 x 3 mL

#### EQUIPMENTS AND OPERATIONAL INPUTS

Thermostated spectrophotometer, pipettes, watch or stopwatches, test tubes, Biocontrol N and Biocontrol P Bioclin. They can be found at markets specialized on Laboratories of Clinical Analysis.

#### TRANSPORTATION AND STORAGE CONDITIONS

The storage temperature should be between 2 to 8°C. The transport at temperatures up to 30°C should not exceed 5 days. Protect from light and avoid moisture. **Do not freeze.**

#### SPECIAL CARE

- 1- For *in vitro* diagnostic use only.**
- Strictly follow the methodology proposed to obtain exact results.
- Water used in material cleaning must to be recent and free of contaminants.
- Saturated deionizer columns release alkaline water, many ions, oxidizing agents and reducers that may alter the results significantly.
- Handle Reagents N° 1 and N° 2 with care, for it contains Sodium Azide which can be irritating to skin and mucous membranes.
- Determinate the factor periodically and each batch of product.
- Heavy metals and fluoride (above 2 mg/dL) are inhibitors of Urease.
- No smoking near the place where they carry out the measurement, because the smoke contains ammonia fumes that contaminate the sample and glassware, leading to falsely increase.
- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.
- To obtain information related to biosafety or in case of accidents with the product, consult the SDS (Safety Data Sheet) available on the website [www.bioclin.com.br](http://www.bioclin.com.br) or upon request by the SAC (Customer Advisory Service) of Quibasa.
- Do not use the product in case of damaged packaging.
- It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

#### SAMPLES

Serum, plasma (EDTA, Heparin or Citrate) and Urine.

Do not use anticoagulants containing ammonia.

The Serum or Plasma sample is stable for 7 days at 2 to 8°C, and 3 months at -20°C.

Urine sample is stable for 7 days at 2 to 8°C and 4 weeks at -20°C.

For determination of Urea in the urine sample collected 24 hours in a vial containing 2.0 mL of HCl 50% (v/v). Centrifuge the sample before starting the technique and to conduct analysis within a few hours, because the urea excreted in urine is easily decomposed by bacterial action.

#### PROCESS DESCRIPTION

The calibration stability of the Urea UV kit installed on refrigerated equipment is up to 07 days. This stability may vary depending on the conditions of the test, equipment and environment. Therefore, it is suggested to follow the product performance using control serum.

#### PREPARATION OF WORKING REAGENT

**The reagents are ready to use in bireagent format and on automated equipment that uses dedicated presentations.**

For use in monoreagent format, prepare the working reagent as follows

Mix 4 parts of Reagent N° 1 (Buffer) with 1 part Reagent N° 2 (Enzymatic Reagent). The Working Reagent is stable for 05 days between 15 and 30°C and 30 days between 2 and 8°C.

#### REACTION CONDITIONS

Is indispensable conditions the use a thermostated cuvette at 37°C (25 to 30°C), 1 cm light path and read at 340 nm (334 - 365 nm).

#### TECHNIQUE

We recommend using the Bioclin Biocontrol N and Biocontrol P kits as control serums.

#### Monoreagent:

Place in a test tube 1.0 mL of Working Reagent and add 10 µL of Sample or Standard, mix and transfer to a thermostated cuvette at 37°C. Trigger the stopwatch simultaneously and measure the absorbance at 30 and 90 seconds at 340 nm (334 - 365 nm).

Differences in absorbance (ΔA) between two times, the Standard and the Sample, will be used to calculate the results.

For determination of Urea in urine, following the above technique use sample diluted 1:50 (0.1 mL urine + 4.9 mL distilled or deionized water). Multiply the result by 50.

To calibrate the bireagent kit, Bioclin recommends using the BIOCAL Bioclin calibrator serum.

**Check the equipment's programming by calling Customer Service at 0800 031 5454**

#### CALCULATIONS

$$\text{Urea (mg/dL)} = \frac{\Delta A \text{ of Sample}}{\Delta A \text{ of Standard}} \times 70$$

As the reaction follows the Beer-Lambert Law, the calibration factor can be used.

$$\text{Calibration Factor} = \frac{\text{Standard Concentration (70 mg/dL)}}{\Delta A \text{ of Standard}}$$

$$\text{Urea (mg/dL)} = \Delta A \text{ of sample} \times \text{Calibration Factor}$$

#### Urine

$$\text{Urea (g/24h)} = \frac{\Delta A \text{ of Sample} \times \text{Factor} \times 50 \times \text{Volume (L)}}{100}$$

The results for serum or plasma are expressed in mg/dL and urine g/24h.

#### PROCEDURE LIMITATIONS

Contamination of water, glassware, environment and/or the sample with ammonia, may lead to falsely elevated results.

#### INTERFERENT

No interference was observed by Bilirubin until 20 mg/dL, Hemoglobin until 300 mg/dL and Triglycerides until 1800 mg/dL. Some drugs such as Aminoglycosides, Cephalosporins, Purinol, Methylidopa, Furosemide, Propranolol can produce interference in the results, increasing the values of Urea.

#### INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

**TRACEABILITY**

The kit standard and the BIOCAL calibrator are traceable to the NIST (National Institute of Standards and Technology) SRM 912 reference material.

**REFERENCE VALUES**

The reference values in mg/dL, for this method were obtained through the determination of urea in healthy populations of male and female.

Serum or plasma	
1 day - 12 month	2 - 34 mg/dL
1 - 13 years	8 - 36 mg/dL
Adults	15 - 40 mg/dL
Urine	
All ages	26 - 43 g/24h

To convert the values to mg/dL to mmol/L (SI), multiply by 0.166.

These values should be used as guidelines, each laboratory should establish its range of reference values, according to the population served.

The results provided by this kit must be interpreted by the medical professional responsible, not being the only criterion for the determination of diagnosis and/or treatment of the patient.

**PRODUCT PERFORMANCE****ACCURACY****Comparison of Methods**

The Urea UV kit was compared with another method for dosage of Urea, commercially available. 42 analyzes were conducted and the results were evaluated. The linear equation obtained was  $Y = 0.9655X + 2.2447$  with correlation coefficient equal to 0.9975. With these results, we can conclude that the kit shows good specificity methodology.

**PRECISION****Repeatability**

The repeatability was calculated from 40 successive determinations, using 3 samples with different concentrations, obtaining the following results:

Repeatability	Sample 1	Sample 2	Sample 3
Average Concentration (mg/dL)	33.75	24.00	41.03
Standard Deviation (mg/dL)	0.95	0.60	1.46
Coefficient of Variation (%)	2.83	2.50	3.56

**Reproducibility**

The reproducibility was calculated from 40 successive determinations for 3 consecutive days, using 3 samples with different concentrations, obtaining the following results:

Reproducibility	Sample 1	Sample 2	Sample 3
Average Concentration (mg/dL)	33.36	23.89	40.12
Standard Deviation (mg/dL)	0.34	0.13	0.79
Coefficient of Variation (%)	1.02	0.54	1.96

**SENSITIVITY**

The sensitivity was calculated from 40 determinations of

a sample free of Urea. The average found that indicates the method detection limit was 0.175 mg/dL, with standard deviation of 0.446 mg/dL. The sensibility, which indicates the quantification limit of the method, corresponds the average plus 3 times the standard deviation, and is equal to 1.514 mg/dL.

**LINEARITY**

The reaction is linear up to 300 mg/dL. For samples with values above 300 mg/dL, dilute it with distilled or deionized water, repeat the dosage and multiply the result by the dilution factor.

**DIAGNOSTIC SIGNIFICANCE**

Urea, the main product of the catabolism of proteins and amino acids, has its serum concentration affected by diet and hydration status, providing a rough indication of the state of renal function.

Elevations values of plasma Urea are classified as:

Pre-Cause kidney: excretion resulting from defects observed in cardiac decompensation, hemorrhagic shock, acute dehydration, high protein catabolism (burns, fever).

Kidney Cause: as a consequence of acute or chronic kidney disease (glomerulonephritis, pyelonephritis, tubular necrosis) with levels of Urea plasma of 300 mg/dL or more.

Post-Cause Kidney: usually resulting from a urinary tract obstruction, can occur in kidney lithiasis tumors by compression of the bladder.

The decrease in serum Urea occurs only in a few situations such as acute liver failure, starvation, in the last trimester of pregnancy.

**BIBLIOGRAPHIC REFERENCES**

- 1 - BALLETER, W. D. E.; BUSITMAN, C. S.; TIDWELL, P.W., Anal. Chim., 33-59.
- 2 - WINDMANN, F. K.; TURNER, K., Clin. Chem, 1987, 21:1754-1770.
- 3 - BERGMAYER, HU., Methods of Enzymatic Analysis, vol. 9, VCH Publishers, 1985, 449-453.
- 4 - TALKER, H.; SCHUBERT, G. E., Klin, 1965, 43-174.
- 5 - CARL, A. B. and EDWARD R. A., Tietz Textbook of Clinical Chem. 2a Ed., 1994, 1528-1531.
- 6 - WHO. Use of anticoagulants in diagnostic laboratory investigations and stability of blood, plasma and serum samples. WHO/DIL/LAB/99.1 rev. 2, 2002:31.
- 7- QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

**QUALITY ASSURANCE**

Before being released for consumption, all Bioclin reagents are tested by the Quality Control Department. The quality of the reagents is ensured until the expiration date mentioned on the presentation packaging, provided they are stored and transported under the appropriate conditions.

**QUIBASA QUÍMICA BÁSICA Ltda**

Rua Teles de Menezes, 92 - Santa Branca  
CEP 31565-130 - Belo Horizonte/MG - Brasil  
Tel.: +55 31 3439 5454 | E-mail: bioclin@bioclin.com.br  
CNPJ: 19.400.787/0001-07 - Indústria Brasileira

**CUSTOMER SERVICE**

Customer Advisory Service  
Phone: 0800 0315454 | E-mail: customerservice@bioclin.com.br

ANVISA registration for Urea UV kit: 10269360068

Review: October/2025

**UNIVERSAL SYMBOLOGY**