LEISHMANIOSE IgG VET FAST

REF| VET004

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FINALIDADE

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USE INSTRUCTIONS

PURPOSE

Fast immunochromatographic test for the qualitative determination of anti-Leishmania antibodies in serum, plasma or whole blood samples.

PRINCIPLE OF ACTION

Methodology: Immunochromatography.
Leishmaniasis IgG VET FAST is an *in vitro* diagnostic test for the rapid and qualitative detection of anti-*Leishmania* antibodies in blood, serum or plasma samples. The cassette consists of a nitrocellulose membrane impregnated with the Leishmania recombinant antigen in test line. The conjugate consists of a mouse antibody anti-dog IgG conjugated to colloidal gold particles. When the sample is added to the cassette, it migrates by capillary action until it reaches the conjugate. In this region, the anti-*Leishmania* IgG antibodies present in the sample interact with the conjugate, forming the complex: dog IgG – conjugate. The sample continues to migrate through the immunochromatographic strip until it reaches the test line, where the dog IgG antibody also interacts with the immobilized antigen of *Leishmania* in the test line, forming the complex: *Leishmania* antigen – dog IgG – conjugate. This shows the formation of a pink-red line on the test line; indicating positive result. Sample excess continues to migrate through the test strip until it reaches the control line, which is composed of impregnated anti-conjugate antibody. Thus, the formation of a second line (control line) is observed; indicating that the test yielded a valid result.

REAGENTS

Reagent Number 1: Cassette with pipette in aluminized sachet - Store between 15 and 30°C. **Do not freeze.** It contains:

01 plastic cassette with 01 test strip containing

- Control region (C): Goat antibody anti-mouse antibody immobilized.
- Test region (T): Immobilized Leishmania recombinant antigen.
- Conjugate: Mouse antibody anti-dog IgG, colloidal gold conjugate.

Reagent Number 2: Diluent - Store between 15 and 30°C. **Do not freeze.** Contains: Buffer solution and preservative.

PRESENTATIONS

Reagent	Presentation	Number of Cassettes for Package	Diluent
Cassette + Diluent	VET004-1	10	1 x 2,5 mL
	VET004-2	25	1 x 5,0 mL
	VET004-3	10	10 x 0,5 mL
	VET004-4	25	25 x 0,5 mL

EQUIPMENTS AND OPERATIONAL INPUTS

Clock or stopwatch. They are in the specialized market of articles for Clinical Analysis Laboratories.

TRANSPORTATION AND STORAGE CONDITIONS

The storage and transport temperature of the product must be between 15 and 30°C. Keep away from direct sunlight. **Do not freeze.**

The cassette is sensitive to air humidity. Therefore, it should always be kept in its inviolate primary packaging (aluminized sachet) until the time of use.

SPECIAL CARE

1- Product for *in vitro* diagnosis, for veterinary use only.
2 - Strictly follow the methodology proposed for obtaining results.
3- Cassettes (test devices) do not contain contaminants; but since the samples are biological, the professional must strictly follow the safety rules and routines. The use of disposable gloves and other personal protective equipment is essential.
4- Use clean, dry containers to package biological samples.
5- Open the aluminum envelop to remove the cassette only at the time of use.
6- The kit cannot be frozen or exposed to high temperatures, since this will cause deterioration of its components.
7- The reagents, as well as the samples, should reach room temperature before starting the test.
8- Do not use reagents from different lots or kits.
9- Do not use kit beyond its expiration date. The expiry date is stated on the product label.
10- The cassette (test device) and the capillary tube are disposable and cannot be reused.
11- The kit materials that had contact with the sample of the patient must have a special destination. The capillary tube the cassettes (test devices) and diluent vials should be disposed of in organic waste.

12- Secondary packaging (box), Instruction for use, aluminum envelope and silica sachet can be disposed of in trash common.

13- We recommend to apply the local, state and federal environmental protection so that the disposal of reagents and biological material is done in accordance with the current legislation.
14- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or through request by the Customer Service Department (SAC) of the Quibasa.

SAMPLES

Perform the test with serum, plasma or whole blood.
Serum: Allow blood to clot and then centrifuge to separate serum.
Plasma: Collect whole blood in a tube containing anticoagulants such as heparin, citrate, or EDTA. Centrifuge the blood and separate the plasma.
Whole blood: Collect the blood sample in a tube containing EDTA, heparin, or citrate using standard laboratory procedures. Whole blood samples should be analyzed within 24 hours, can be stored at 2-8°C and never be frozen. Serum and plasma samples should be analyzed immediately after collection, but if necessary, they can be stored at 2-8°C for up to 24 hours or at -20°C for longer periods.

PROCESS FOR CANINE LEISHMANIASIS FAST TEST

1- The sample should be at room temperature before starting the test.
2- Remove the cassette (test device) from the protective packaging and place it on a clean, dry and level surface
3- Collect the sample using the disposable capillary tube that accompanies the kit.
4- Tighten the capillary tube above the marked trace and touch the sample. Relieve the pressure so the sample is aspirated. Aspirate until the trace marked to obtain the volume of 10 µL.
5- Dispense the sample by pressing the pipette again, touching the test membrane at the sample well.
6- Use 1 drop of serum / plasma (10 µL) or 2 drops of blood (20 µL) to perform the test.
Note: Alternatively, an automatic pipette can be used to collect and add the sample.
7- Hold the Diluent bottle (Reagent Nº 2) vertically and apply 2 drops (70 µL) in sample well.

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
8- Wait for the lines to form. Interpret the results between 15 and 20 minutes. Do not interpret after 20 minutes.

INTERPRETATION OF RESULTS

Reactive Test: Formation of two lines: one in the control region (C) and one in the test region (T). Do not interpret after 20 minutes.	
Not Reactive Test: Line formation only in control region (C). Do not interpret after 20 minutes.	
Invalid Test: The absence of line formation in the control region (C), with or without line formation in the test region (T), indicates procedure error or deterioration of the cassette. In this case, repeat the test using a new cassette.	

TEST LIMITATIONS

This product is intended for screening for Leishmaniasis in dogs. Interpretation should not be based on a single essay. Other confirmatory tests must be included before a sample is considered positive. Therefore, for the diagnosis of the disease, the professional must take into account not only the result of this test; but also, clinical signs, animal health conditions, and other diagnostic results that may be available.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal program of quality control, where procedures, standards, limits and tolerance for variations are clearly established. It is important to note that all measuring systems characteristic analytical variability, which should be monitored by the laboratories themselves. Therefore, it is the use of controls, which allow the assessment of the precision and accuracy of dosages.

PRODUCT PERFORMANCE

QUALITY CONTROL

CLINICAL SENSITIVITY AND SPECIFICITY

This product has been tested comparing with other methods. A total of 105 samples characterized by canine Leishmaniasis were analyzed; 51 samples were diagnosed

as positive and 54 were diagnosed as negative. The results show that the clinical sensitivity of the product is >99% and the clinical specificity is 96,30%.

METHOD		REFERENCE		Total
Cassette Leishmaniasis IgG VET FAST	Result	Reactive	Not reactive	
	Reactive	51	2	
	Not reactive	0	52	
Total result		51	54	105

Clinical sensitivity: (51/51): >99%
Clinical Specificity: (52/54): 96,30%
Precision: (51+52/51+54): 98,09%

REPEATABILITY

Repeatability was calculated from 10 successive determinations, using 2 different samples, obtaining the following results:

Samples	Nº of Repetitions	Expected Result	Found Result
Sample 1	10	Reactive	Reactive
Sample 2	10	Reactive	Reactive

REPRODUCIBILITY

Reproducibility was calculated from 10 successive determinations over 3 consecutive days, using 2 different samples, obtaining the following results:

Samples	Nº of Days	Nº of Repetitions	Expected Result	Found Result
Sample 1	03	10	Reactive	Reactive
Sample 2	03	10	Reactive	Reactive

BIBLIOGRAPHIC REFERENCES

1- QUIBASA: Data from the Research and Development Department.

QUALITY GUARANTEE

Before being released for consumption, all Quibasa reagents are tested by the Quality Control Department. The quality of the reagents is assured until the validity date mentioned on the packaged, as long as they are stored and transported under the appropriate conditions.

MANUFACTURER INFORMATION

QUIBASA QUÍMICA BÁSICA Ltda
Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Tel.: (31) 3439.5454
E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Industria Brasileira

CONSUMER SERVICE

Customer Support Service
Tel.: 0800 0315454
E-mail: sac@bioclin.com.br

Product licensed in the Ministry of Agriculture, Livestock and Supply since 30/12/2019 with the number 10.299/2019.

Technical manager: Dr. Hugo Vieira Fróes (CRMV/MG 8963).

Review: February/2020

UNIVERSAL SYMBOLOGY

	CATALOG NUMBER		MANUFACTURED BY
	BATCH CODE		CONTROL
	DATE OF MANUFACTURE		POSITIVE CONTROL
	USED BY (last day of month)		NEGATIVE CONTROL
	TEMPERATURE LIMITATION (store at)		BIOLOGICAL RISK
	CONTAINS SUFFICIENT FOR <N> TESTS		INFLAMMABLE
	CONSULT INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC DEVICE		POISON
	EUROPEAN AUTHORIZED REPRESENTATIVE		CE MARK
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED