VETLISA AIE IgG

REF VET037

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FINALIDADE

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USE INSTRUCTIONS

FUNCTION

Test for the qualitative determination of IgG antibodies against Equine Infectious Anemia Virus in equine serum or plasma, by immunoenzymatic assay, in microplate.

PRINCIPLE OF ACTION

Methodology: Enzyme immunoassay or immunoenzymatic The VETLISA AIE IgG is a solid phase immunoenzymatic assay based on the immunocapture principle for the qualitative detection of IgG antibodies against Equine Infectious Anemia Virus (EIAV) in equine serum or plasma. Antibodies against Equine Infectious Anemia Virus present in the sample bind to the purified recombinant antigen that covers the microplate, forming the immobilized complex: antigen - anti-EIAV antibodies. After the initial incubation, the microplate is washed to remove unbound materials. The Antibody-Enzyme Conjugate binds to the immobilized antigen – anti-EIAV antibody complex on the plate. After the second incubation, the microplate is washed and the Substrate is added and incubated. The intensity of the blue color produced by the addition of the Substrate is proportional to the amount of antibodies against the Equine Infectious Anemia Virus present in the samples. The Stop Solution is added to stop the reaction, promoting a color change to yellow, measured in a microplate reader.

REAGENTS

- 1- Sensitized Plate** - Store between 2 and 8°C. Contains: Microplate coated with Equine Infectious Anemia virus (EIAV) purified recombinant antigen.
- 2- Concentrated Conjugate** - Store between 2 and 8°C. Contains: Peroxidase-linked anti-equine IgG antibody.
- 3- Concentrated Wash** - Store between 2 and 8°C. Contains: Buffered solution and preservative.
- 4- Diluent** - Store between 2 and 8°C. Contains: Buffered solution and preservative.
- 5- TMB Substrate** - Store between 2 and 8°C. Contains: Solution containing tetramethylbenzidine (TMB) and preservative.
- 6- Stop Solution** - Store between 2 and 8°C. Contains: Hydrochloric Acid.
- 7- Negative Control** - Store between 2 and 8°C. Contains: IgG antibodies not reactive for AIE virus and preservative. **Potentially infectious.**
- 8- Positive Control** - Store between 2 and 8°C. Contains: Anti-VAIE Antibodies and preservative. **Potentially infectious.**
- 9- Plate Sealers**

PRESENTATIONS

Components	Presentation		
	1	2	3
	96 wells	192 wells	480 wells
1- Sensitized Plate	1 unit (96 wells)	2 units (2 x 96 wells)	5 units (5 x 96 wells)
2- Concentrated Conjugate	1 x 300 µL	2 x 300 µL	5 x 300 µL
3- Concentrated Wash Solution	1 x 20mL	2 x 20mL	5 x 20mL
4- Diluent	1 x 60 mL	2 x 60 mL	5 x 60 mL
5- Substrate TMB	1 x 12 mL	2 x 12 mL	5 x 12 mL
6- Stop Solution	1 x 12 mL	2 x 12 mL	5 x 12 mL
7- Negative Control	1 x 1 mL	2 x 1 mL	5 x 1 mL
8- Positive Control	1 x 1 mL	2 x 1 mL	5 x 1 mL
9- Plate Sealers	3 units	6 units	15 units

EQUIPMENTS AND OPERATIONAL INPUTS

Materials contained in the kit:

- Reagents described in the previous table.
- Instructions for use (manual).

Materials needed, but not contained in the kit:

- 1-** Set of calibrated pipettes, capable of dispensing volumes between 5 and 1000 µL, with a coefficient of variation of less than 1.5%.
- 2-** Repeater or multichannel pipette calibrated at 100 µL and 300 µL volume, with a coeffi cient of variation of less than 1.5% (optional).
- 3-** Microplate washing machine (optional).
- 4-** ELISA reader with absorbance capacity at 450 and 630 nm wavelength.
- 5-** Absorbent paper to dry the microcavities.
- 6-** Stopwatch or clock.
- 7-** Bottle to store the Washing Solution, after dilution.
- 8-** Distilled or deionized water.
- 9-** Quality Control Tools

TRANSPORTATION AND STORAGE CONDITIONS

The product storage temperature should be 2 to 8 °C. Transport can be carried out at room temperature (30°C) for up to 5 days. Keep away from light and avoid humidity. **Do not freeze.**

SPECIAL CARE

- 1- Product for *in vitro* diagnosis, for veterinary use only.**
- 2-** Strictly follow the methodology proposed to obtain accurate results.
- 3-** The professional must strictly follow safety rules and routines when handling biological samples. The use of disposable gloves and other personal protective equipment is essential.
- 4- IMPORTANT:** Before starting the test, allow all kit components to reach room temperature. Open the envelope containing the wells only after reaching room temperature.

- 5-** Do not mix reagents from different kits or lots. Do not use expired kit components.
- 6-** Do not eat, drink, or smoke at the test site.
- 7-** Do not pipette reagents or samples by mouth. Do not use the same tip to pipette different samples.
- 8-** Positive and negative controls should be retested for each new test performed.
- 9-** The Stop Solution contains Sulfuric Acid, which is a strong acid. Therefore, handle it with care.
- 10-** Handle the kit components with care to avoid contamination. Take special care with the substrate, which is a light sensitive solution. Use new or appropriately clean tips and containers for handling and do not allow them to be exposed to strong light during storage or during the incubation periods of the assay.
- 11-** Use contaminant-free freshly distilled or deionized water when preparing solutions.
- 12-** Always add the reagents in the same order to minimize the difference in reaction time between the microcavities.
- 13-** Before reading the plate, ensure that the bottom of the microcavities are clean and dry. Make sure there are no bubbles on the surface of the liquid.
- 14-** We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.
- 15-** To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa.
- 16-** Do not use the product in case of damaged packaging.
- 17-** It is essential that the instruments and equipment used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Use serum or plasma (EDTA or Heparin). Hemolyzed or highly lipemic samples should not be used. Samples may be refrigerated at between 2 and 8°C for a maximum of 5 days. If samples cannot be analyzed within 5 days, they can be stored for up to 30 days at -20°C (freezer).

PROCESS DESCRIPTION

ATTENTION: The Positive and Negative Controls are ready to use.

PREPARATION OF WORK REAGENTS AND SAMPLES

1- Washing Solution
Dilute the content of the Concentrated Wash Solution (Reagent N° 3) in 1000 mL of freshly distilled or deionized water. Store at 2 to 8°C until the validity date printed on the original bottle. It can be stored at room temperature. If crystallization occurs, heat to 37°C until dissolution.

2- Conjugate Dilution

Dilute the Concentrated Conjugate (Reagent N° 2) in the ratio of 1/50 in Diluent (Reagent N° 4). Prepare the solution at the moment of performing the test. To perform an assay using all wells in the kit, mix 220 µL of Concentrated Conjugate in 10,8 mL of Diluent. To perform an assay using 8 wells (1 strip), mix 20 µL of the Concentrated Conjugate in 980 µL of Diluent. **IMPORTANT:** The diluted conjugate solution cannot be stored. Therefore, prepare only the amount necessary to perform the test.

3- Sample Dilution

In a test tube, dilute 10 µl of the sample in 490 µl of diluent. Cap the tube and gently vortex or manually mix by inversion. Dilutions cannot be stored.

TECHNIQUE

- Before starting the assay, place all reagents, Controls and Samples to stabilize at room temperature (15 - 30°C) for at least 30 minutes. Return the unused strips to the original sealed packaging.
- 1-** Separate the wells to be used considering: Controls and Samples (it is recommended to test in duplicate).
- 2-** Separate the first well for the Blank (OPTIONAL).
- 3-** Pipette 100 µL of the Negative Control and the Positive Control into the previously determined wells. **Note: The controls are ready to use, there is no need to dilute them.**
- 4-** Pipette 100 µL of the previously diluted samples into the previously determined wells. In the Blank cavity (OPTIONAL), if you have chosen, pipette only 100 µL of the diluent.
- 5-** Homogenize gently for ± 10 seconds. Cover the cavities with plate sealer.
- 6-** Incubate for 30 minutes at room temperature.
- 7-** Discard the content of the cavities by aspiration (washing machine) or by decantation (manual). Use approximately 300µL/well of previously diluted Wash Solution, for a total of three (3) wash cycles. Shake for three seconds after each wash cycle. To guarantee the drying of the plate, at the end of the wash, tap the plate for a few seconds on absorbent paper.
- Note:** Poor washing/drying can cause inadequate results.
- 8-** Pipette 100 µL of the Conjugate previously diluted into each well, even in the Blank cavity.
- 9-** Homogenize gently for ± 10 seconds. Cover the cavities with plate sealer.
- 10-** Incubate for 30 minutes at room temperature.
- 11-** Repeat item 7.
- 12-** Add 100 µL of Substrate TMB into each well.
- 13-** Homogenize gently for ± 3 seconds. Cover the cavities with plate sealer.
- 14-** Incubate for 10 minutes at room temperature, protected from light.
- 15-** Remove the plate sealer from the cavities.
- 16-** Pipette 100 µL of Stop Solution into all wells.

17- Homogenize gently for ± 3 seconds.
18- Read the absorbance in an ELISA reader on a dual filter 450 nm (primary filter) / 630nm (secondary filter) within a maximum of 10 minutes after adding the Stop Solution.

TECHNIQUE VERIFICATION

Verify if the results obtained from reading the Blank and the Controls are compatible with the values presented above:

ITEM	ABSORBANCE (DUAL FILTER)
Blank	< 0,150
Negative Control	> 0,100 e < 0,350
Positive Control	> 1,000

If the values are out the expected ones, the test should be repeated.

CUT OFF CALCULATION

If the results of the controls are valid, calculate the cut off value with the following formula:

Cut Off = Negative Control Average Absorbance + 0,300

CALCULATION OF THE SAMPLE INDEX

Calculate the index by dividing the absorbance of the sample by the cut off value. Example:

ITEM	ABSORBANCE
Sample	1,1
Negative Control	0,210
Cut Off Value	0,210 + 0,300 = 0,510
Index: Sample Abs. / Cut Off Value	1,1 / 0,510 = 2,156

INTERPRETATION OF THE RESULT - QUALITATIVE

RESULTS	INDEX
Negative	< 0,9
Positive	≥ 0,9 e ≤ 1,1
Undetermined	> 1,1

PROCEDURE LIMITATIONS

The interpretation of a diagnostic test should not be based on a single trial. Other confirmatory tests must be included before a sample is considered positive. A negative result does not exclude the possibility of exposure. However, all results should be interpreted in conjunction with vaccination history, available clinical and laboratory information.

INTERNAL QUALIT CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present an analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

PRODUCT PERFORMANCE
QUALITY CONTROL

CLINICAL SENSITIVITY AND SPECIFICITY

This product was tested comparing to other enzyme immunoassay methods. A total of 92 samples were analyzed, with 46 samples diagnosed as positive and 46 diagnosed as negative. The results show that the clinical sensitivity of the product is 91,3% and the clinical specificity is 95,6%.

METHOD		REFERENCE		TOTAL
VETLISA AIE IgG	Result	Positive	Negative	
	Positive	42	2	44
	Negative	4	44	48
TOTAL RESULT		46	46	92

Clinical sensitivity: 91,3% (42/46)
Clinical specificity: 95,6% (44/46)
Precision: 93,4% [(42+44) / (46+46)]

Precision

REPEATABILITY

Repeatability was calculated from 10 successive determinations, using 3 samples with different values, obtaining the following results:

REPEATABILITY	SAMPLE		
	1	2	3
Mean	1,009	1,373	0,053
Standard Deviation	0,061	0,082	0,006
Coefficient of Variation (%)	6,08	6,00	11,14

REPRODUCIBILITY

Reproducibility was calculated from 10 successive determinations over 3 consecutive days, using 3 samples with different values, obtaining the following results:

REPRODUCIBILITY	SAMPLE		
	1	2	3
Mean	1,019	1,364	0,051
Standard Deviation	0,056	0,069	0,006
Coefficient of Variation (%)	5,49	5,06	10,98

CLINICAL SIGNIFICANCE

Equine Infectious Anemia (EIA) is a cosmopolitan disease that affects the equines (horses, mules and donkeys). EIA is caused by a Lentivirus of the Retroviridae Family, mechanically transmitted by flies and horseflies such as *Stomoxys calcitrans*, *Chrysops* sp. and *Tabanus* sp., or through contaminated fomites. EIA virus infection causes fever, anemia, hemorrhage, weight loss, depression, changes varied neurological and immunosuppression that may appear acutely, chronically or intermittently. Asymptomatic animals, which act as reservoirs and spreaders of the virus, are common. The AIE is part of the National Equine Health Program (PNSE) in that in the occurrence of suspicious cases the Agricultural Defense Service must be notified.

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QUALITY ASSURANCE

Before being released for consumption, all Quibasa reagents are tested by the Quality Control Department. The quality of the reagents is assured until the validity date mentioned on the packaged, as long as they are stored and transported under the appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda.

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Tel.: (31) 3439.5454
E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Industria Brasileira

CONSUMER SERVICE

Customer Support Service
Tel.: 0800 0315454
E-mail: sac@bioclin.com.br

Product licensed in the Ministry of Agriculture, Livestock and Supply since 11/11/2019 with the number 10.280/2019.

Technical manager:

Dra. Camila Eckstein (CRMV/MG 20611)

Review: September/2021

UNIVERSAL SYMBOLOGY

CATALOG NUMBER

BATCH CODE

DATE OF MANUFACTURE

USED BY
(last day of month)

TEMPERATURE LIMITATION
(store at)

CONTAINS SUFFICIENT
FOR <N> TESTS

CONSULT INSTRUCTIONS
FOR USE

IN VITRO DIAGNOSTIC DEVICE

EUROPEAN AUTHORIZED
REPRESENTATIVE

KEEP AWAY
FROM SUNLIGHT

MANUFACTURED BY

CONTROL

POSITIVE CONTROL

NEGATIVE CONTROL

BIOLOGICAL RISK

INFLAMMABLE

CORROSIVE

POISON

CE MARK

DO NOT USE IF
PACKAGE IS
DAMAGED