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**Bioclin<sup>+</sup>Vet**

### VETLISA FeLV Ag

REF VET041

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##### FINALIDADE

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- 1- American Animal Hospital Association; American Association of Feline Practitioners. 2020 AAHA/AAFP Feline Vaccination Guidelines. Journal of the American Animal Hospital Association, 56(4), 249–265, 2020.
- 2- American Association of Feline Practitioners. 2020 AAFP Feline Retrovirus Testing and Management Guidelines. Journal of Feline Medicine and Surgery, 22, 5–30, 2020.
- 3- Hofmann-Lehmann, R.; Hartmann, K. Feline Leukaemia Virus Infection – ABCD Recommendations and Review of the Literature. European Advisory Board on Cat Diseases (ABCD), 2021.
- 4- QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

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1- American Animal Hospital Association; American Association of Feline Practitioners. 2020 AAHA/AAFP Feline Vaccination Guidelines. Journal of the American Animal Hospital Association, 56(4), 249–265, 2020.

2- American Association of Feline Practitioners. 2020 AAFP Feline Retrovirus Testing and Management Guidelines. Journal of Feline Medicine and Surgery, 22, 5–30, 2020.

3- Hofmann-Lehmann, R.; Hartmann, K. Feline Leukaemia Virus Infection – ABCD Recommendations and Review of the Literature. European Advisory Board on Cat Diseases (ABCD), 2021.

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### INSTRUCTIONS FOR USE

#### FUNCTION

Test for the qualitative determination of the p27 antigen of the feline leukemia virus in serum or plasma, by microplate enzyme immunoassay. For *in vitro* diagnostic use only.

#### PRINCIPLE OF THE ASSAY

**Methodology:** Enzyme immunoassay or immunoenzymatic assay. The VETLISA FeLV Ag is a solid-phase enzyme immunoassay based on the "sandwich" principle for the qualitative detection of the p27 antigen of the feline leukemia virus (FeLV) in serum or plasma. The microplate is coated with highly specific monoclonal antibodies against the FeLV p27 antigen. Initially, the sample is incubated with the Conjugate, which consists of monoclonal anti-FeLV antigen antibodies conjugated to peroxidase. FeLV antigens present in the sample bind to the monoclonal antibodies specific for the FeLV p27 antigen immobilized on the plate and to the monoclonal anti-FeLV antigen antibodies conjugated to peroxidase, forming an immunocomplex: anti-FeLV antigen – FeLV antigen – anti-FeLV antigen conjugated to peroxidase. After the initial incubation, the microplate is washed to remove unbound materials. Then, the Substrate is added and incubated, producing a blue color that indicates the presence of FeLV p27 antigen in the sample. The Stop Solution is added to stop the reaction, promoting a change of color from blue to yellow, measured in a microplate reader.

#### REAGENTS

- 1 – Sensitized Plate** – Store between 2 and 8°C. Contains: Polystyrene plate, divided into 12 strips of 8 wells each, impregnated with monoclonal antibody against feline leukemia virus (FeLV).
- 2 – Conjugate** – Store between 2 and 8°C. Contains: Monoclonal antibody solution against FeLV antigen, conjugated to Peroxidase.
- 3 – Concentrated Wash (20x)** – Store between 2 and 8°C. Contains: Phosphate Buffer Solution, surfactant and preservative.
- 4 – Substrate TMB** – Store between 2 and 8°C. Contains: Solution containing Tetramethylbenzidine (TMB), Citric Acid Solution and Urea Peroxide.
- 5 – Stop Solution** – Store between 2 and 8°C. Contains: 1N Sulfuric Acid Solution.
- 6 – Negative Control** – Store between 2 and 8°C. Contains: Negative serum for feline viral leukemia and preservative.
- 7 – Positive Control** – Store between 2 and 8°C. Contains: Positive serum for feline viral leukemia and preservative.
- 8 – Plate Sealers**

#### PRESENTATIONS

Components	Presentation		
	1	2	3
	96 cavities	192 cavities	480 cavities
<b>1- Sensitized Plate</b>	1 unit	2 units	5 units
<b>2- Conjugate</b>	1 x 12 mL	2 x 12 mL	5 x 12 mL
<b>3- Concentrated Wash</b>	1 x 50 mL	2 x 50 mL	5 x 50 mL
<b>4- Substrate TMB</b>	1 x 12 mL	2 x 12 mL	5 x 12 mL
<b>5- Stop Solution</b>	1 x 12 mL	2 x 12 mL	5 x 12 mL
<b>6- Negative Control</b>	1 x 1 mL	2 x 1 mL	5 x 1 mL
<b>7- Positive Control</b>	1 x 1 mL	2 x 1 mL	5 x 1 mL
<b>8- Plate Sealers</b>	3 units	6 units	15 units

#### EQUIPMENTS AND OPERATIONAL INPUTS

##### Materials contained in the kit:

- Reagents described in the previous table.
- Instructions for use (manual).

##### Materials needed, but not contained in the kit:

- 1 – Pipettes capable of dispensing volumes from 5 to 500 µL with a coefficient of variation lower than 1.5%.
- 2 – Repipettor for repetitive pipetting of volumes of 500 µL with coefficient of variation less than 1.5% or multichannel pipette (optional).
- 3 – Microplate washer (optional).
- 4 – ELISA reader with absorbance capacity at 450 and 630 nm wavelength.
- 5 – Absorbent paper for drying microwells.
- 6 – Timer or stopwatch.
- 7 – Container for storing the diluted Wash Solution.
- 8 – Distilled or deionized water.
- 9 – Quality Control tools.
- 10 – Incubator at 37°C ± 2°C.

#### STORAGE AND TRANSPORT CONDITIONS

The product storage temperature should be 2 to 8 °C. Transport can be carried out at room temperature (30°C) for up to 5 days. Keep away from light and avoid humidity. **Do not freeze.**

#### SPECIAL CARE

- 1– For professional veterinary *in vitro* diagnostic use only.
- 2– Follow the recommended methodology strictly to obtain accurate results.
- 3– The sachet containing the microplate must be opened only after it reaches room temperature. Replace strips from unused microwells in the sachet, seal and store at 2 to 8 °C.
- 4– The water used for cleaning materials must be fresh and free of contaminants.
- 5– Saturated deionizer columns release alkaline water, different ions, and oxidizing and reducing agents that can significantly alter the results.
- 6– Do not mix reagents from different kits or lots. Do not use expired components.
- 7– Do not eat, drink, or smoke in the testing area.
- 8– Do not pipette reagents or sample(s) by mouth. Do not use the same pipette tip for different samples.
- 9– The Positive and Negative Controls must be retested with each

new assay performed.

- 10– The Stop Solution contains Sulfuric Acid which is a strong acid. Therefore, handle it with due care.
- 11– The professional must strictly follow the safety rules and routines when handling biological samples. The use of disposable gloves and other personal protective equipment is essential.
- 12– Pipette the reagents in the same order each time to minimize reaction time differences between microwells.
- 13– As a safety measure, cover the plate during the reaction.
- 14– Make sure that the bottom of the cavity is clean and dry and that there are no bubbles on the liquid surface before reading the plate. Do not allow the wells to dry out during the test.
- 15– Handle the kit components with due care, in order to avoid contamination. Use new or properly cleaned tips and containers for handling and do not expose reagents, especially Substrate, to strong light or Hypochlorite vapors during storage or incubation steps.
- 16– Products containing serum are potentially capable of transmitting diseases. Appropriate biosafety precautions must be taken when handling these materials.
- 17– We recommend applying local, state and federal environmental protection standards so that reagents and biological material are disposed of in accordance with current legislation.
- 18– For biosafety information or in case of accidents with the product, consult the SDS (Safety Data Sheet) available on the website [www.bioclin.com.br](http://www.bioclin.com.br) or request it via Quibasa's Customer Support Service.
- 19– Do not use the product if the packaging is damaged.
- 20– It is essential that all instruments and equipment used are properly calibrated and undergo regular maintenance.

#### SAMPLES

##### Serum or Plasma (EDTA)

Use serum or plasma. Hemolyzed or highly lipemic samples should not be used. Samples may be stored refrigerated between 2 and 8°C, for a maximum period of 5 days. If the samples cannot be analyzed within 5 days, they can be stored for up to 30 days at –20°C.

#### PROCESS DESCRIPTON

##### Stability After Open

The results of the stability test show that the VETLISA FeLV Ag kit is stable for at least 30 days after opening. This stability may vary depending on test and environmental conditions. Therefore, it is recommended to monitor product performance using the internal kit controls and the assay validation criteria.

#### PREPARATION OF WORK REAGENTS

##### Samples and Controls

Samples should be added directly to the plate; pre-dilution is not required.

**ATTENTION: The Positive and Negative Controls are ready to use.**

##### Washing Solution

Dilute the contents of the Concentrated Wash (Reagent No. 3) in 1000 mL of distilled or deionized water. After preparation, the solution can be stored at 2 to 30 °C up to 30 days. If crystallization occurs, heat to 37°C until dissolved.

##### Conjugate

The Conjugate is ready to use.

##### Substrate

The Substrate is ready to use.

## VETERINARY USE

#### TECHNIQUE

**For use in automatic equipment, consult the Customer Service Department (SAC).**

Before starting the test, place all reagents, samples and controls to stabilize at room temperature (15 - 30 °C) for at least 40 minutes.

- 1– Separate the cavities to be used considering: Controls and Samples (it is recommended to test in duplicate). Return the unused strips of the microplate to the original sealed package.
  - 2– Separate the first cavity for Blank (OPTIONAL).
  - 3– Pipette 50 µL of the Negative Control and Positive Control in the previously determined wells. **Note: Controls are ready to use and do not require dilution.**
  - 4– Pipette 50 µL of the samples in the previously determined wells.
  - 5– Pipette 50 µL of Conjugate into each well, including the Blank well.
  - 6– Gently homogenize for ± 10 seconds. Cover the wells with plate sealers.
  - 7– Incubate for 60 minutes ± 2 minutes in an incubator at 37°C ± 2°C.
  - 8– Remove the plate sealer from the wells.
  - 9– Discard the contents of the wells by aspiration (washer) or decanting (manual).
- Use approximately 300 µL of Washing Solution, previously prepared, and perform a total of five (5) wash cycles. Shake for 3 seconds with each wash. To guarantee the drying of the plate, at the end of the washing, beat the plate for a few seconds on absorbent paper.
- Note:** Poor washing / drying can cause inadequate results.
- 10– Pipette 100 µL of Substrate into all wells.
  - 11– Gently homogenize for ± 10 seconds. Cover the wells with a plate sealer.
  - 12– Incubate for 10 minutes ± 2 minutes in an incubator at 37°C ± 2°C.
  - 13– Remove the plate sealer from the wells.
  - 14– Pipette 100 µL of Stop Solution into all wells.
  - 15– Mix gently for approximately 10 seconds.
  - 16– Read the plate using dual-wavelength filters: 450 nm / 630 nm, within 15 minutes (maximum).

**TECHNIQUE VERIFICATION**

Verify if the results obtained for reading the Blank and the Controls are compatible with the values presented below:

ITEM	ABSORBANCE (DUAL FILTER)
Blank	< 0.100
Negative Control	0.100 to 0.300
Positive Control	> 0.800

If the values are outside the expected values, the technique should be repeated.

**Calculations**

Calculate Cut-off according to the following formula:

$$\text{Cut-off} = \text{Average absorbance of the Negative Control} + 0.200$$

Example:

ITEM	ABSORBANCE
Negative Control	0.152
	0.148
<b>Cut-off = Average absorbance of the Negative Control + 0.200</b>	$((0.152 + 0.148) / 2) + 0.200 = 0.350$

Calculate the Index by dividing the sample absorbance by the Cut-off value. Example:

ITEM	ABSORBANCE
Sample	1.538
Cut-off Value	0.350
<b>Index: Sample Abs. / Cut-off Value</b>	$1.538 / 0.350 = 4.39$

**Note:** The data presented in the examples are for illustration only and cannot be used to calculate the results.

**INTERPRETATION OF THE RESULT**

After calculating the sample indices, consider the indices below to determine the results.

RESULTS	INDEX
Negative	< 0.9
Positive	> 1.1
Undetermined	$\geq 0.9 - \leq 1.1$

**Note:** An indeterminate result indicates that the sample showed a response in the assay's gray zone, between the cut-off limits, preventing confirmation or exclusion of FeLV antigenemia. Intermediate values may occur in very early stages of infection, in cases of very low antigenic load, or due to pre-analytical variations. It is recommended to repeat the test, preferably using a new sample, to rule out analytical variability or interfering factors. In cats with recent exposure risk or compatible clinical suspicion, PCR for FeLV proviral DNA detection may be required as a confirmatory method. When recent exposure cannot be excluded and the animal is clinically stable, repeating p27 testing after approximately 4 to 6 weeks aligns with the antigenemia detection dynamics described in the literature. Until the infectious status is clarified, the cat should be considered potentially infectious to other felines. Final interpretation should be performed by the Veterinarian, integrating history, clinical signs, and complementary tests.

**PROCEDURE LIMITATIONS**

The interpretation of a diagnostic test should not be established on the basis of a single test. Other confirmatory tests must be included before a sample is considered positive. A negative result does not exclude the possibility of exposure. All results must be interpreted in conjunction with other clinical information available before the definitive diagnosis of the disease.

**INTERFERENTS**

No interference was observed for Triglycerides 1500 mg/dL, Ascorbic Acid 2 g/dL, Creatine 200 mg/dL, Bilirubin 1 g/dL, Albumin 2 g/dL, Hemoglobin 1000 mg/dL, Oxalic Acid 60 mg/dL, Factor Rheumatoid 980 IU/mL, C Reactive Protein 41.2 mg/dL and Antistreptolysin O 1023 IU/mL.

**CROSS REACTIVITY**

A cross-reactivity study was conducted, evaluating 88 samples that were negative for Feline Leukemia Virus (FeLV) but positive for other infections. The samples included 10 positive for Feline Infectious Peritonitis (FIP), 7 positive for Feline Herpesvirus, 8 positive for Feline Calicivirus, 8 positive for Feline Panleukopenia (Parvovirus), 9 positive for *Mycoplasma haemofelis* (hemoplasmosis), 7 positive for *Bartonella henselae*, 10 positive for Cryptococcosis, 10 positive for Histoplasmosis, 10 positive for Toxoplasmosis, and 10 positive for Feline Immunodeficiency Virus (FIV). No cross-reactivity was observed with any of the samples. Despite these results, the possibility of cross-reactivity cannot be completely excluded. Final diagnosis should consider the patient's clinical data along with other laboratory and clinical information.

**INTERNAL QUALIT CONTROL**

The Clinical Laboratory must have an internal quality control program, where procedures, standards, limits and tolerance for variations are clearly established. It is important to note that all measurement systems have a characteristic analytical variability, which must be monitored by the laboratories themselves. For that, it is recommended to use controls, which allow to evaluate the precision and accuracy of the measurements.

**PRODUCT PERFORMANCE****ACCURACY****Repeatability**

Repeatability was calculated from 10 successive determinations, using 3 samples with different values. The following results were obtained:

REPEATABILITY	SAMPLE		
	1	2	3
Average	1.550	0.750	0.100
Standard Deviation	0.108	0.044	0.005
Coefficient of Variation (%)	6.954	5.913	5.119

**Reproducibility**

Reproducibility was calculated from 10 successive determinations over 3 consecutive days, using 3 samples with different values. The following results were obtained:

REPRODUCIBILITY	SAMPLE		
	1	2	3
Average	1.500	0.723	0.102
Standard Deviation	0.093	0.043	0.006
Coefficient of Variation (%)	6.191	5.967	5.398

**CLINICAL SENSITIVITY AND SPECIFICITY**

This product was tested in comparison with another immunoassay method. A total of 131 samples were analyzed, including 38 samples diagnosed as positive and 93 diagnosed as negative. The results show that the clinical sensitivity of the product is >99.9% and the clinical specificity is 96.8%.

Method	Reference		Total	
	Result	Positive		Negative
VELLISA FeLV Ag	Positive	90	0	90
	Negative	3	38	41
<b>Total Result</b>		93	38	131

Clinical Sensitivity: >99.9% (38 / 38) – 95% CI: 95%–100%  
 Clinical Specificity: 96.8% (90 / 93) – 95% CI: 93.24%–100%  
 Accuracy: 97.71% [(90 + 38) / (93 + 38)]

**CLINICAL SIGNIFICANCE**

Feline leukemia virus (FeLV) infection is caused by a retrovirus of the genus *Gammaretrovirus*, belonging to the subfamily *Orthoretrovirinae*. Following oronasal infection, the virus replicates in lymphoid tissues of the oropharynx and, in the absence of an effective immune response, reaches the bone marrow, establishing viremia and systemic dissemination. Transmission occurs primarily via saliva, during social or aggressive interactions between cats. The clinical outcome depends on the host's immune response. In abortive infection, the virus is eliminated before viremia occurs; in regressive infection, initial viral replication and transient viremia may occur, followed by immune control that halts active viral particle production, resulting in the persistence of viral DNA; and in progressive infection, persistent viremia and continuous shedding of viral particles are observed. Cats with progressive infection are the main spreaders and have an increased risk of anemia, immunosuppression, secondary infections, and hematopoietic neoplasms such as lymphomas. The p27 antigen, the main structural capsid protein of FeLV, represents the most reliable marker of active infection and is detectable in the blood of viremic cats. Although it may be transiently present in cases of regressive infection, its persistence indicates progressive infection and correlates with viral transmission potential. In this context, VETLISA FeLV Ag stands out as a reliable and accurate diagnostic tool for detecting FeLV p27 antigen in serum, plasma, or whole blood. The test offers high sensitivity, specificity, and reproducibility, providing consistent results that enable early diagnosis and assist Veterinarians in clinical decision-making, monitoring of infected cats, and implementation of effective FeLV control and prevention measures.

**BIBLIOGRAPHIC REFERENCES**

- American Animal Hospital Association; American Association of Feline Practitioners. 2020 AAHA/AAFP Feline Vaccination Guidelines. Journal of the American Animal Hospital Association, 56(4), 249–265, 2020.
- American Association of Feline Practitioners. 2020 AAFP Feline Retrovirus Testing and Management Guidelines. Journal of Feline Medicine and Surgery, 22, 5–30, 2020.
- Hofmann-Lehmann, R.; Hartmann, K. Feline Leukaemia Virus Infection – ABCD Recommendations and Review of the Literature. European Advisory Board on Cat Diseases (ABCD), 2021.
- QUIBASA: Data from the Research and Development Department.

**QUALITY GUARANTEE**

Before being released for consumption, all Quibasa reagents are tested by the Quality Control Department. The quality of the reagents is assured until the validity date mentioned on the packaged, as long as they are stored and transported under the appropriate conditions.

**QUIBASA QUÍMICA BÁSICA Ltda.**

Rua Teles de Menezes, 92 - Santa Branca  
 CEP 31565-130 - Belo Horizonte - MG - Brasil  
 Tel.: +55 31 3439.5454 | E-mail: bioclin@bioclin.com.br  
 CNPJ: 19.400.787/0001-07 - Made in Brazil

**CONSUMER SERVICE**

Customer Support Service  
 Tel.: 0800 0315454 | E-mail: sac@bioclin.com.br

**Product Licensed at the Ministry of Agriculture, Livestock and Supply from 10/12/2019 under the number 10.290/2019.**

**Technical manager:** Dra. Camila Eckstein (CRMV/MG 20611)

**Review:** March/2026

**UNIVERSAL SYMBOLOGY**

	CATALOG NUMBER		MADE BY
	LOT NUMBER		CONTROL
	MANUFACTURING DATE		POSITIVE CONTROL
	VALIDITY DATE (last day of the month)		NEGATIVE CONTROL
	TEMPERATURE LIMIT (store)		BIOLOGICAL RISK
	CONTENT IS SUFFICIENT FOR <N> TEST		FLAMMABLE
	SEE INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC PRODUCT		TOXIC
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED
	DO NOT REUSE		PRODUCT STERILIZED
	CAUTION		DANGER