VETLISA ERLIQUIOSE IgM

REF
VET066

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FINALIDADE

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USAGE INSTRUCTIONS

PURPOSE

Test for qualitative determination of IgM antibodies against *Ehrlichia canis* in serum or plasma of dogs, through immunoenzymatic assay, in microplate. For *in vitro* diagnostics only.

PRINCIPLE OF ACTION

Methodology: Enzyme immunoassay or immunoenzymatic assay.

VETLISA Ehrlichiosis IgM is a solid phase immunoenzymatic assay, based on the immunocapture principle for the qualitative detection of IgM antibodies against *Ehrlichia canis* in serum or plasma of dogs. Antibodies against Ehrlichia canis present in the sample bind to the recombinant antigens of *Ehrlichia canis* coated on the microplate, forming immobilized antigen-antibody immunocomplexes. After the initial incubation, the microplate is washed to remove unbound materials. The Conjugate, formed by anti-dog IgM antibody linked to Peroxidase, is added and binds to IgM antibodies linked to the antigens immobilized on the plate. After a second incubation, the microplate is washed to remove unbound materials. Then, the Substrate is added and incubated, producing a blue color, whose intensity is proportional to the amount of IgM antibodies against *Ehrlichia canis* present in the sample. The Stop Solution is added to stop the reaction, promoting a change of color from blue to yellow, measured in a microplate reader.

REAGENTS

1- Sensitized Plate - Store between 2 and 8°C. Contains: Polystyrene plate, divided into 12 strips of 8 wells each, impregnated with *Ehrlichia canis* recombinant antigen.

2- Concentrated Conjugate (100X) - Store at 2 to 8°C. Contains: Dog anti-IgM antibody conjugated to peroxidase and protein stabilizer.

3- Concentrated Wash (20X) - Store at 2 to 8°C. Contains: Buffer solution, surfactant and preservative.

4- Diluent - Store at 2 to 8°C. Contains: Buffer solution, stabilizers, dye, surfactant and preservative.

5- Substrate TMB - Store between 2 and 8°C. Contains: Solution containing 3, 3', 5, 5' - Tetramethylbenzidine (TMB), Urea Peroxide and preservative.

6- Stop Solution - Store at 2 to 8°C. Contains: Sulfuric Acid Solution 1 mol / L.

7- Negative Control - Store between 2 and 8°C. Contains: Serum containing negative serum for IgM antibodies against *Ehrlichia canis* and protein stabilizer. **Potentially infective.**

8- Positive Control - Store between 2 and 8°C. Contains serum containing positive serum for IgM antibodies against *Ehrlichia canis* and protein stabilizer. **Potentially infective.**

9- Plate Sealers

PRESENTATIONS

COMPONENTS	PRESENTATION			
	1	2	3	4
	96 wells	192 wells	480 wells	960 wells
1- Sensitized Plate	1 unit	2 units	5 units	10 units
2- Concentrated Conjugate (100x)	1 x 300 µL	2 x 300 µL	5 x 300 µL	10 x 300 µL
3- Concentrated Wash (20x)	1 x 50 mL	2 x 50 mL	5 x 50 mL	10 x 50 mL
4- Diluent	1 x 60 mL	2 x 60 mL	5 x 60 mL	10 x 60 mL
5- Substrate TMB	1 x 15 mL	2 x 15 mL	5 x 15 mL	10 x 15 mL
6- Stop Solution	1 x 15 mL	2 x 15 mL	5 x 15 mL	10 x 15 mL
7- Negative Control	1 x 1 mL	2 x 1 mL	5 x 1 mL	10 x 1 mL
8- Positive Control	1 x 1 mL	2 x 1 mL	5 x 1 mL	10 x 1 mL
9- Plate Sealers	3 units	6 units	15 units	30 units

EQUIPMENT AND OPERATIONAL INPUTS

Materials contained in the kit:

- Reagents described in the previous table

- Instructions for use (manual)

Necessary materials not contained in the kit:

1- Pipettes capable of dispensing volumes of 5 to 500 µL with a variation coefficient of less than 1.5%.

2- Repipettor for repetitive pipetting of volumes of 500 µL with coefficient of variation less than 1.5% or multichannel pipette (optional).

3- Microplate washer (optional).

4- ELISA reader with absorbance capacity at 450 and 630 nm wavelength.

5- Absorbent paper to dry the wells.

6- Stopwatch or watch.

7- Bottle to store the Washing Solution after dilution.

8- Distilled or deionized water.

9- Quality Control Tools.

10- Incubator at 37 ° C ± 2 ° C.

STORAGE ANS TRANSPORT CONDITIONS

The product storage temperature is 2 to 8 °C. Transport can be done at room temperature (up to 30 ° C) for up to 5 days. Keep away from light and avoid humidity. **Do not freeze.**

SPECIAL CARE

1- For professional veterinary *in vitro* diagnostic use only.

2- Strictly follow the proposed methodology to obtain accurate results.

3- The sachet containing the microplate must be opened only after reaching room temperature. Replace strips from unused microwells in the sachet, seal and store at 2 to 8 °C.

4- The water used to clean the material must be fresh and free of contaminants.

5- Saturated deionizer columns release alkaline water, different ions, and oxidizing and reducing agents that can significantly alter the results.

6- The Stop Solution contains Sulfuric Acid which is a strong acid. So, handle it with due care.

7- The handling of any product that contains serum is potentially capable of transmitting diseases. Therefore, it is necessary to take proper biosafety precautions when handling these products.

8- Pipette the reagents always in the same order to minimize the difference in reaction time between the microwells.

9- For protection measure, the plate must be covered during the reaction.

10- Ensure that the bottom of the cavity is clean and dry and that there are no bubbles on the surface of the liquid before reading the plate. Do not allow wells to dry during the assay.

11- Do not expose reagents, especially Substrate, to strong light or Hypochlorite vapors during storage or incubation steps.

12- We recommend applying local, state and federal regulations for environmental protection so that the disposal of reagents and biological material is done in accordance with current legislation.

13- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Safety Information Sheet for Chemical Products) available on the website www.bioclin.com.br or upon request by the SAC (Service of Customer Service) by Quibasa.

14- Do not use the product in case of damage to the packaging.

15- It is essential that the instruments and equipment used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Serum or Plasma (EDTA or Heparin).

Hemolyzed or highly lipemic samples should not be used.

The samples can be kept refrigerated, between 2 and 8 °C, for a maximum period of 5 days. If the samples cannot be analyzed within 5 days, they can be stored for up to 30 days at -20 °C.

PROCESS DESCRIPTON

Stability After Open

The results of the stability test show that the VETLISA Ehrlichiosis IgM kit is stable after being opened up to 30 days. This stability may vary according to the test conditions and the environment. Therefore, it is suggested to monitor the product's performance using internal kit controls and the technique validation criteria.

ATTENTION: The Positive and Negative Controls are ready to use.

PREPARATION OF WORK REAGENTS

Washing Solution

Dilute the contents of the Concentrated Wash (Reagent N° 3) in 1000 mL of distilled or deionized water. After preparation, the solution can be stored at 2 to 30 °C up to 30 days. If crystallization occurs, heat to 37°C until dissolved.

Conjugate Solution

Dilute the Concentrated Conjugate (Reagent N° 2) in the proportion of 1: 101 in Diluent (Reagent N° 4). Prepare the solution when performing the test.

To perform a test using all the wells of the kit, mix 110 µL of the Concentrated Conjugate in 11 mL of Diluent.

To perform an assay using 8 wells (1 strip), mix 10 µL of the Concentrated Conjugate in 1 mL of Diluent.

IMPORTANT: The diluted conjugate solution cannot be stored. Therefore, prepare only the amount necessary to perform the test.

Dilution of Samples

In a test tube, dilute 5 µL of the sample in 500 µL of Diluent (Reagent N° 4), if the test is going to be performed in duplicate. Cap the tube and vortex gently or mix manually by inversion. Dilutions cannot be stored.

SUBSTRATE

The Substrate is ready to use.

TECHNIQUE

For use in automatic equipment, consult the Customer Service Department (SAC).

Before starting the test, place all reagents, samples and controls to stabilize at room temperature (15 - 30 °C) for at least 40 minutes.

1- Separate the cavities to be used considering: Controls and Samples (it is recommended to test in duplicate). Return the unused strips of the microplate to the original sealed package.

2- Separate the first cavity for Blank (OPTIONAL).

3- Pipette 100 µL of the Negative Control and Positive Control in the previously determined wells. **Note: The controls are ready for use, there is no need to dilute them.**

4- Pipette 100 µL of the samples previously diluted in the previously determined wells. In the Blank cavity (OPTIONAL), if you have made the option, pipette only 100 µL of the diluent.

5- Gently mix for ± 10 seconds. Cover the wells with plate sealer.

6- Incubate for 30 minutes ± 2 minutes in an incubator at 37°C ± 2°C.

7- Remove the sealer from the cavities.

8- Discard the contents of the cavities by aspiration (Washer) or by decantation (Manual). Use approximately 300 µL of Washing Solution, **previously prepared**, and perform a total of five (5) wash cycles. Shake for 3 seconds with each wash. To guarantee the drying of the plate, at the end of the washing, beat the plate for a few seconds on absorbent paper.

IMPORTANT: Poor washing / drying can cause inadequate results.

9- Pipette 100 µL of Conjugate previously diluted in all wells.

10- Gently homogenize for ± 10 seconds. Cover the cavities with the plate sealer.

11- Incubate for 30 minutes ± 2 minutes in an incubator at 37°C ± 2°C.

12- Remove the plate sealer from the cavities.

13- Repeat item 8.

14- Pipette 100 µL of Substrate into all wells.

15- Gently homogenize for ± 10 seconds. Cover the cavities with the plate sealer.

16- Incubate for 10 minutes ± 2 minutes in an incubator at 37°C ± 2°C.

17- Remove the plate sealer from the cavities.

18- Pipette 100 µL of Stop Solution into all wells.

19- Gently mix for ± 10 seconds.

20- Read using a double filter: 450 nm / 630 nm in up to 10 minutes (maximum).

TECHNIQUE VERIFICATION

Verify if the results obtained from reading the Blank and the Controls are compatible with the specification above:

ITEM	ABSORBANCE (DUAL FILTER)
Blank	< 0.100
Negative Control	0.050 a 0.250
Positive Control	> 0.800

If the values are out the expected values, the technique should be repeated.

CALCULATIONS

Calculate the Cut Off according to the following formula:
Cut Off = Average Negative Control Absorbance + 0.500
Example:

ITEM	ABSORBANCIA
	A1 = 0.153
Negative Control	A2 = 0.157
Cut Off = Negative Control Average Absorbance + 0.500	(0.153 + 0.157) / 2 + 0.500 = 0.655

Calculate the index by dividing the absorbance of the sample by the cut off value.
Example:

ITEM	ABSORBANCE
Sample	1.680
Cut Off Value	0.655
Index: Sample Abs. / Cut Off Value	1.680 / 0.655 = 2.34

Note: The data presented in the examples are for illustration only and cannot be used to calculate the results.

INTERPRETATION OF THE RESULT

After calculating the sample index, consider the indices below to determine the results.

RESULTS	INDEX
Negative	< 0.9
Positive	> 1.1
Undetermined	0.9 y 1.1

Note: In the case of an undetermined result, the sample must be re-analyzed. Samples that obtain repeatedly indeterminate results must be retested using an alternative method. If the results remain indeterminate, a new sample should be collected in two weeks. The result of the last sample collected must prevail. The results provided by this kit must be interpreted by the responsible medical professional, and it is not the only criterion for determining the diagnosis and / or treatment of the patient.

PROCESS LIMITATIONS

The interpretation of a diagnostic test should not be established on the basis of a single test. Other confirmatory tests must be included before a sample is considered positive. A negative result does not exclude the possibility of exposure. All results must be interpreted in conjunction with other clinical information available before the definitive diagnosis of the disease.

INTERFERENTS

No interference was observed for Triglycerides 1500 mg/dL, Ascorbic Acid 2 g/dL, Creatine 200 mg/dL, Bilirubin 1 g/dL, Albumin 2 g/dL, Hemoglobin 1000 mg/dL, Oxalic Acid 60 mg/dL, Factor Rheumatoid 980 IU/mL. C Reactive Protein 41.2 mg/dL and Anti-streptolysin O 1023 IU/mL.

CROSS REACTIVITY

A cross-reactivity study was carried out, evaluating 35 samples negative for IgM antibodies against *Ehrlichia canis*, but positive for other infections. 6 positive samples for Brucellosis, 5 positive samples for Distemper, 1 positive sample for Parvovirus, 12 positive samples for Babesiosis, 4 positive samples for Heartworm, 4 positive samples for Leptospirosis and 3 positive samples for Leishmaniasis were tested. No cross reactivity was observed with any of the samples. Despite the results found, the possibility of cross-reactivity cannot be completely ruled out. The final diagnosis should consider the patient's clinical data along with other laboratory data.

INTERNAL QUALIT CONTROL

The Clinical Laboratory must have an internal quality control program, where procedures, standards, limits and tolerance for variations are clearly established. It is important to note that all measurement systems have a characteristic analytical variability, which must be monitored by the laboratories themselves. For that, it is recommended to use controls, which allow to evaluate the precision and accuracy of the measurements.

PRODUCT PERFORMANCE

ACCURACY

Repeatability

Repeatability was calculated from 10 successive determinations, using 3 samples with different values. The following results were obtained:

REPEATABILITY	ABSORBANCE		
	1	2	3
Mean	1.492	0.775	0.227
Standard Deviation	0.077	0.012	0.007
Coefficient of Variation (%)	5.18	1.51	3.28

Reproducibility

Reproducibility was calculated from 10 successive determinations over 3 consecutive days, using 3 samples with different values. The following results were obtained:

REPRODUCIBILITY	ABSORBANCE		
	1	2	3
Mean	1.473	0.793	0.197
Standard Deviation	0.077	0.020	0.032
Coefficient of Variation (%)	5.26	2.46	16.31

CLINICAL SENSITIVITY AND SPECIFICITY

In the VETLISA Ehrlichiosis IgM, clinical samples were analyzed in comparison with another method of Enzyme immunoassay. The results show that the clinical sensitivity of the kit is 97.8% and the clinical specificity is 96.4%.

VETLISA Ehrlichiosis IgM	Method	Reference		Total
	Result	Positive	Negative	
	Positive	44	2	46
	Negative	1	53	54
Total Result		45	55	100

Clinical Sensitivity: 97.8% (44/45) - IC 95%: 91.5% - 100%

Clinical Specificity: 96.4% (53/55) - IC 95%: 93.5% - 100%

Precision: 97.0% (44+53) / (55+45)

CLINICAL SIGNIFICANCE

Ehrlichiosis is a disease caused by bacteria of the genus *Ehrlichia*, which parasitize blood cells such as leukocytes, erythrocytes, endothelial cells and platelets. Transmission occurs through the bite of the brown tick (*Rhipicephalus sanguineus*), which is contaminated with the bacteria after the blood meal in dogs with high parasitemia. The incubation period for ehrlichiosis in the dog varies from 8 to 20 days, and the clinical signs observed include petechiae on the skin and mucous membranes, enlarged lymph nodes, hepatomegaly and splenomegaly, apathy, weight loss, lack of appetite, anorexia, fever, runny nose and eyepiece. When treatment does not result in the total elimination of the organism, subsequent relapses may occur. The dog's immune response against *Ehrlichia* is essential to fight infection and also allows diagnosis by detecting antibodies in serological tests. Immunoglobulin M (IgM), is the first class of antibodies produced by the immune system, and its production indicates the occurrence of recent infection. Thus, the detection of IgM favors early diagnosis and favors early treatment.

BIBLIOGRAPHIC REFERENCES

- 1- COUTO CG. Doenças Rickettsiais In: BIRCHARD, SHERDING, Manual Saunders: Clínica de pequenos animais. Ed. Roca: 139-42, 1998.
- 2- ETTINGER SJ, FELDMAN EC. Tratado de Medicina Interna Veterinária. 4 ed. São Paulo: Editora Manole. 1992.
- 3- ROSEZ KV, ALVES FR, BLEICH I. Eriiquiose canina. Revista Cães & Gatos, n.96, p. 25-28. 2001.
- 4- WOLDEHIWET Z, RISTIC M. Rickettsial and chlamydial diseases of domestic animals, 1ed, Oxford, Pergamon Press, 1993, 427 p.
- 5- WOODY BJ, HOSKINS JD. Erlichial Diseases of dogs. Veterinary Clinics of North America:Small Animal Practice. n.1. v. 21, p.129-140. Janeiro1991.
- 6- QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

QUALITY GUARANTEE

Before being released for consumption, all Quibasa reagents are tested by the Quality Control Department. The quality of the reagents is assured until the validity date mentioned on the packaged, as long as they are stored and transported under the appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda.

Rua Teles de Menezes, 92 - Santa Branca

CEP 31565-130 - Belo Horizonte - MG - Brasil

Tel.: +55 31 3439.5454 | E-mail: bioclin@bioclin.com.br

CNPJ: 19.400.787/0001-07 - Industria Brasileira

CONSUMER SERVICE

Customer Support Service

Tel.: 0800 0315454 | E-mail: sac@bioclin.com.br

Product licensed in the Ministry of Agriculture, Livestock and Supply since 04/09/2022 with the number 10.531/2022.

Technical manager: Dra. Camila Eckstein (CRMV/MG 20611)

Review: October/2022

UNIVERSAL SYMBOLOGY

CATALOG NUMBER

LOT NUMBER

MANUFACTURING DATE

VALIDITY DATE
(last day of the month)

TEMPERATURE LIMIT
(store)

CONTENT IS SUFFICIENT
FOR <N> TEST

SEE INSTRUCTIONS
FOR USE

IN VITRO DIAGNOSTIC
PRODUCT

KEEP AWAY
FROM SUNLIGHT

DO NOT REUSE

CAUTION

MADE BY

CONTROL

POSITIVE CONTROL

NEGATIVE CONTROL

BIOLOGICAL RISK

FLAMMABLE

CORROSIVE

TOXIC

DO NOT USE IF
PACKAGE IS
DAMAGED

PRODUCT
STERILIZED

DANGER