VETLISA LEISHMANIOSE IgG

REF VET043

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FINALIDADE

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PELIGRO

FUNCTION

Test for qualitative determination of antibodies against *Leishmania infantum* using canine serum or plasma, by enzyme immunoassay, in microplate. For *in vitro* diagnostic use only.

PRINCIPLE OF ACTION

Methodology: Enzyme immunoassay or immunoenzymatic

This product is an immunoenzymatic assay in solid phase, based on the principle of immunocapture for the qualitative detection of IgG antibodies against *Leishmania* in dog serum or plasma. Antibodies against *Leishmania* present in the sample, bind to the recombinant k39 protein of *L. infantum* coated on the microplate, forming immobilized antigen-antibody complexes. After the initial incubation, the microplate is washed to eliminate unbound materials. The Conjugate, formed by the anti-dog IgG antibody conjugated to peroxidase, binds the immobilized antigen-antibody complex on the plate. Then, the microplate is washed and incubated with Substrate. The intensity of the blue color produced by the addition of the Substrate is proportional to the amount of antibodies against *Leishmania* present in the sample. The Stop Solution is added to stop the reaction, promoting a color change to yellow, measured in a microplate reader.

REAGENTS

1- Sensitized Plate - Store between 2 and 8°C. Contains: Polystyrene plate, divided into 12 strips of 8 cavities each, impregnated with recombinant *L. infantum* antigen.

2- Concentrated Conjugate (100x) - Store between 2 and 8°C. Contains: Peroxidase-linked anti-dog IgG antibody solution.

3- Concentrated Wash Solution (20x) - Store between 2 and 8°C. Contains: Buffer Solution, surfactant and preservative.

4- Diluent - Store between 2 and 8°C. Contains: Buffer solution, surfactant, stabilizer and preservative.

5- Substrate TMB - Store between 2 and 8°C. Contains: Solution containing Tetramethylbenzidine (TMB <1.0mg/mL), Citric Acid Solution <5% and Urea Peroxide <1%.

6- Stop Solution - Store between 2 and 8°C. Contains: 1M Sulfuric Acid.

7- Negative Control - Store between 2 and 8°C. Contains: Non-reactive serum for leishmaniasis and preservative. **Potentially infectious.**

8- Positive Control - Store between 2 and 8°C. Contains: Reactive serum for leishmaniasis and preservative. **Potentially infectious.**

9- Plate Sealers

REAGENTS	Presentation		
	1	2	3
1- Sensitized Plate	1 unit (96 cavities)	2 units (192 cavities)	5 units (480 cavities)
2- Concentrated Conjugate (100x)	1 x 350 µL	2 x 350 µL	5 x 350 µL
3- Concentrated Wash Solution (20x)	1 x 50 mL	2 x 50 mL	5 x 50 mL
4- Diluent	1 x 60 mL	2 x 60 mL	5 x 60 mL
5- Substrate TMB	1 x 15 mL	2 x 15 mL	5 x 15 mL
6- Stop Solution	1 x 15 mL	2 x 15 mL	5 x 15 mL
7- Negative Control	1 x 1.0 mL	2 x 1.0 mL	5 x 1.0 mL
8- Positive Control	1 x 1.0 mL	2 x 1.0 mL	5 x 1.0 mL
9- Plate Sealers	3 unit	6 units	15 units

EQUIPMENT AND OPERATIONAL INPUTS

Materials contained in the kit:

- Reagents described in the previous table.
- Instructions for use (manual).

Materials needed, but not contained in the kit:

- 1- Set of calibrated pipettes, capable of dispensing volumes between 5 and 1000 µL, with a coefficient of variation of less than 1.5%.
- 2- Repeater or multichannel pipette calibrated at 100 µL volume, with a coefficient of variation of less than 1.5% (optional).
- 3- Microplate washing machine (optional).
- 4- ELISA reader with absorbance capacity at 450 and 630 nm wavelength.
- 5- Absorbent paper to dry the microcavities.
- 6- Stopwatch or clock.
- 7- Bottle to store the Washing Solution, after dilution.
- 8- Distilled or deionized water.
- 9- Quality Control Tools.

STORAGE AND TRANSPORT CONDITIONS

The product storage temperature should be 2 to 8 °C. Transport can be carried out at temperatures up to 30°C for up to 5 days. Keep away from light and avoid humidity. **Do not freeze.**

SPECIAL CARE

- 1- **Veterinary use only, for *in vitro* diagnostic.**
- 2- Strictly follow the methodology proposed to obtain accurate results.
- 3- The professional must strictly follow safety rules and routines when handling biological samples. The use of disposable gloves and other personal protective equipment is essential.
- 4- **IMPORTANT:** Before starting the test, allow all kit components to reach room temperature. Open the envelope containing the cavities only after reaching room temperature.
- 5- Do not mix reagents from different kits or lots. Do not use expired kit components.

- 6- Do not eat, drink, or smoke at the test site.
- 7- Do not pipette reagents or samples by mouth. Do not use the same tip to pipette different samples.
- 8- Positive and negative controls should be retested for each new test performed.
- 9- The Stop Solution contains Sulfuric Acid, which is a strong acid. Therefore, handle it with care.
- 10- Handle the kit components with care to avoid contamination. Take special care with the substrate, which is a light sensitive solution. Use new or appropriately clean tips and containers for handling and do not allow them to be exposed to strong light during storage or during the incubation periods of the assay.
- 11- Use contaminant-free freshly distilled or deionized water when preparing solutions.
- 12- Always add the reagents in the same order to minimize the difference in reaction time between the microcavities.
- 13- Before reading the plate, ensure that the bottom of the microcavities are clean and dry. Make sure there are no bubbles on the surface of the liquid.
- 14- We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.
- 15- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Customer Advisory Service) of Quibasa
- 16- Do not use the product in case of damaged packaging
- 17- It is essential that the instruments and equipment used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Hemolysis, lipemia and bilirubin (jaundice) can interfere with the test result. The result obtained from samples under these conditions must be evaluated with caution.

Samples may be refrigerated at between 2 and 8°C for a maximum of 5 days.

If samples cannot be analyzed within 5 days, they can be stored for up to 30 days at -20°C (freezer)

Use serum or plasma (EDTA or Heparin).

PROCESS DESCRIPTION

Stability After Open

The results of the stability test show that the VETLISA LEISHMANIASIS IgG kit is stable after being opened for up to 30 days. This stability may vary according to the test conditions and the environment. Therefore, it is suggested to monitor the product's performance using internal kit controls and the technique validation criteria.

ATTENTION: The Positive and Negative Controls are ready to use.

PREPARATION OF WORK REAGENTS AND SAMPLES

1- Washing Solution

Dilute the content of the Concentrated Wash Solution (Reagent N° 3) in 1000 mL of freshly distilled or deionized water. Store at 2 to 8°C until the validity date printed on the original bottle. If crystallization occurs, heat to 37°C until dissolution. It can be stored at room temperature.

2- Conjugate Solution

Dilute the Concentrated Conjugate (Reagent N° 2) in the ratio of 1:101 in Diluent (Reagent N° 4). Prepare the solution at the moment of performing the test.

To perform an assay using all cavities in the kit, mix 110 µL of Concentrated Conjugate in 11 mL of Diluent.

To perform an assay using 8 cavities (1 strip), mix 10 µL of the Concentrated Conjugate in 1 mL of Diluent.

IMPORTANT: The diluted conjugate solution cannot be stored. Therefore, prepare only the amount necessary to perform the test.

3- Sample Dilution

In a test tube, dilute 15 µL of the sample in 300 µL of Diluent (Reagent N° 4), if the test is going to be performed in duplicate. Cap the tube and gently vortex or manually mix by inversion. Dilutions cannot be stored.

TECHNIQUE

Before starting the assay, place all reagents, Controls and Samples to stabilize at room temperature (15 - 30°C) for at least 30 minutes. Return the unused strips to the original sealed packaging.

- 1- Separate the cavities to be used considering: Controls and Samples ((it is recommended to test in duplicate).
- 2- Separate the first cavity for the Blank (OPTIONAL).
- 3- Pipette 100 µL of the Negative Control and the Positive Control into the previously determined cavities. **Note: The controls are ready to use, there is no need to dilute them.**
- 4- Pipette 100 µL of the previously diluted samples into the previously determined cavities. In the Blank cavity (OPTIONAL), if you have chosen, pipette only 100 µL of the diluent.
- 5- Homogenize gently for ± 10 seconds. Cover the cavities with plate sealer.
- 6- Incubate for 30 minutes at room temperature.
- 7- Discard the content of the cavities by aspiration (washing machine) or by decantation (manual). Use approximately 300 µL/ cavity of previously diluted Wash Solution, for a total of three (3) wash cycles. Shake for three seconds after each wash cycle. To guarantee the drying of the plate, at the end of the wash, tap the plate for a few seconds on absorbent paper.

IMPORTANT: Poor washing/drying can cause inadequate results.

- 8- Pipette 100 µL of the Conjugate previously diluted into each cavity, even in the Blank cavity.

- 9- Homogenize gently for ± 10 seconds. Cover the cavities with plate sealer.
- 10- Incubate for 30 minutes at room temperature.
- 11- Repeat item 7.
- 12- Add 100 µL of Substrate TMB into each cavitie.
- 13- Homogenize gently for ± 3 seconds. Cover the cavities with plate sealer.
- 14- Incubate for 10 minutes at room temperature, protected from light.
- 15- Remove the plate sealer from the cavities.
- 16- Pipette 100 µL of Stop Solution into all cavities.
- 17- Homogenize gently for ± 3 seconds
- 18- Read the absorbance in an ELISA reader on a dual filter 450 nm (primary filter) / 630nm (secondary filter) within a maximum of 10 minutes after adding the Stop Solution.

TECHNIQUE VERIFICATION

Verify if the results obtained from reading the Controls and the Blank are compatible with the specification above:

ITEM	ABSORBANCE (DUAL FILTER)
Blank	< 0.050
Negative Control	< 0.450
Positive Control	> 1.000

If the values are out the expected ones, the test should be repeated.

CUT OFF CALCULATION

If the results of the controls are valid, calculate the cut off value with the following formula:

Cut Off = Negative Control Average Absorbance + 0.500

CALCULATION OF THE SAMPLE INDEX

Calculate the index by dividing the absorbance of the sample by the cut off value. Example:

ITEM	ABSORBANCE
Sample	1.093
Negative Control	0.201
Cut Off Value	0.201 + 0.500 = 0.701
Index: Sample Abs. / Cut Off Value	1.093 / 0.701 = 1.559

INTERPRETATION OF THE RESULT - QUALITATIVE

RESULTS	INDEX
Negative	< 0.9
Positive	> 1.1
Undetermined	≥ 0.9 e ≤ 1.1

IMPORTANT: Samples with indeterminate results must be retested. Samples that repeatedly show indeterminate results should be analyzed using an alternative method. If this sample shows a positive result, then it should be considered positive.

PROCEDURE LIMITATIONS

The interpretation of a diagnostic test should not be based on a single trial. Other confirmatory tests must be included before a sample is considered positive. A negative result does not exclude the possibility of exposure. However, all results should be interpreted in conjunction with vaccination history, available clinical and laboratory information.

INTERNAL QUALIT CONTROL

The Clinical Laboratory must have an internal quality control, where all procedures, rules, limits and tolerance to variations be clearly established. It is important to mention that all measurement systems present a analytical variety, and it must be monitor by the laboratory. Therefore, it is recommendable the use of controls, allowing the precision and accuracy of the dosages.

PRODUCT PERFORMANCE

PRECISION

Repeatability

Repeatability was calculated from 10 successive determinations, using 3 samples with different values, obtaining the following results:

Repeatability	Sample		
	1	2	3
Mean	1.794	1.029	0.040
Standard Deviation	0.053	0.035	0.003
Coefficient of Variation (%)	2.977	3.370	7.764

Reproducibility

Reproducibility was calculated from 10 successive determinations over 3 consecutive days, using 3 samples with different values, obtaining the following results:

Reproducibility	Sample		
	1	2	3
Mean	1.807	1.026	0.040
Standard Deviation	0.050	0.035	0.003
Coefficient of Variation (%)	2.760	3.422	6.895

CLINICAL SENSITIVITY AND SPECIFICITY

The VETLISA LEISHMANIASIS IgG analyzed clinical samples compared to other enzyme immunoassay methods. The results show that the clinical sensitivity of the product is > 99.9% and the clinical specificity is 96%.

METHOD		REFERENCE		Total
VETLISA LEISHMANIASIS IgG	Result	Positive	Negative	
	Positive	144	0	144
	Negative	6	150	156
Total Result		150	150	300

Clinical sensitivity: >99.9% (150/150) - IC 95%: 97.47 - 100%
Clinical specificity: 96% (144/150) - IC 95%: 91.82 - 98.58%
Precision: 98% [(144 + 150) / (150 + 150)]

CROSS REACTION

This product was tested with serum samples from animals diagnosed as negative for canine leishmaniasis and positive for the following diseases: ehrlichiosis, babesiosis, leptospirosis, anaplasmosis, distemper, parvovirus, and *Trypanosoma cruzi* infection. The product showed a negative result under all tested conditions.

REACTION WITH VACCINATED DOGS

This product was tested with serum samples from dogs vaccinated for Leishmaniasis, with a vaccine available on the market. The product was negative for all tested samples.

CLINICAL SIGNIFICANCE

Leishmaniasis is an infectious and zoonotic disease caused by protozoa of the genus *Leishmania*, and transmitted by sandfly bites. The incubation period of Leishmaniasis can range from a few weeks to several years, and the main clinical manifestations are generalized lymphadenopathy, exfoliative dermatitis, progressive weight loss and keratoconjunctivitis. Due to the high frequency of asymptomatic cases, the diagnosis serology is an effective tool and indicated by the Ministry of Health. According to official bodies, the diagnosis of Leishmaniasis should consider screening the dog using immunochromatographic assays, followed by the confirmatory immunoenzymatic assay.

BIBLIOGRAPHIC REFERENCES

- 1- QUIBASA: Data from the Research and Development Department.
- 2- KILLICK-KENDRICK R, KILLICK-KENDRICK M, PINELLI E, DEL REAL G, MOLINA R, VITUTIA MM, CANAVATE MC, NIETO J. A laboratory model of canine leishmaniasis: the inoculation of dogs with Leishmania infantun promastigotes from midguts of experimentally infected phlebotomine sandfl ies. Parasite, 1994, 7, 311-318.
- 3- RIBEIRO VM. Leishmanioses. In: Associação Nacional de Clínicos Veterinários de Pequenos Animais;
- 4- DENARDI AB, ROZA MR, organizadores. PROMOVET Pequenos Animais: Programa de Atualização em Medicina Veterinária: Ciclo 1. Porto Alegre: Artmed Panamericana; 2016.p.107- 50.
- 5- Ministério da Saúde. Manual de Vigilância e Controle da Leishmaniose Visceral. 1ª ed. Brasília- DF, 2006.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Quality Control Department. The quality of the reagents is assured until the validity date mentioned on the packaged, as long as they are stored and transported under the appropriate conditions.

 **QUIBASA QUÍMICA BÁSICA Ltda.**
Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Tel.: +55 31 3439.5454 | E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

CONSUMER SERVICE

Customer Support Service
Tel.: 0800 0315454 | E-mail: sac@bioclin.com.br

Product licensed in the Ministry of Agriculture, Livestock and Supply since 02/08/2019 with the number 10.264/2019.

Technical manager: Dra. Camila Eckstein (CRMV/MG 20611)

Review: April/2023

CATALOG NUMBER

LOT NUMBER

MANUFACTURING DATE

VALIDITY DATE
(last day of the month)

TEMPERATURE LIMIT
(store)

CONTENT IS SUFFICIENT
FOR <N> TEST

SEE INSTRUCTIONS
FOR USE

IN VITRO DIAGNOSTIC
PRODUCT

KEEP AWAY
FROM SUNLIGHT

DO NOT REUSE

CAUTION

MADE BY

CONTROL

POSITIVE CONTROL

NEGATIVE CONTROL

BIOLOGICAL RISK

FLAMMABLE

CORROSIVE

TOXIC

DO NOT USE IF
PACKAGE IS
DAMAGED

PRODUCT
STERILIZED

DANGER