Bioclin ^{+Vet}

VETLISA FIV Ac

REF **VET042**

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FINALIDADE

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- Hosie MJ, Addie DD, Belák S, et al. ABCD guidelines on feline immunodeficiency virus (FIV) infection. *Journal of Feline Medicine and Surgery*. 2009;11(7):575–584. Atualizado em 2017.
- Little S, Levy J, Hartmann K, et al. 2020 AAFP Feline Retrovirus Testing and Management Guidelines. *Journal of Feline Medicine and Surgery*. 2020;22(1):5–30.
- Spriffler F, Jongwattanapisan P, Luengyosuechakul S, et al. Prevalence and Risk Factors of Feline Immunodeficiency Virus and Feline Leukemia Virus Infection in Healthy Cats in Thailand. *Frontiers in Veterinary Science*. 2022;8:764217.
- Westman ME, Coggins SJ, van Dorselaer M, et al. Feline immunodeficiency virus (FIV) infection in domestic pet cats in Australia and New Zealand: Guidelines for diagnosis, prevention and management. *Australian Veterinary Journal*. 2022;100(8):345–359.
- Westman ME, Malik R, Norris JM. Diagnosing feline immunodeficiency virus (FIV) and feline leukaemia virus (FeLV) infection: an update for clinicians. *Australian Veterinary Journal*. 2019;97(3):47–55.
- QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

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FUNCTION

Test for the qualitative determination of total antibodies against feline immunodeficiency virus (FIV) in feline serum or plasma samples, using a microplate enzyme immunoassay. For *in vitro* diagnostic use only.

PRINCIPLE OF ACTION

Methodology: Enzyme immunoassay or immunoenzymatic assay. The VETLISA FIV Ac is a solid-phase enzyme immunoassay based on the principle of immunocapture for the qualitative detection of total antibodies against feline immunodeficiency virus (FIV) in serum or plasma. Antibodies against feline immunodeficiency virus present in the sample bind to the recombinant viral antigens that coat the microplate, forming immobilized complexes: FIV antigen–anti-FIV antibodies. After the initial incubation, the microplate is washed to remove unbound materials. The Conjugate, composed of anti-feline antibodies conjugated to peroxidase, is added and binds to the immobilized antigen–antibody complex on the plate. After the second incubation, the microplate is washed to remove excess material. Following this step, the Substrate is added and incubated, producing a blue color that indicates the presence of antibodies against FIV. The Stop Solution is added to interrupt the reaction, promoting a color change from blue to yellow, which is measured in a microplate reader.

REAGENTS

- 1- Sensitized Plate** – Store between 2 and 8 °C. Contains: Polystyrene plate, divided into 12 strips of 8 wells each, impregnated with recombinant antigen of feline immunodeficiency virus (FIV).
- 2- Concentrated Conjugate (100x)** – Store between 2 and 8 °C. Contains: Solution with anti-feline antibodies linked to Peroxidase.
- 3- Concentrated Wash (20x)** – Store between 2 and 8 °C. Contains: Phosphate Buffer Solution, surfactant and preservative.
- 4- Diluent** – Store between 2 and 8 °C. Contains: Buffer Solution, surfactant, stabilizer and preservative.
- 5- Substrate TMB** – Store between 2 and 8 °C. Contains: Solution containing Tetramethylbenzidine (TMB), Citric Acid Solution and Urea Peroxide.
- 6- Stop Solution** – Store between 2 and 8 °C. Contains: 1N Sulfuric Acid Solution.
- 7- Negative Control** – Store between 2 and 8 °C. Contains: Negative serum for feline viral immunodeficiency and preservative.
- 8- Positive Control** – Store between 2 and 8 °C. Contains: Positive serum for feline viral immunodeficiency and preservative.
- 9- Plate Sealers**

PRESENTATIONS

Components	Presentation		
	1	2	3
	96 cavities	192 cavities	480 cavities
1- Sensitized Plate	1 unit	2 units	5 units
2- Concentrated Conjugate	1 x 250 µL	2 x 250 µL	5 x 250 µL
3- Concentrated Wash	1 x 50 mL	2 x 50 mL	5 x 50 mL
4- Diluent	1 x 60 mL	2 x 60 mL	5 x 60 mL
5- Substrate TMB	1 x 12 mL	2 x 12 mL	5 x 12 mL
6- Stop Solution	1 x 12 mL	2 x 12 mL	5 x 12 mL
7- Negative Control	1 x 1 mL	2 x 1 mL	5 x 1 mL
8- Positive Control	1 x 1 mL	2 x 1 mL	5 x 1 mL
9- Plate Sealers	3 units	6 units	15 units

EQUIPMENTS AND OPERATIONAL INPUTS
Materials contained in the kit:

- Reagents described in the previous table.
- Instructions for use (manual).

Materials needed, but not contained in the kit:

- 1 – Pipettes capable of dispensing volumes from 5 to 500 µL with a coefficient of variation lower than 1.5%.
- 2 – Re pipettor for repetitive pipetting of volumes of 500 µL with coefficient of variation less than 1.5% or multichannel pipette (optional).
- 3 – Microplate washer (optional).
- 4 – ELISA reader with absorbance capacity at 450 and 630 nm wavelength.
- 5 – Absorbent paper for drying microwells.
- 6 – Timer or stopwatch.
- 7 – Container for storing the diluted Wash Solution.
- 8 – Distilled or deionized water.
- 9 – Quality Control tools.
- 10 – Incubator at 37°C ± 2°C

STORAGE AND TRANSPORT CONDITIONS

The product storage temperature should be 2 to 8 °C. Transport can be carried out at room temperature (30°C) for up to 5 days. Keep away from light and avoid humidity. **Do not freeze.**

SPECIAL CARE

- 1– For professional veterinary *in vitro* diagnostic use only.
- 2– Follow the recommended methodology strictly to obtain accurate results.
- 3– The sachet containing the microplate must be opened only after it reaches room temperature. Replace strips from unused microwells in the sachet, seal and store at 2 to 8 °C.
- 4– The water used for cleaning materials must be fresh and free of contaminants.
- 5– Saturated deionizer columns release alkaline water, different ions, and oxidizing and reducing agents that can significantly alter the results.

6– Do not mix reagents from different kits or lots. Do not use expired components.

7– Do not eat, drink, or smoke in the testing area.

8– Do not pipette reagents or sample(s) by mouth. Do not use the same pipette tip for different samples.

9– The Positive and Negative Controls must be retested with each new assay performed.

10– The Stop Solution contains Sulfuric Acid which is a strong acid. Therefore, handle it with due care.

11– The professional must strictly follow the safety rules and routines when handling biological samples. The use of disposable gloves and other personal protective equipment is essential.

12– Pipette the reagents in the same order each time to minimize reaction time differences between microwells.

13– As a safety measure, cover the plate during the reaction.

14– Make sure that the bottom of the cavity is clean and dry and that there are no bubbles on the liquid surface before reading the plate. Do not allow the wells to dry out during the test.

15– Handle the kit components with due care, in order to avoid contamination. Use new or properly cleaned tips and containers for handling and do not expose reagents, especially Substrate, to strong light or Hypochlorite vapors during storage or incubation steps.

16– Products containing serum are potentially capable of transmitting diseases. Appropriate biosafety precautions must be taken when handling these materials.

17– We recommend applying local, state and federal environmental protection standards so that reagents and biological material are disposed of in accordance with current legislation.

18– For biosafety information or in case of accidents with the product, consult the SDS (Safety Data Sheet) available on the website www.bioclin.com.br or request it via Quibasa's Customer Support Service.

19– Do not use the product if the packaging is damaged.

20– It is essential that all instruments and equipment used are properly calibrated and undergo regular maintenance.

SAMPLES
Serum or Plasma (EDTA)

Use serum or plasma. Hemolyzed or highly lipemic samples should not be used. Samples may be stored refrigerated between 2 and 8°C, for a maximum period of 5 days. If the samples cannot be analyzed within 5 days, they can be stored for up to 30 days at –20°C..

PROCESS DESCRIPTION
Stability After Open

The results of the stability test show that the VETLISA FIV Ac kit is stable for at least 30 days after opening. This stability may vary depending on test and environmental conditions. Therefore, it is recommended to monitor product performance using the internal kit controls and the assay validation criteria.

PREPARATION OF WORK REAGENTS
Samples and Controls

In a test tube, dilute 15 µL of the sample in 300 µL of Diluent (Reagent No. 4), if performing the assay in duplicate. Cap the tube and gently vortex or manually mix by inversion. The dilutions cannot be stored.

ATTENTION: The Positive and Negative Controls are ready to use.

Washing Solution

Dilute the contents of the Concentrated Wash (Reagent No. 3) in 1000 mL of distilled or deionized water. After preparation, the solution can be stored at 2 to 30 °C up to 30 days. If crystallization occurs, heat to 37°C until dissolved.

VETERINARY USE
Conjugate Solution

Dilute the Concentrated Conjugate (Reagent N° 2) in the proportion of 1: 101 in Diluent (Reagent N° 4). Prepare the solution when performing the test.

To perform a test using all the wells of the kit, mix 110 µL of the Concentrated Conjugate in 11 mL of Diluent.

To perform an assay using 8 wells (1 strip), mix 10 µL of the Concentrated Conjugate in 1 mL of Diluent.

IMPORTANT: The diluted conjugate solution cannot be stored. Therefore, prepare only the amount necessary to perform the test.

Substrate

The Substrate is ready to use.

TECHNIQUE
For use in automatic equipment, consult the Customer Service Department (SAC).

Before starting the test, place all reagents, samples and controls to stabilize at room temperature (15 - 30 °C) for at least 40 minutes.

1- Separate the cavities to be used considering: Controls and Samples (it is recommended to test in duplicate). Return the unused strips of the microplate to the original sealed package.

2- Separate the first cavity for Blank (OPTIONAL).

3- Pipette 100 µL of the Negative Control and Positive Control in the previously determined wells. **Note: Controls are ready to use and do not require dilution.**

4- Pipette 100 µL of the samples previously diluted in the previously determined wells. In the Blank cavity (OPTIONAL), if you have made the option, pipette only 100 µL of the diluent.

5- Gently homogenize for ± 10 seconds. Cover the wells with plate sealers.

6- Incubate for 30 minutes ± 2 minutes in an incubator at 37°C ± 2°C.

7- Remove the plate sealer from the wells.

8- Discard the contents of the wells by aspiration (washer) or decanting (manual).

Use approximately 300 µL of Washing Solution, previously prepared, and perform a total of five (5) wash cycles. Shake for 3 seconds with each wash. To guarantee the drying of the plate, at the end of the washing, beat the plate for a few seconds on absorbent paper.

Note: Poor washing / drying can cause inadequate results.

9- Pipette 100 µL of Conjugate previously diluted in all wells.

10- Gently homogenize for ± 10 seconds. Cover the wells with plate sealers.

11- Incubate for 30 minutes ± 2 minutes in an incubator at 37°C ± 2°C.

12- Remove the plate sealer from the wells.

13- Repeat item 8.

14- Pipette 100 µL of Substrate into all wells.

15- Gently homogenize for ± 10 seconds. Cover the wells with a plate sealer.

16- Incubate for 10 minutes ± 2 minutes in an incubator at 37°C ± 2°C.

17- Remove the plate sealer from the wells.

18- Pipette 100 µL of Stop Solution into all wells.

19- Mix gently for approximately 10 seconds.

20- Read the plate using dual-wavelength filters: 450 nm / 630 nm, within 15 minutes (maximum).

TECHNIQUE VERIFICATION

Verify if the results obtained for reading the Blanck and the Controls are compatible with the values presented below:

ITEM	ABSORBANCE (DOUBLE FILTER)
Blank	< 0.100
Negative Control	0.100 to 0.300
Positive Control	> 0.800

If the values are outside the expected values, the technique should be repeated.

Calculations

Calculate Cut-off according to the following formula:

Cut-off = Average absorbance of the Negative Control + 0.250
--

Example:

ITEM	ABSORBANCIA
Negative Control	0.151
	0.149
Cut-off = Negative Control Average Absorbance + 0.250	$((0.151+0.149) / 2) + 0.250 = 0.400$

Calculate the Index by dividing the sample absorbance by the Cut-off value. Example

ITEM	ABSORBANCE
Sample	1.626
Cut-off Value	0.400
Index: Sample Abs. / Cut-off Value	$1.626 / 0.400 = 4.1$

Note: The data presented in the examples are for illustration only and cannot be used to calculate the results.

INTERPRETATION OF THE RESULT

After calculating the sample indices, consider the indices below to determine the results.

RESULTS	INDEX
Negative	< 0.9
Positive	> 1.1
Undetermined	$\geq 0.9 - \leq 1.1$

Note: An indeterminate result indicates that the sample showed absorbance within the assay's gray zone, between the cutoff limits established for the microplate ELISA used to detect anti-FIV antibodies, which does not allow confirmation or exclusion of exposure to, or infection with, feline immunodeficiency virus. Intermediate values may occur during the early phase of seroconversion (especially in the first weeks after exposure, during which antibody production may take up to approximately 60 days to become clearly detectable), in the presence of maternal antibodies in kittens, or due to pre-analytical interferences and analytical variations. It is recommended to repeat the test, preferably using a new sample, in order to rule out analytical variability or interfering factors. In cats with recent exposure risk or compatible clinical suspicion, confirmatory testing—such as PCR for detection of FIV provirus—may be indicated according to current guidelines. When recent exposure cannot be excluded, it is recommended to repeat serology after an interval that allows full development of the immune response, generally around 60 days after the last possible exposure. In kittens, the presence of maternal antibodies may cause indeterminate results or weak reactions, and in such cases, reevaluation should be performed when the cat reaches an age at which these antibodies have normally disappeared, ideally after six months of age, allowing determination of whether the detected antibodies reflect true infection or merely passive transfer. Until the infectious status is clarified, the cat should be considered potentially infectious to other felines. Final interpretation must be performed by the Veterinarian, integrating history, clinical signs, and complementary tests.

PROCEDURE LIMITATIONS

The interpretation of a diagnostic test should not be established on the basis of a single test. Other confirmatory tests must be included before a sample is considered positive. A negative result does not exclude the possibility of exposure. All results must be interpreted in conjunction with other clinical information available before the definitive diagnosis of the disease.

INTERFERENTS

No interference was observed for Triglycerides 1500 mg/dL, Ascorbic Acid 2 g/dL, Creatine 200 mg/dL, Bilirubin 1 g/dL, Albumin 2 g/dL, Hemoglobin 1000 mg/dL, Oxalic Acid 60 mg/dL, Factor Rheumatoid 980 IU/mL, C Reactive Protein 41.2 mg/dL and Antistreptolysin O 1023 IU/mL.

GROSS REACTIVITY

A cross-reactivity study was conducted, evaluating 88 samples that were negative for feline immunodeficiency virus (FIV) but positive for other infections. The following samples were tested: 10 samples positive for Feline Infectious Peritonitis (FIP), 7 samples positive for Feline Herpesvirus, 8 samples positive for feline calicivirus, 8 samples positive for Feline Panleukopenia (Parvovirus), 9 samples positive for *Mycoplasma haemofelis* (hemoplasmosis), 7 samples positive for *Bartonella henselae*, 10 samples positive for Cryptococcosis, 10 samples positive for Histoplasmosis, 10 samples positive for Toxoplasmosis, and 10 samples positive for feline leukemia virus (FeLV). No cross-reactivity was observed with any of the samples. Despite these results, the possibility of cross-reactivity cannot be completely ruled out. The final diagnosis should consider the patient's clinical data along with other laboratory and clinical findings.

INTERNAL QUALIT CONTROL

The Clinical Laboratory must have an internal quality control program, where procedures, standards, limits and tolerance for variations are clearly established. It is important to note that all measurement systems have a characteristic analytical variability, which must be monitored by the laboratories themselves. For that, it is recommended to use controls, which allow to evaluate the precision and accuracy of the measurements.

PRODUCT PERFORMANCE**ACCURACY****Repeatability**

Repeatability was calculated from 10 successive determinations, using 3 samples with different values. The following results were obtained:

REPEATABILITY	SAMPLE		
	1	2	3
Average	0.420	0.820	0.050
Standard Deviation	0.033	0.053	0.004
Coefficient of Variation (%)	7.866	6.465	8.411

Reproducibility

Reproducibility was calculated from 10 successive determinations over 3 consecutive days, using 3 samples with different values. The following results were obtained:

REPRODUCIBILITY	SAMPLE		
	1	2	3
Average	0.428	0.817	0.051
Standard Deviation	0.033	0.038	0.005
Coefficient of Variation (%)	7.739	4.669	9.174

CLINICAL SENSITIVITY AND SPECIFICITY

This product was tested in comparison with another immunoassay method. A total of 120 samples were analyzed, with 40 samples diagnosed as positive and 80 diagnosed as negative. The results show that the clinical sensitivity of the product is > 99.9% and the clinical specificity is 91.2%.

VETLISA FIV Ac	Method	Reference		Total
	Result	Negative	Positive	
	Negative	73	0	
Positive	7	40	47	
TOTAL RESULT		80	40	120

Clinical Sensitivity: > 99.9% (40 / 40) – 95% CI: 95%–100%

Clinical Specificity: 91.2% (73 / 80) – 95% CI: 85.31%–97.19%

Accuracy: 94.17% $[(73 + 40) / (80 + 40)]$

CLINICAL SIGNIFICANCE

Feline immunodeficiency virus (FIV) infection is caused by a retrovirus of the genus *Lentivirus*, which integrates into the host genome and establishes a persistent infection. The clinical course is generally divided into three phases: an acute phase, characterized by fever and mild lymphadenopathy; a prolonged asymptomatic phase, which may last for years; and a chronic phase, marked by progressive immunosuppression resulting from a decrease in CD4+ lymphocytes and the inversion of the CD4:CD8 ratio. In this phase, the animal becomes more susceptible to opportunistic infections, chronic inflammatory diseases, and neoplasms such as lymphomas. Diagnosis is traditionally based on the detection of specific antibodies against FIV. These antibodies typically appear between two and eight weeks after infection and remain detectable throughout the cat's life, making serology a sensitive and effective tool for identifying infected cats. However, negative results may occur in very early stages of infection, before seroconversion, or in cases of advanced immunosuppression. In young kittens, the presence of maternal antibodies may interfere with interpretation of results. Identifying infected cats is essential for appropriate clinical management, controlling the spread of the virus, and establishing preventive measures. The ELISA method offers high sensitivity, specificity, and the ability to process multiple samples simultaneously, allowing standardization, traceability, and application in laboratory routines and epidemiological studies. In this context, VETLISA FIV Ac stands out as a reliable and accurate diagnostic tool, enabling rapid and safe detection of anti-FIV antibodies and assisting veterinarians in diagnosis, monitoring, and clinical decision-making related to FIV infection.

BIBLIOGRAPHIC REFERENCES

- Hosie MJ, Addie DD, Belák S, et al. ABCD guidelines on feline immunodeficiency virus (FIV) infection. Journal of Feline Medicine and Surgery. 2009;11(7):575–584. Atualizado em 2017.
- Little S, Levy J, Hartmann K, et al. 2020 AAEP Feline Retrovirus Testing and Management Guidelines. Journal of Feline Medicine and Surgery. 2020;22(1):5–30.
- Sprößler F, Jongwattanapisan P, Luengyosuechakul S, et al. Prevalence and Risk Factors of Feline Immunodeficiency Virus and Feline Leukemia Virus Infection in Healthy Cats in Thailand. Frontiers in Veterinary Science. 2022;8:764217.
- Westman ME, Coggins SJ, van Dorsseleer M, et al. Feline immunodeficiency virus (FIV) infection in domestic pet cats in Australia and New Zealand: Guidelines for diagnosis, prevention and management. Australian Veterinary Journal. 2022;100(8):345–359.
- Westman ME, Malik R, Norris JM. Diagnosing feline immunodeficiency virus (FIV) and feline leukaemia virus (FeLV) infection: an update for clinicians. Australian Veterinary Journal. 2019;97(3):47–55.
- QUIBASA: Data from the Research and Development Department.

QUALITY GUARANTEE

Before being released for consumption, all Quibasa reagents are tested by the Quality Control Department. The quality of the reagents is assured until the validity date mentioned on the packaged, as long as they are stored and transported under the appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda.
 Rua Teles de Menezes, 92 - Santa Branca
 CEP 31565-130 - Belo Horizonte - MG - Brasil
 Tel.: +55 31 3439.5454 | E-mail: bioclin@bioclin.com.br
 CNPJ: 19.400.787/0001-07 - Industria Brasileira

CONSUMER SERVICE

Customer Support Service

Tel.: 0800 0315454 | E-mail: sac@bioclin.com.br

Product Licensed at the Ministry of Agriculture, Livestock and Supply from 11/11/2019 under the number 10.281/2019.

Technical manager: Dra. Camila Eckstein (CRMV/MG 20.611)

Review: March/2026

UNIVERSAL SYMBOLOGY