VITAMINA D

REF K300 / K301 / K302

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FINALIDADE

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- Garland, C. F. et al. The Role of Vitamin D in Cancer Prevention. Am J Public Health. 2006, 96(2): 252–261.
- Giovannucci, E. The epidemiology of vitamin D and cancer incidence and mortality: a review (United States). Cancer Causes Control. 2005, 16(2):83-95.
- Van den Bemd, G. J., Chang, GT. Vitamin D and vitamin D analogs in cancer treatment. Curr Drug Tar-gets. 2002, 3(1):85-94.
- Danik, J. S., Manson J. E, Vitamin D and Cardiovascular Disease. Curr Treat Options Cardiovasc Med. 2012, 14(4): 414–424.
- Reid I. R.; Bolland M. J. Role of Vitamin D Deficiency in Cardiovascular Disease. Heart. 2012, 98(8):609-614.
- Lavie, C. J.; Lee, J. H.; Milani, R. V. Vitamin D and Cardiovascular Disease Will It Live Up to its Hype? J Am Coll Cardiol. 2011, 58(15):1547-1556
- Holick, MF. Vitamin D Status: Measurement, Interpretation and Clinical Application. Ann Epidemiol. 2009, 19(2): 73–78.
- Morris H. A. Vitamin D: A Hormone for All Seasons-How Much is enough? Clin. Biochem. Rev., 2005, 26, 21-32.
- Bikle D. D. Vitamin D and the skin. J. Bone Miner. Metab. 2010, 28, 117-30.
- Zerwekh J. E. Blood biomarkers of vitamin D status. Am. J. Clin. Nutr. 2008, 87:1087S-91S.
- Malabanan A, Veronikis IE, Holick MF. Redefining Vitamin D insufficiency. Lancet, 1998, 351, 805-806.
- Holick MF et al. Evaluation, treatment, and prevention of vitamin D deficiency: an Endocrine Society clinical practice guideline. J Clin Endocrinol Metab. 2011;96 (7):1911-30.
- Borai, Anwar et al. “The pre-analytical stability of 25-hydroxyvitamin D: Storage and mixing effects.”Journal of clinical laboratory analysis, vol. 34.2 (2020) .
- QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

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- Garland, C. F. et al. The Role of Vitamin D in Cancer Prevention. Am J Public Health. 2006, 96(2): 252–261.
- Giovannucci, E. The epidemiology of vitamin D and cancer incidence and mortality: a review (United States). Cancer Causes Control. 2005, 16(2):83-95.
- Van den Bemd, G. J., Chang, GT. Vitamin D and vitamin D analogs in cancer treatment. Curr Drug Tar-gets. 2002, 3(1):85-94.
- Danik, J. S., Manson J. E, Vitamin D and Cardiovascular Disease. Curr Treat Options Cardiovasc Med. 2012, 14(4): 414–424.
- Reid I. R.; Bolland M. J. Role of Vitamin D Deficiency in Cardiovascular Disease. Heart. 2012, 98(8):609-614.
- Lavie, C. J.; Lee, J. H.; Milani, R. V. Vitamin D and Cardiovascular Disease Will It Live Up to its Hype? J Am Coll Cardiol. 2011, 58(15):1547-1556
- Holick, MF. Vitamin D Status: Measurement, Interpretation and Clinical Application. Ann Epidemiol. 2009, 19(2): 73–78.
- Morris H. A. Vitamin D: A Hormone for All Seasons-How Much is enough? Clin. Biochem. Rev., 2005, 26, 21-32.
- Bikle D. D. Vitamin D and the skin. J. Bone Miner. Metab. 2010, 28, 117-30.
- Zerwekh J. E. Blood biomarkers of vitamin D status. Am. J. Clin. Nutr. 2008, 87:1087S-91S.
- Malabanan A, Veronikis IE, Holick MF. Redefining Vitamin D insufficiency. Lancet, 1998, 351, 805-806.
- Holick MF et al. Evaluation, treatment, and prevention of vitamin D deficiency: an Endocrine Society clinical practice guideline. J Clin Endocrinol Metab. 201196 (7):1911-30.
- Borai, Anwar et al. “The pre-analytical stability of 25-hydroxyvitamin D: Storage and mixing effects.”Journal of clinical laboratory analysis, vol. 34.2 (2020) .
- QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

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INSTRUCTIONS FOR USE

FUNCTION

Test for quantitative determination of total Vitamin D (25-OH), in biological samples (plasma or serum) using immunoturbidimetric assay. Only for *in vitro* diagnostic use.

PRINCIPLE OF ACTION

Methodology: Immunoturbidimetry.

The vitamin D in the sample reacts with anti-Vitamin D antibodies conjugated to latex particles, forming the immunocomplex, which can be measured at 700 nm. The intensity of the turbidity of the solution is directly proportional to the amount of vitamin D present in the sample.

This reaction allows the concentration of vitamin D present in the sample to be quantified using a turbidimetric method. The suspended latex particles, covered with Anti-Vitamin D antibodies, are bound by the vitamin D present in the sample. The agglutination process causes an increase in particle size and, consequently, an increase in absorbance, which is measured by comparison with the calibrator of known concentration.

REAGENTS

Number 1 (R1) - Buffer - Store between 2 and 8 °C. **Do not freeze.** Contains: Buffer < 200 mmol/L, stabilizers, surfactant and preservative.

Number 2 (R2) - Vitamin D Latex - Store between 2 and 8 °C. **Do not freeze.** Contains: Suspension of Latex Particles sensitized with anti-Vitamin D antibody, stabilizers and preservative.

Number 3 (R3) - Calibrator Level 1 - Store between 2 and 8°C. **Do not freeze.** Contains: Vitamin D, stabilizers and preservative. **Potentially infectious.**

Number 4 (R4) - Calibrator Level 2 - Store between 2 and 8°C. **Do not freeze.** Contains: Vitamin D, stabilizers and preservative. **Potentially infectious.**

Number 5 (R5) - Calibrator Level 3 - Store between 2 and 8°C. **Do not freeze.** Contains: Vitamin D, stabilizers and preservative. **Potentially infectious.**

Number 6 (R6) - Calibrator Level 4 - Store between 2 and 8°C. **Do not freeze.** Contains: Vitamin D, stabilizers and preservative. **Potentially infectious.**

Number 7(R7) - Calibrator Level 5 - Store between 2 and 8°C. **Do not freeze.** Contains: Vitamin D, stabilizers and preservative. **Potentially infectious.**

Number 8 (R8) - Control Level 1 - Store between 2 and 8°C. **Do not freeze.** Contains: Vitamin D, stabilizers and preservative. **Potentially infectious.**

Number 9 (R9) - Control Level 2 - Store between 2 and 8°C. **Do not freeze.** Contains: Vitamin D, stabilizers and preservative. **Potentially infectious.**

Attention: The concentration of Vitamin D for calibrators and controls varies according to the batch - See bottle label.

PRESENTATION

Reagent	Presentation				
Buffer	1 x 20 mL	1 x 20 mL	1 x 40 mL	-	-
Vitamin D Latex	1 x 5 mL	1 x 5 mL	1 x 10 mL	-	-
Calibrator Level 1	1 x 1 mL	-	-	1 x 1 mL	-

Calibrator Level 2	1 x 1 mL	-	-	1 x 1 mL	-
Calibrator Level 3	1 x 1 mL	-	-	1 x 1 mL	-
Calibrator Level 4	1 x 1 mL	-	-	1 x 1 mL	-
Calibrator Level 5	1 x 1 mL	-	-	1 x 1 mL	-
Control Level 1	1 x 1 mL	-	-	-	1 x 1 mL
Control Level 2	1 x 1 mL	-	-	-	1 x 1 mL

OPERATIONAL EQUIPMENT AND CONSUMPTIONS

Automatic biochemical equipment, Bioclin Vitamin D Calibrator, Bioclin Vitamin D Controls. These items are found in the specialized market for articles for Clinical Analysis Laboratories.

STORAGE AND TRANSPORT CONDITIONS

The storage temperature should be 2 to 8°C. Transport at temperatures up to 30°C should not exceed 5 days. Keep away from light and avoid moisture. **Do not freeze.**

SPECIAL CARES

1- For *in vitro* diagnostic use only.

2- Strictly follow the proposed methodology to obtain exact results.

3- The water used to clean the material must be recent and free from contaminating agents.

4- Saturated deionizing columns release alkaline water, various ions and oxidizing and reducing agents, which can significantly alter the results.

5- The calibrators and controls were tested for anti-HIV, anti-HCV antibodies and HBs antigen using state-of-the-art methods and present negative results. The risk of infection cannot be excluded and the reagent should be handled with the same care as for patient plasma. **Potentially infectious.**

6- Perform a new calibration curve whenever necessary and for each batch of the product.

7- We recommend applying local, state and federal environmental protection standards so that reagents and biological material are disposed of in accordance with current legislation.

8- To obtain information related to biosafety or in case of accidents with the product, consult the SDS (Safety Data Sheet) available on the website www.bioclin.com.br or through a request through SAC (Customer Advisory Service) from Quibasa.

9- Do not use the product if the packaging is damaged.

10- It is essential that the instruments and equipment used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Only serum or plasma samples should be used for analysis. Samples can be stored under refrigeration, between 2 and 8 °C, for a maximum period of 1 day. If samples cannot be analyzed during this period, they can be stored for up to 3 months at -20 °C (freezer).

To obtain plasma, collect it using lytic heparin or EDTA as an anticoagulant.

PROCESS DESCRIPTION

Stability After Opening

The stability of the VITAMIN D product after opening (In Use / On Board – Product installed on equipment) is 28 days. This stability may vary according to the conditions of use of the product and environment. Therefore, it is suggested to monitor the product's performance using control sera.

TECHNIQUE

The Buffer, Vitamin D Latex, Calibrators and Controls reagents are ready to use.

To use the Vitamin D Buffer and Latex reagents, use the 5 VITAMIN D calibration levels as a calibrator and the two VITAMIN D control levels as a control.

The kit is only indicated for use in automatic biochemical analyzers. **Check the programming for the equipment on the website www.bioclin.com.br or via SAC 0800 031 5454.**

PROCESS LIMITATIONS

Product performance directly depends on the performance of the automatic equipment used. Therefore, to guarantee and obtain better test results, periodic checking of the equipment's optical and pipetting system is of fundamental importance.

INTERFERENCES

No interference was observed up to 40 mg/dL of unconjugated Bilirubin, 40 mg/dL of conjugated Bilirubin, 600 mg/dL of Hemoglobin, 1000 mg/dL of Triglycerides, 12 g/dL of Total Proteins and 200 IU/mL of Factor rheumatoid.

CROSS REACTIVITY

The VITAMIN D kit has 101% specificity for measuring 25-(OH)₂ Vitamin D₃ and 98% specificity for 25-(OH)₂ Vitamin D₂. No falsely elevated results were observed for the other analytes above, since the results obtained present a concentration lower than 2.7 ng/mL, which is the sensitivity of the assay, for concentrations of up to 0.58 ng/mL of 1.25-(OH)₂ Vitamin D₃, 0.58 ng/mL of 1,25-(OH)₂ Vitamin D₂, 100 ng/mL Vitamin D₃ and 100 ng/mL of Vitamin D₂.

Despite the results found, the possibility of other cross-reactivities cannot be completely ruled out. The final diagnosis must consider the patient's clinical data with other laboratory data.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control program, where procedures, standards, limits and tolerance for variations are clearly established. It is important to highlight that all measurement systems present a characteristic analytical variability, which must be monitored by the laboratories themselves. To this end, it is recommended to use controls, which allow evaluating the precision and accuracy of the dosages.

TRACEABILITY

The kit standard is normalized using a purified vitamin D preparation in mass and concentration.

REFERENCE VALUES

The reference values, in ng/mL, for the present method, were obtained through the determination of vitamin D in healthy male and female populations.

Reference values	
LEVEL	INTERVAL (ng/mL)
Deficiency	< 20
Insufficient	20-30
Ideal	30-100
Risk of toxicity	>100

These values should be used as a guideline, and each laboratory should create its range of reference values, according to the population served.
The results provided by this kit must be interpreted by the responsible medical professional, and are not the only criteria for determining the patient's diagnosis and/or treatment.

PRODUCT PERFORMANCE

ACCURACY

Comparison of methods and methodological specificity

The VITAMIN D kit was compared with another commercially available method for measuring Vitamin D. 171 analyzes were performed and the results were evaluated. The linear equation obtained was Y=1.04X-1.86 and the correlation coefficient was 0.96. With these results, it can be concluded that the kit presents good methodological specificity.

PRECISION

Repeatability

Repeatability was calculated from 40 successive determinations, using 2 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2
Concentration (ng/mL)	22.3	43.8
Standard deviation (ng/mL)	1.15	1.18
Coefficient of variation (%)	5	3

The product presents results within the expected variation range (coefficient of variation ≤ 10%).

Reproducibility

Reproducibility was calculated from 40 successive determinations over 3 consecutive days, using 2 samples with different concentrations, obtaining the following results:

	Sample 1	Sample 2
Concentration (ng/mL)	21.9	43.0
Standard deviation (ng/mL)	1.4	1.7
Coefficient of variation (%)	6	4

The product presents results within the expected variation range (coefficient of variation ≤ 10%).

SENSITIVITY

Sensitivity was calculated from 40 determinations of a sample free of Vitamin D. The average found was 0.5 ng/mL with a standard deviation of 0.7 ng/mL. The sensitivity, which indicates the detection limit of the method, corresponds to the mean plus 3 times the standard deviation and is equal to 2.71 ng/mL.

LINEARITY

The reaction is linear up to the last highest point of the calibration curve (around 141 ng/mL).

DIAGNOSTIC MEANING

The term Vitamin D is used to designate the compounds vitamin D₃ (cholecalciferol) and vitamin D₂ (ergocalciferol). Vitamin D₂ is obtained from dairy products, while vitamin D₃ is produced in the skin after exposure to ultraviolet radiation (sunlight).
In the liver, vitamin D is hydroxylated at its 25th carbon to form 25-OH vitamin D. Vitamin D deficiency can result in osteoporosis, bone fracture, rickets and osteomalacia. The sum of the concentrations of 25-OH Vitamin D₂ and 25-OH Vitamin D₃, in serum or plasma, is referred to as "Total 25-OH Vitamin D".

BIBLIOGRAPHIC REFERENCES

1. Garland, C. F. et al. The Role of Vitamin D in Cancer Prevention. Am J Public Health. 2006, 96(2): 252–261.
2. Giovannucci, E. The epidemiology of vitamin D and cancer incidence and mortality: a review (United States). Cancer Causes Control. 2005, 16(2):83-95.
3. Van den Bemd, G. J., Chang, GT. Vitamin D and vitamin D analogs in cancer treatment. Curr Drug Tar-gets. 2002, 3(1):85-94.
4. Danik, J. S., Manson J. E, Vitamin D and Cardiovascular Disease. Curr Treat Options Cardiovasc Med. 2012, 14(4): 414–424.
5. Reid I. R.; Bolland M. J. Role of Vitamin D Deficiency in Cardiovascular Disease. Heart. 2012, 98(8):609-614.
6. Lavie, C. J.; Lee, J. H.; Milani, R. V. Vitamin D and Cardiovascular Disease Will It Live Up to its Hype? J Am Coll Cardiol. 2011, 58(15):1547-1556
7. Holick, MF. Vitamin D Status: Measurement, Interpretation and Clinical Application. Ann Epidemiol. 2009, 19(2): 73–78.
8. Morris H. A. Vitamin D: A Hormone for All Seasons-How Much is enough? Clin. Biochem. Rev., 2005, 26, 21-32.
9. Bikle D. D. Vitamin D and the skin. J. Bone Miner. Metab. 2010, 28, 117-30.
10. Zerwekh J. E. Blood biomarkers of vitamin D status. Am. J. Clin. Nutr. 2008, 87:1087S-91S.
11. Malabanan A, Veronikis IE, Holick MF. Redefining Vitamin D insufficiency. Lancet, 1998, 351, 805-806.
12. Holick MF et al. Evaluation, treatment, and prevention of vitamin D deficiency: an Endocrine Society clinical practice guideline. J Clin Endocrinol Metab. 201196 (7):1911-30.
13. Borai, Anwar et al. "The pre-analytical stability of 25-hydroxyvitamin D: Storage and mixing effects."Journal of clinical laboratory analysis, vol. 34,2 (2020) .
14. QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

QUALITY WARRANTY

Before being released for consumption, all **Bioclin** reagents are tested by the Quality Control Department. The quality of the reagents is assured until the expiration date mentioned on the packaging, as long as they are stored and transported under appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda
Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Tel.: (31) 3439.5454 | E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

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