TT BIOCLIN

 **K219**

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FINALIDADE

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USAGE INSTRUCTIONS

FUNCTION

Method for manual or automated determination of Thrombin Time in citrated plasma. For *in vitro* diagnostic use only.

PRINCIPLE OF ACTION

It consists of measuring the plasma clotting time after the addition of thrombin (activated factor II), promoting the conversion of fibrinogen to fibrin.

REAGENT

Number 1 - Thrombin Solution Store at 2 to 8°C. Contains Bovine Thrombin 4.0 NIH/mL, Calcium Chloride (10 - 50 mmol/L), stabilizers and preservative.

Presentation ready to use.

PRESENTATION

Reagent 1
2 x 2 mL
6 x 2 mL
12 x 2 mL
6 x 5 mL

EQUIPMENT AND OPERATING INPUTS

Tips, pipettes, tubes, semi-automatic or automatic coagulation equipment (coagulometer), or timer (manual technique). These articles are available in the specialized market for Clinical Analysis Laboratories.

STORAGE AND TRANSPORT CONDITIONS

The storage temperature should be 2 to 8°C. Transport at temperatures up to 30°C should not exceed 5 days. Keep out of the light. **Do not freeze.**

SPECIAL CARE

- 1- **For professional *in vitro* diagnostic use only.**
- 2- Strictly follow the proposed methodology to obtain exact results.
- 3- The water used to clean the material must be fresh and free from contaminants.
- 4- The water level in the water bath must be higher than the level of the samples in the test tubes.
- 5- Reagent N° 1 contains Sodium Azide, irritating to skin and mucous membranes. Handle with care.
- 6- Avoid heating of reagent N° 1.
- 7- We recommend applying local, state and federal environmental protection standards so that the disposal of reagents and biological material is carried out in accordance with current legislation.

- 8- To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Chemical Safety Information Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Service of Chemicals) Customer Service) of Quibasa.
- 9- Do not use the product in case of damage to the packaging.

SAMPLES

Material: Plasma
Anticoagulant: Use the Bioclin A001 Citrate kit as anticoagulant according to the instructions for use, or 3.8% Sodium Citrate (0.130 M), in the proportion of 9 parts of blood to 1 of anticoagulant (as an example: 4.5 mL of blood + 0.5 mL of anticoagulant). Do not use EDTA or Heparin.
Collection: Blood must be obtained by venipuncture, avoiding hemolysis, prolonged tourniquet, bubble formation and tissue fluid aspiration.
Preparation: The blood should be mixed with the anticoagulant immediately after collection. Centrifuge immediately at 3000 rpm for 15 minutes. Remove the plasma without pipetting red cells or the yellow layer.

Sample Stability

Plasma: 4 hours stored at 15 - 25°C.
Plasma of a patient on heparin therapy: 2 hours stored at 15 - 25°C.
Frozen Plasma for a maximum of 2 weeks at -20°C.
Do not use clotted samples.

DESCRIPTION OF THE PROCESS

The stability of the TT Bioclin kit after opening is at least 14 days. This stability may vary according to the test, equipment and environment conditions. Therefore, it is suggested to monitor the product's performance using control plasmas.

TECHNIQUE

Perform tests and controls in duplicate.

Procedure

- Allow reagent N° 1 to reach room temperature before use.
- 1-Incubate 0.2 mL of plasma in a water bath at 37°C for 2 minutes.
- 2- Add 0.2 mL of reagent N° 1 unheated, simultaneously start the timer and record clot formation times.

RESULTS

The results of the Thrombin Time test are expressed in seconds. When the difference in the clotting times of the same sample is greater than 5%, it is recommended to repeat the test discarding the previous measurements.

INTERFERENTS

According to the scientific literature, hemolysis and lipemia are interferents that, in excess, can unpredictably alter the results of thrombin time dosing, which can result in falsely elevated or decreased results.

TRACEABILITY

The TT Bioclin kit is traceable to the reference material NIBSC WHO 2nd International Standard for Thrombin 01/580 US FDA/CBER Thrombin Standard Lot K.

INTERNAL QUALITY CONTROL

The Clinical Laboratory must have an internal quality control program, where procedures, limit norms and tolerance for variations are clearly established. It is important to note that all measurement systems have a characteristic analytical variability, which must be monitored by the laboratories themselves. For that, it is recommended to use controls, which allow to evaluate the precision and accuracy of the dosages.

REFERENCE

VALUES Normal values: <21 seconds.

The range should be used as a guide only. It is recommended that each Laboratory establish, in the population served, its own reference intervals. As an example, 40 plasmas from healthy donors were tested using the TT Bioclin reagent. An average clotting time of 17.5 seconds was obtained, with a standard deviation of 0.65 seconds. (97.5% of the tests obtained values below 19.0 seconds).

Aspects such as: quality of collection, conservation of the sample and equipment used can impact the sample's clotting time. It is recommended that each laboratory define its reference values for its population served.

PRODUCT PERFORMANCE

Accuracy

COMPARISON OF METHODS AND SPECIFICITY METHODOLOGICAL

The TT Bioclin kit was compared with another method for evaluating commercially available Thrombin Time. 48 analyzes were performed and the results were evaluated. The correlation coefficient is 0.9868. With these results, it can be concluded that the kit has good methodological specificity.

SPECIFICITY

The specificity of the kit is guaranteed by the addition of an exclusive source of thrombin, a specific activator of factor I (Fibrinogen), which eliminates the possibility of direct activation of other factors in the coagulation cascade.

REPETIBILITY

Repeatability was calculated from 10 consecutive determinations of three samples with different clotting times, obtaining the following results.

Repeatability	Sample 1	Sample 2	Sample 3
Clotting time (sec)	18.48	22.32	19.87
Standard deviation (sec)	0.17	0.20	0.26
Coefficient of variation (%)	0.92	0.90	1.31

REPRODUCIBILITY

Reproducibility was calculated from 10 consecutive determinations of three samples with different clotting times, for 3 consecutive days, obtaining the following results.

Repeatability	Sample 1	Sample 2	Sample 3
Clotting time (sec)	15.68	18.54	55.84
Standard deviation (sec)	0.29	0.30	1.98
Coefficient of variation (%)	1.82	1.64	3.55

ANALYTICAL SENSITIVITY

The sensitivity of the kit is defined by the concentration of thrombin coagulation factor activity in the kit, which is 4.0 NIH/mL. GAFFNEY, Patrick J .; EDGELL, Tracey A. The International and "NIH" Units for Thrombin-How Do They Compare? .Thrombosis and haemostasis, v. 73, n. 03, p. 900-903, 1995.

DIAGNOSTIC MEANING

Thrombin Time determines the time that fibrinogen is converted to fibrin in the presence of thrombin. Abnormality of clotting time in TT can be caused when fibrinogen concentrations are below 70 to 100 mg/dL, dysfibrinogenemia, disseminated intravascular coagulation (DIC), liver disease, in some paraproteinemias. Clotting times are heavily influenced by heparin and hirudin, or any treatment or disease that interferes with the conversion of fibrinogen to fibrin.

BIBLIOGRAPHICAL REFERENCES

- 1-SARAIVA, A. S. L. et al. Manual de diagnóstico laboratorial das Coagulopatias Hereditárias e Plaquetopatias. **Ministério da Saúde-Brasília-DF**, 2012.
- 2-HARRISON, Robert L.; BIROTTE, Regeina. The Thrombin Clotting Time. **American journal of clinical pathology**, v. 89, n. 1, p. 81-87, 1988.
- 3-FLANDERS, Michele M.; CRIST, Ronda; RODGERS, George M. Comparison of five thrombin time reagents. **Clinical chemistry**, v. 49, n. 1, p. 169-172, 2003.
- 4- QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

QUALITY ASSURANCE

Before being released for consumption, all Bioclin reagents are tested by the Quality Control Department. The quality of the reagents is ensured until the expiration date mentioned on the presentation packaging, as long as they are stored and transported under the appropriate conditions.

QUIBASA QUÍMICA BÁSICA Ltda

Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte/MG - Brasil

Phone.: +55 (31) 3439.5454

E-mail: bioclin@bioclin.com.br

CNPJ: 19.400.787/0001-07 – Made In Brazil

CONSUMER SERVICE

Customer Service

Tel .: 0800 0315454

E-mail: sac@bioclin.com.br

Registration number of TT Bioclin kit at ANVISA: 10269360327

Revision: November/2020

UNIVERSAL SYMBOLOGY

	CATALOG NUMBER		MADE BY
	LOT NUMBER		CONTROL
	MANUFACTURING DATE		POSITIVE CONTROL
	VALIDITY DATE (last day of the month)		NEGATIVE CONTROL
	TEMPERATURE LIMIT (store)		BIOLOGICAL RISK
	CONTENT IS SUFFICIENT FOR <N> TEST		FLAMMABLE
	SEE INSTRUCTIONS FOR USE		CORROSIVE
	IN VITRO DIAGNOSTIC PRODUCT		TOXIC
	EUROPEAN AUTHORIZED REPRESENTATIVE		CE MARK
	KEEP AWAY FROM SUNLIGHT		DO NOT USE IF PACKAGE IS DAMAGED
	DO NOT REUSE		PRODUCT STERILIZED
	CAUTION		DANGER