Bioclin

WRIGHT CORANTE HEMATOLÓGICO

REF

C006

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FINALIDADE

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USAGE INSTRUCTIONS

FUNCTION

Staining method for peripheral blood cells, bone marrow or cytological study of cellular elements. Only for *in vitro* diagnostic use.

PRINCIPLE OF ACTION

Methodology: Dye according to Wright
Wright Dye is a mixture of dye with neutral characteristics, which stain the nuclear and cytoplasmic components of cells.

REAGENTS

Number 1 - Wright Dye - Store between 15 and 30°C.
Contains: Eosin Methylene Blue and solubilizer.

PRESENTATION

Reagent	Volume
Nº1	1000 mL

EQUIPMENTS AND OPERATIONAL INPUTS

Support for staining, watch or stopwatch, microscope, pipettes and blades. These items can be found in the specialized market of articles for Clinical Analysis Laboratories.

TRANSPORTATION AND STORAGE CONDITIONS

The storage and transportation temperature should be between 15 to 30°C. Protect from light and avoid moisture.

SPECIAL CARE

- 1- For *in vitro* diagnostic use only.**
- 2-** Strictly follow the methodology proposed.
- 3-** Water used in material cleaning must be recent and free of contaminants.
- 4-** Handle reagent with care, avoid contact with skin and mucous.
- 5-** The slides for the fabrication of smears should be clean and fat-free.
- 6-** A drop of blood used for making the smear should not be too large because the larger the droplet the more thicker the smear.
- 7-** We recommend applying the local, state and federal rules for environmental protection, so that disposal of reagents and biological material can be made in accordance with current legislation.
- 8-** To obtain information related to biosafety or in case of accidents with the product, consult the MSDS (Material Safety Data Sheet) available on the website www.bioclin.com.br or upon request by the SAC (Advisory Service Customer) of Quibasa.

- 9-** Do not use the product in case of damaged packaging.
- 10-** It is essential that the instruments and equipments used are properly calibrated and subjected to periodic maintenance.

SAMPLES

Peripheral blood collected by digital puncture or venous collection, without anticoagulant or with EDTA Bioclin. The smears should be done preferably without the use of anticoagulants, because the anticoagulant may alter the leukocyte morphology. Therefore, slides made from blood collected with anticoagulants should be made up to 30 minutes after collection to avoid cell deformations.⁴

PROCESS DESCRIPTION

TECHNIQUE

Place slide on staining rack.

Note: Due to the pH variation of tap water, Bioclin recommends the use of deionized water in steps 2 and 3 of the staining technique.

- 1-** Cover the smear with 20 drops of dye. Wait one to three minutes.
- 2-** Add 20 drops of deionized water. Homogenize the dye and water mixture. Wait for the coloration from three to five minutes.
- 3-** Rinse the slide with deionized water and carefully wipe the part opposite the smear for better viewing of the slide under the microscope.
- 4-** Allow it to dry in an upright position.
- 5-** Examine with an immersion objective.

RESULTS

Erythrocytes: Pinkish.
Platelets: Bluish color.
Neutrophiles: Dark blue nucleus and pale pink cytoplasm with granules ranging from pink to blue.
Lymphocytes: Violet blue nucleus and blue cytoplasm.
Basophiles: purple to dark blue nucleus and cytoplasm with bulky dark blue granules.
Monocytes: Violet blue nucleus and pale blue cytoplasm.
Eosinophiles: Blue nucleus and pale pink cytoplasm with bulky granules ranging from red to orange.

PROCEDURE LIMITATIONS

Smears that present a blue or green color are alkaline due to the prolonged action of the dye, thick smears or insufficient washing of the slide.

Smears that are present red color are overly acidic due to weak staining, washing time or excessive acidity of the dye.
Precipitation may appear when the slides are unwashed, dust or drying of the dye occurs during the staining process.
Eventually the water used by the laboratory does not provide a good coloring. In this case use a Phosphate Tampon (pH 7,2), using the following formula:

1 - Solution A	
KH ₂ PO ₄	9,1 g
Distilled H ₂ O q.s.p.	1.000 mL

2 - Solution B	
Na ₂ HPO ₄	9,5 g
Distilled H ₂ O q.s.p.	1.000 mL

Usage Tampon

Solution A	28 mL
Solution B	72 mL

PRODUCT PERFORMANCE

ACCURACY
Comparison of Methods and Methodological Specificity
Wright Hemotologic Dye kit was compared with other smears coloring methods commercially available. According to the coloring comparison of 20 clinical samples, the method shows an excellent correlation.

PRECISION

Repeatability

The repeatability was calculated from 10 successive colorations, using 1 sample, obtaining the following results:

Number	Sample 1
1	Standard Coloring
2	Standard Coloring
3	Standard Coloring
4	Standard Coloring
5	Standard Coloring
6	Standard Coloring
7	Standard Coloring
8	Standard Coloring
9	Standard Coloring
10	Standard Coloring

Reproducibility

The reproducibility was calculated from 10 successive colorations for 3 consecutive days, using 1 sample, obtaining the following results:

Days	N° of Repetitions	Expected Results	Found Results
1° Day	10	Standard Coloring	Standard Coloring
2° Day	10	Standard Coloring	Standard Coloring
3° Day	10	Standard Coloring	Standard Coloring

BIBLIOGRAPHIC REFERENCES

1 - WALF., Pratical Clinical Hematology Interpretation and Techniques, 375-376,1973.
2 - LIMA, O. A.; SOARES, J.B.; GALIZZI, J.; GRECO, J. B.; CANÇADO J. R., 12-21, 1992.
3 - CARVALHO, M. G.; SILVA, M. B. S., Hematologia, Técnicas Laboratoriais e Interpretação, 28-32, 1988.
4 - NARASIMHA, A., KUMAR, H., PRASAD, CSBR. Anticoagulant induced artefacts in peripheral blood smears. Indian Journal of Hematology & Blood Transfusion. 2008;24(2):43-48.
5 - QUIBASA: Dados do Departamento de Pesquisa e Desenvolvimento.

QUALITY ASSURANCE

Before being released for consumption, all **Bioclin** reagents are tested by the Department of Quality Control. The quality of reagents is assured until expiration date stated on the presentation packaging, when stored and transported under appropriate conditions.

 **QUIBASA QUÍMICA BÁSICA Ltda**
Rua Teles de Menezes, 92 - Santa Branca
CEP 31565-130 - Belo Horizonte - MG - Brasil
Phone.: +55 31 3439.5454 | E-mail: bioclin@bioclin.com.br
CNPJ: 19.400.787/0001-07 - Made in Brazil

CUSTOMER SERVICE

Customer Advisory Service
Phone.: 0800 0315454 | E-mail: sac@bioclin.com.br

ANVISA registration for Wright Hematologic Dye kit: 10269360095

Review: August/2022

UNIVERSAL SYMBOLOGY

CATALOG NUMBER

LOT NUMBER

MANUFACTURING DATE

VALIDITY DATE
(last day of the month)

TEMPERATURE LIMIT
(store)

CONTENT IS SUFFICIENT
FOR <N> TEST

SEE INSTRUCTIONS
FOR USE

IN VITRO DIAGNOSTIC
PRODUCT

KEEP AWAY
FROM SUNLIGHT

DO NOT REUSE

CAUTION

MADE BY

CONTROL

POSITIVE CONTROL

NEGATIVE CONTROL

BIOLOGICAL RISK

FLAMMABLE

CORROSIVE

TOXIC

DO NOT USE IF
PACKAGE IS
DAMAGED

PRODUCT
STERILIZED

DANGER