Catálogo: K040-1

ALBUMINA MONOREAGENTE

Nº de Testes: 1250

Versão: 08/08/2018

Cód.: 00

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**K040-1**

Reagente Nº 1 – Reagente de Cor – 1 x 250 mL

Reagente Nº 2 – Padrão – 1 x 2 mL

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<b>TEST CODE</b>	<b>0 000</b>	<b>SAMPLE BLANK ABSORBANCE CHECKS</b>	
Name	ALBUMINA	Absorbance Limit	-0,1   3
Abbreviation	ALB	Reagent Abs Limit	-0,1   3
Mode	Ponto de Desfecho	Substrate Depletion	0 Abs
Reagent Onboard Stability	1 day	Result Type	Quantitative
Reagent Low Vol Warming	5 Tests	Primary Units	g/dL
Service Assay		Unit Conversion Factor A	-
Suppress Test Result		Unit Conversion Factor B	-
		Secondary Units	-
Serum	✓	Decimals	1
Urine	-		
WholeBlood (disabled)	-	<b>REFERENCE RANGES</b>	
CSF (disabled)	-	<b>Serum</b>	
(disabled)	-	Panic	0   0
Replicates	1	Male	3.5   3.5
Sample Blank	On Request	Female	3.5   3.5
Sample Volume	2.0 µL	Pediatric	0   0
Reagent Volume -2.25 min (R1)	200 µL	Urine	-   -
Reagent Volume 2.90 min (R2)	00 µL	Whole Blood	-   -
Reagent Volume 4.70 min (R3)	00 µL	<b>CALIBRATION TYPE</b>	Calibrator
<b>WAVELENGTH</b>	620 nm	Correlation	Factor Offset
			1.000 0.000
<b>Incubation Time</b>	4.5 min	Calibrator Name	<a href="#">Biocal – K072</a>
<b>ABSORBANCE CHECKS</b>		Replicates	1
Absorbance Limit (Abs)	-0,1   3	Expiration Time	7 days
Reagent Abs Limit (Abs)	-0,1   3	Auto Predilution	-
Reaction Direction	Increasing	Reagent Blank	✓
Reagent Deviation Abs	-	Reag. Blank Expiration	1 day
Substrate Depletion	0 Abs	Reag. Blank Replicates	1
Alinearity Limit	10 %	Reag. Blank Abs Limits	-0,1   3
		Number of Points	1
		Algorithm	Linear
		Algorithm Inaccuracy	5%
		Extended Abs Limits	1   1
		Extended Concent Limits	0   10

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