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Reagente Nº 1 – Tampão – 2 x 24 mL

Reagente Nº 2 – Substrato – 2 x 6 mL

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TEST CODE	0 000	SAMPLE BLANK ABSORBANCE CHECKS
Name	GAMA GT	Absorbance Limit -0,1 3
Abbreviation	GGT	Reagent Abs Limit -0,1 3
Mode	Kinetic	Substrate Depletion 0 Abs
Reagent Onboard Stability	1 day	Result Type Quantitative
Reagent Low Vol Warming	5 Tests	Primary Units U/L
Service Assay		Unit Conversion Factor A -
Suppress Test Result		Unit Conversion Factor B -
Serum	✓	Secondary Units -
Urine	-	Decimals 1
WholeBlood (disabled)	-	
CSF (disabled)	-	REFERENCE RANGES
(disabled)	-	Serum
Replicates	1	Panic 0.0 0.0
Sample Blank	No	Male 0.0 55.0
Sample Volume	10 µL	Female 0.0 38.0
Reagent Volume -2.25 min (R1)	160 µL	Pediatric 0.0 0.0
Reagent Volume 2.90 min (R2)	40 µL	Urine - -
Reagent Volume 4.70 min (R3)	00 µL	Whole Blood - -
WAVELENGTH	405 nm	CALIBRATION TYPE
Kinetic	Delay 50 sec	Correlation
Min Time 159 sec		Factor 1.000
ABSORBANCE CHECKS		Offset 0.000
Absorbance Limit (Abs)	-0,1 3	Calibrator Name Biocal – K072
Reagent Abs Limit (Abs)	-0,1 3	Replicates 1
Reaction Direction	Increasing	Expiration Time 7 days
Reagent Deviation Abs	-	Auto Predilution -
Substrate Depletion	0 Abs	Reagent Blank ✓
Alinearity Limit	10 %	Reag. Blank Expiration 1 day
		Reag. Blank Replicates 1
		Reag. Blank Abs Limits -0,1 3
		Number of Points 1
		Algorithm Linear
		Algorithm Inaccuracy 5%
		Extended Abs Limits 1 1
		Extended Concent Limits 0 10

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