Catálogo: K084-2

LACTATO

Nº de Testes: 166

Versão: 08/08/2018

Cód.: 00

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K084-2

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TEST CODE	0 000	SAMPLE BLANK ABSORBANCE CHECKS	
Name	LACTATO	Absorbance Limit	-0,1 3
Abbreviation	LAC	Reagent Abs Limit	-0,1 3
Mode	Ponto de Desfecho	Substrate Depletion	0 Abs
Reagent Onboard Stability	1 day	Result Type	Quantitative
Reagent Low Vol Warming	5 Tests	Primary Units	mg/dL
Service Assay		Unit Conversion Factor A	-
Suppress Test Result		Unit Conversion Factor B	-
Serum	✓	Secondary Units	-
Urine	-	Decimals	1
WholeBlood (disabled)	-	REFERENCE RANGES	
CSF (disabled)	-	Serum	
(disabled)	-	Panic	0 0
Replicates	1	Male	4.5 19.8
Sample Blank	Yes	Female	4.5 19.8
Sample Volume	3,0 µL	Pediatric	0 0
Reagent Volume -2.25 min (R1)	240 µL	Urine	- -
Reagent Volume 2.90 min (R2)	60 µL	Whole Blood	- -
Reagent Volume 4.70 min (R3)	00 µL	CALIBRATION TYPE	
WAVELENGTH	340 nm	Correlation	Calibrator
Incubation Time	4.5 min	Factor	1.000
ABSORBANCE CHECKS		Offset	0.000
Absorbance Limit (Abs)	-0,1 3	Calibrator Name	Biocal – K072
Reagent Abs Limit (Abs)	-0,1 3	Replicates	1
Reaction Direction	Increasing	Expiration Time	7 days
Reagent Deviation Abs	-	Auto Predilution	-
Substrate Depletion	0 Abs	Reagent Blank	✓
Alinearity Limit	10 %	Reag. Blank Expiration	1 day
		Reag. Blank Replicates	1
		Reag. Blank Abs Limits	-0,1 3
		Number of Points	1
		Algorithm	Linear
		Algorithm Inaccuracy	5%
		Extended Abs Limits	1 1
		Extended Concent Limits	0 10

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