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Reagente Nº 1 – Reagente de Cor – 2 x 100 mL

Reagente Nº 2 – Padrão – 1 x 3 mL

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TEST CODE	0 000	SAMPLE BLANK ABSORBANCE CHECKS	
Name	MAGNÉSIO	Absorbance Limit	-0,1 3
Abbreviation	MAG	Reagent Abs Limit	-0,1 3
Mode	Ponto de Desfecho	Substrate Depletion	0 Abs
Reagent Onboard Stability	1 day	Result Type	Quantitative
Reagent Low Vol Warming	5 Tests	Primary Units	mg/dL
Service Assay		Unit Conversion Factor A	-
Suppress Test Result		Unit Conversion Factor B	-
Serum	✓	Secondary Units	-
Urine	-	Decimals	1
WholeBlood (disabled)	-	REFERENCE RANGES	
CSF (disabled)	-	Serum	
(disabled)	-	Panic	0.0 0.0
Replicates	1	Male	1.6 2.4
Sample Blank	On Request	Female	1.6 2.4
Sample Volume	2,0 µL	Pediatric	0.0 0.0
Reagent Volume -2.25 min (R1)	200 µL	Urine	- -
Reagent Volume 2.90 min (R2)	00 µL	Whole Blood	- -
Reagent Volume 4.70 min (R3)	00 µL	CALIBRATION TYPE	Calibrator
WAVELENGTH	505 nm	Correlation	Factor Offset
Incubation Time	1.5 min	Calibrator Name	1.000 0.000
ABSORBANCE CHECKS		Replicates	Biocal – K072
Absorbance Limit (Abs)	-0,1 3	Expiration Time	1
Reagent Abs Limit (Abs)	-0,1 3	Auto Predilution	7 days
Reaction Direction	Increasing	Reagent Blank	-
Reagent Deviation Abs	-	Reagent Blank Expiration	✓
Substrate Depletion	0 Abs	Reagent Blank Replicates	1 day
Alinearity Limit	10 %	Reagent Blank Abs Limits	1
		Reagent Blank Abs Limits	-0,1 3
		Number of Points	1
		Algorithm	Linear
		Algorithm Inaccuracy	5%
		Extended Abs Limits	1 1
		Extended Concent Limits	0 10

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