Catálogo: K049-6

TRANSAMINASE ALT (TGP) CINÉTICA

Nº de Testes: 500

Versão: 08/08/2018

Cód.: 00

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**K049-6**

Reagente Nº 1 – Substrato – 2 x 40 mL

Reagente Nº 2 – Coenzima – 2 x 10 mL

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<b>TEST CODE</b>	<b>0 000</b>	<b>SAMPLE BLANK ABSORBANCE CHECKS</b>	
Name	ALT_TGP	Absorbance Limit	-0,1   3
Abbreviation	ALT_TGP	Reagent Abs Limit	-0,1   3
Mode	Kinetic	Substrate Depletion	0 Abs
Reagent Onboard Stability	1 day	Result Type	Quantitative
Reagent Low Vol Warming	5 Tests	Primary Units	U/L
Service Assay		Unit Conversion Factor A	-
Suppress Test Result		Unit Conversion Factor B	-
Serum	✓	Secondary Units	-
Urine	-	Decimals	1
WholeBlood (disabled)	-	<b>REFERENCE RANGES</b>	
CSF (disabled)	-	<b>Serum</b>	
(disabled)	-	Panic	0.0   0.0
Replicates	1	Male	5.0   38.0
Sample Blank	No	Female	5.0   38.0
Sample Volume	20 µL	Pediatric	0.0   0.0
Reagent Volume -2.25 min (R1)	200 µL	Urine	-   -
Reagent Volume 2.90 min (R2)	00 µL	Whole Blood	-   -
Reagent Volume 4.70 min (R3)	00 µL	<b>CALIBRATION TYPE</b>	Calibrator
<b>WAVELENGTH</b>	340 nm	Correlation	Factor 1.000 Offset 0.000
<b>Kinetic</b>	Delay 50 sec Min Time 186 sec	Calibrator Name	<a href="#">Biocal – K072</a>
<b>ABSORBANCE CHECKS</b>		Replicates	1
Absorbance Limit (Abs)	-0,1   3	Expiration Time	7 days
Reagent Abs Limit (Abs)	-0,1   3	Auto Predilution	-
Reaction Direction	Decreasing	Reagent Blank	✓
Reagent Deviation Abs	-	Reag. Blank Expiration	1 day
Substrate Depletion	0 Abs	Reag. Blank Replicates	1
Alinearity Limit	10 %	Reag. Blank Abs Limits	-0,1   3
		Number of Points	1
		Algorithm	Linear
		Algorithm Inaccuracy	5%
		Extended Abs Limits	1   1
		Extended Concent Limits	0   10

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