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Reagente Nº 1 – Substrato – 2 x 40 mL

Reagente Nº 2 – Coenzima – 2 x 10 mL

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TEST CODE	0 000	SAMPLE BLANK ABSORBANCE CHECKS
Name	ALT_TGP	Absorbance Limit -0,1 3
Abbreviation	ALT_TGP	Reagent Abs Limit -0,1 3
Mode	Kinetic	Substrate Depletion 0 Abs
Reagent Onboard Stability	1 day	Result Type Quantitative
Reagent Low Vol Warming	5 Tests	Primary Units U/L
Service Assay		Unit Conversion Factor A -
Suppress Test Result		Unit Conversion Factor B -
Serum	✓	Secondary Units -
Urine	-	Decimals 1
WholeBlood (disabled)	-	REFERENCE RANGES
CSF (disabled)	-	Serum
(disabled)	-	Panic 0.0 0.0
Replicates	1	Male 5.0 38.0
Sample Blank	No	Female 5.0 38.0
Sample Volume	20 µL	Pediatric 0.0 0.0
Reagent Volume -2.25 min (R1)	160 µL	Urine - -
Reagent Volume 2.90 min (R2)	40 µL	Whole Blood - -
Reagent Volume 4.70 min (R3)	00 µL	CALIBRATION TYPE
WAVELENGTH	340 nm	Correlation Calibrator
Kinetic	Delay 50 sec	Factor 1.000
	Min Time 186 sec	Offset 0.000
ABSORBANCE CHECKS		Calibrator Name Biocal – K072
Absorbance Limit (Abs)	-0,1 3	Replicates 1
Reagent Abs Limit (Abs)	-0,1 3	Expiration Time 7 days
Reaction Direction	Decreasing	Auto Predilution -
Reagent Deviation Abs	-	Reagent Blank ✓
Substrate Depletion	0 Abs	Reag. Blank Expiration 1 day
Alinearity Limit	10 %	Reag. Blank Replicates 1
		Reag. Blank Abs Limits -0,1 3
		Number of Points 1
		Algorithm Linear
		Algorithm Inaccuracy 5%
		Extended Abs Limits 1 1
		Extended Concent Limits 0 10

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