

Versão: 17/08/2020

| Test Name | Code | Local Code | Channel | | Release | Modified on |
|-----------|------|------------|---------|--------------------------------------------|---------|-------------|
| CLORETOS | CLOR | CLOR | 00 | <input checked="" type="checkbox"/> Enable | | 17/08/2020 |

Result

| | |
|------|-------|
| Unit | mEq/L |
|------|-------|

| | |
|------------------|---|
| Decimal Position | 0 |
|------------------|---|

Manual patient validation

Correlation

Slope Intercept

| | |
|---------|---------|
| 1.00000 | 0.00000 |
|---------|---------|

| Reagent short name | Reagent number |
|--------------------|----------------|
| CLOR | |

365

Man/Default Woman Child

| | | | |
|--------------------------------------------------------------------------------------------------|-----|-----|--|
|  High check | 112 | 112 | |
|--------------------------------------------------------------------------------------------------|-----|-----|--|

Man/Default Woman Child

☐ Low check

☐ High check ☐ ☐ ☐

Catálogo: K221

CLORETOS MONORREAGENTE

Nº de testes: 250

Versão: 17/08/2020

Observações: Para garantir o melhor desempenho nas análises, a Bioclin recomenda o uso do Biocal - K072 para Calibração e o uso do Biocontrol N – K073 e do Biocontrol P - K074 para Controle de Qualidade.

| | | | |
|-----------|------|------------|---------|
| Test Name | Code | Local Code | Channel |
| CLORETOS | CLOR | CLOR | 00 |

Pre-dilution

| | |
|------|--------------------|
| Type | Calibrador Diluent |
| | |

| | | | |
|----------|----------|----------|----------|
| Factor 1 | Factor 2 | Factor 3 | Factor 4 |
| | | | |
| Factor 5 | Factor 6 | Factor 7 | Factor 8 |
| | | | |

Calibration

| | |
|--------------------|---------------|
| Calibration mode | Slope Average |
| Level | 1 |
| Calibration Factor | |
| Run(s) | 2 |

Validity

| | |
|------------------------------------------------|----------------------------------------|
| <input checked="" type="checkbox"/> On request | <input type="checkbox"/> Time Validity |
| Interval | Time Unit |
| 0 | Days |

Factor calibration

| | |
|----------------------|--|
| Low limit check | |
| High limit check | |
| Relative limit check | |

| | | | |
|-------------|-------------------------------------|-----------|-------------------------------------|
| Dev_Rep (%) | <input checked="" type="checkbox"/> | Dev_C (%) | <input checked="" type="checkbox"/> |
|-------------|-------------------------------------|-----------|-------------------------------------|

| |
|-----------------|
| Calibrator Used |
| Biocal |

Checks

| | |
|-------------------------------------|--------------------------------|
| <input checked="" type="checkbox"/> | Reagent Limit Absorbance Check |
|-------------------------------------|--------------------------------|

| | |
|-------------------|----------|
| Reagent Range Low | -3.00000 |
|-------------------|----------|

| | |
|--------------------|---------|
| Reagent Range High | 3.00000 |
|--------------------|---------|

| | |
|-------------------------------------|--------------------------------------|
| <input checked="" type="checkbox"/> | Reagent Blank Limit Absorbance Check |
|-------------------------------------|--------------------------------------|

| | |
|-------------------------|----------|
| Blank Range - Low limit | -3.00000 |
|-------------------------|----------|

| | |
|--------------------------|---------|
| Blank Range - High limit | 3.00000 |
|--------------------------|---------|

Control required

| | |
|-----------------------------------------------|--------------|
| <input checked="" type="checkbox"/> Control 1 | Biocontrol N |
| <input checked="" type="checkbox"/> Control 2 | Biocontrol P |
| <input type="checkbox"/> Control 3 | |

Validity Backup

| | |
|--------------------------|------------------------------------------------|
| <input type="checkbox"/> | Backup time frame without calibration required |
|--------------------------|------------------------------------------------|

| | |
|----------|-----------|
| Interval | Time Unit |
| 0 | Days |

Versão: 17/08/2020

Observações: Para garantir o melhor desempenho nas análises, a Bioclin recomenda o uso do Biocal - K072 para Calibração e o uso do Biocontrol N – K073 e do Biocontrol P - K074 para Controle de Qualidade.

Channel

00

Blank

Diluent

H20

10

X

Mixing Speed

40

Analysis Sequence

H2O Vol (μL)

10

Page 10 of 10

Catálogo: K221

CLORETOS MONORREAGENTE

Nº de testes: 250

Versão: 17/08/2020

Observações: Para garantir o melhor desempenho nas análises, a Bioclin recomenda o uso do Biocal - K072 para Calibração e o uso do Biocontrol N – K073 e do Biocontrol P - K074 para Controle de Qualidade.

Test Name

CLORETOS

Code

CLORETOS

Local Code

CLORETOS

Channel

00

Correlation Factor

Slope 1.00000

Intercept 0.000

Reaction Direction

X

Reaction Direction Check

Reaction Direction Increase

Sample Limit Check

Sample Limit

Sample Limit Cycle

1

Definition

Calculation Type End Point

OD Deviation check

Linear regression

r2 restored SD

First point

First point threshold

Last point

SD factor

Antigen excess activation

Antigen excess limit (%)

Antigen excess point

1

Reaction Limit check

Reaction limit absorbance

Cycle 1

First Reading

Cycle 0

Last Reading

Cycle 10