

Versão: 23/10/2020

Test Name	Code	Local Code	Channel		Release	Modified on
FOSFATASE ALCALINA	FALC	FALC	00	<input checked="" type="checkbox"/> Enable		23/10/2020

Result

Unit	U/L
------	-----

Decimal Position	0
------------------	---

☐ Manual patient validation

Correlation

Slope	Intercept
1.00000	0.00000

Delta Check Validity	Absolut Variation	Relative Variation
365		

Reference Range (μLmol/L)

	Man/Default	Woman	Child
Low check	27	27	75
High check	100	100	390

Man/Default	Woman	Child
<input type="text"/>	<input type="text"/>	<input type="text"/>
<input type="text"/>	<input type="text"/>	<input type="text"/>

☐ Low check

High check			
------------	--	--	--

Catálogo: K224

FOSFATASE ALCALINA IFCC

Nº de testes: 500

Versão: 23/10/2020

Observações: Para garantir o melhor desempenho nas análises, a Bioclin recomenda o uso do Biocal - K072 para Calibração e o uso do Biocontrol N – K073 e do Biocontrol P - K074 para Controle de Qualidade.

Test Name	Code	Local Code	Channel
FOSFATASE ALCALINA	FALC	FALC	00

Pre-dilution

Type	Calibrador Diluent

Factor 1	Factor 2	Factor 3	Factor 4
Factor 5	Factor 6	Factor 7	Factor 8

Calibration

Calibration mode	Slope Average
Level	1
Calibration Factor	
Run(s)	2

Validity

<input checked="" type="checkbox"/> On request	<input type="checkbox"/> Time Validity
Interval	Time Unit
0	Days

Factor calibration

Low limit check	
High limit check	
Relative limit check	

Dev_Rep (%)	<input checked="" type="checkbox"/>	Dev_C (%)	<input checked="" type="checkbox"/>
-------------	-------------------------------------	-----------	-------------------------------------

Calibrator Used
Biocal

Checks

<input checked="" type="checkbox"/>	Reagent Limit Absorbance Check
-------------------------------------	--------------------------------

Reagent Range Low	-3.00000
-------------------	----------

Reagent Range High	3.00000
--------------------	---------

<input checked="" type="checkbox"/>	Reagent Blank Limit Absorbance Check
-------------------------------------	--------------------------------------

Blank Range - Low limit	-3.00000
-------------------------	----------

Blank Range - High limit	3.00000
--------------------------	---------

Control required

<input checked="" type="checkbox"/> Control 1	Biocontrol N
<input checked="" type="checkbox"/> Control 2	Biocontrol P
<input type="checkbox"/> Control 3	

Validity Backup

<input type="checkbox"/>	Backup time frame without calibration required
--------------------------	--

Interval	Time Unit
0	Days

Versão: 23/10/2020

Observações: Para garantir o melhor desempenho nas análises, a Bioclin recomenda o uso do Biocal - K072 para Calibração e o uso do Biocontrol N – K073 e do Biocontrol P - K074 para Controle de Qualidade.

Channel

00

Blank

Diluent

H20

x After

Mixing Speed

40

Analysis Sequence

H2O Vol (μL)

10

Page 10 of 10

Catálogo: K224

FOSFATASE ALCALINA IFCC

Nº de testes: 500

Versão: 23/10/2020

Observações: Para garantir o melhor desempenho nas análises, a Bioclin recomenda o uso do Biocal - K072 para Calibração e o uso do Biocontrol N – K073 e do Biocontrol P - K074 para Controle de Qualidade.

Test Name	Code	Local Code	Channel
FOSFATASE ALCALINA	FALC	FALC	00

Correlation Factor

Slope 1.00000

Intercept 0.000

Reaction Direction

X Reaction Direction Check

Reaction Direction Increase

Sample Limit Check

Sample Limit

Sample Limit Cycle 1

Definition

Calculation Type Kinetic

OD Deviation check

Linear regression

r2 restored SD

First point

First point threshold

Last point

SD factor

Antigen excess activation

Antigen excess limit (%)

Antigen excess point 1

Reaction Limit check

Reaction limit absorbance

Cycle 1

First Reading

Cycle 5

Last Reading

Cycle 20