

Test Name	Code	Local Code	Channel	<input checked="" type="checkbox"/> Enable	Release	Modified on
CREATININA AUT	CRE	CRE	00			19/03/2020
Characteristics			Pre-dilution		Result	
Sample Type	Serum   Plasma		Diluent name		Unit	mg/dL
Number of Reagents	1		Factor		Decimal Position	2
			Incubation time (in cycles)	1	<input type="checkbox"/>	Manual patient validation

## Result

Unit	mg/dL
------	-------

Decimal Position

Manual patient validation

## Correlation

Slope	Intercept
1.0	- 0.3

## Delta Check

Delta Check Validity	Absolut Variation	Relative Variation
365		

**Reference Range (μLmol/L)**

High check      1,2      1,0

**Rerum Range (μLmol/L)**

	Man/Default	Woman	Child
Low check			

☐ High check

Catálogo: K222

CREATININA AUTOMAÇÃO

Nº de testes: 2000

Versão: 19/03/2020

Observações: Para o formato Birreagente recomenda-se a utilização do calibrador multiparâmetro de bioquímica Biocal – K072. Para avaliar a precisão e a exatidão das dosagens, recomendamos o uso dos soros controle Biocontrol N – K073 e Biocontrol P – K074.

Test Name	Code	Local Code	Channel
CREATININA AUT	CRE	CRE	00

Pre-dilution

Type	Calibrador Diluent

Factor 1	Factor 2	Factor 3	Factor 4
Factor 5	Factor 6	Factor 7	Factor 8

Calibration

Calibration mode	Slope Average
Level	1
Calibration Factor	
Run(s)	2

Validity

<input checked="" type="checkbox"/> On request	<input type="checkbox"/> Time Validity
Interval	Time Unit
0	Days

Factor calibration

Low limit check	
High limit check	
Relative limit check	

Dev_Rep (%)	<input checked="" type="checkbox"/>	Dev_C (%)	<input checked="" type="checkbox"/>
-------------	-------------------------------------	-----------	-------------------------------------

Calibrator Used
Biocal

Checks

<input checked="" type="checkbox"/>	Reagent Limit Absorbance Check
-------------------------------------	--------------------------------

Reagent Range Low	-3.00000
-------------------	----------

Reagent Range High	3.00000
--------------------	---------

<input checked="" type="checkbox"/>	Reagent Blank Limit Absorbance Check
-------------------------------------	--------------------------------------

Blank Range - Low limit	-3.00000
-------------------------	----------

Blank Range - High limit	3.00000
--------------------------	---------

Control required

<input checked="" type="checkbox"/> Control 1	Biocontrol N
<input checked="" type="checkbox"/> Control 2	Biocontrol P
<input type="checkbox"/> Control 3	

Validity Backup

<input type="checkbox"/>	Backup time frame without calibration required
--------------------------	--

Interval	Time Unit
0	Days

Catálogo: K222

CREATININA AUTOMAÇÃO

Nº de testes: 2000

Versão: 19/03/2020

Observações: Para o formato Birreagente recomenda-se a utilização do calibrador multiparâmetro de bioquímica Biocal – K072. Para avaliar a precisão e a exatidão das dosagens, recomendamos o uso dos soros controle Biocontrol N – K073 e Biocontrol P – K074.

Test Name	Code	Local Code	Channel
CREATININA AUT	CRE	CRE	00

Cleaner

Cleaner Solution

Before

x

After

Wavelength (nm)

Primary Wavelength

505

Secondary Wavelength

Blank

Reagent Blank

Diluent

H20

Mixing Speed

40

Analysis Sequence

Cycle	Reagent	Volume (µL)	Sample Neede	Volume (µL)	H2O Vol (µL)
1	R1	120	Sample	15	10
10			R2	30	10

Catálogo: K222

CREATININA AUTOMAÇÃO

Nº de testes: 2000

Versão: 19/03/2020

**Observações:** Para o formato Birreagente recomenda-se a utilização do calibrador multiparâmetro de bioquímica Biocal – K072.  
Para avaliar a precisão e a exatidão das dosagens, recomendamos o uso dos soros controle Biocontrol N – K073 e Biocontrol P – K074.

Test Name

Code

Local Code

Channel

CREATININA AUT

CRE

CRE

00

### Correlation Factor

Slope 1.00000

Intercept -0.300

### Reaction Direction

X Reaction Direction Check

Reaction Direction Increase

### Sample Limit Check

Sample Limit

Sample Limit Cycle 1

### Definition

Calculation Type Kinetic

### OD Deviation check

Linear regression

r2 restored SD

First point

First point threshold

Last point

SD factor

### Antigen excess activation

Antigen excess limit (%)

Antigen excess point 1

### Reaction Limit check

Reaction limit absorbance

Cycle 1

### First Reading

Cycle 15

### Last Reading

Cycle 20