

Catálogo: K083

COLESTEROL MONOREAGENTE

Nº de testes: 1000

Versão: 17/01/2017

**Observações:** Para garantir o melhor desempenho nas análises, a Bioclin recomenda o uso do Biocal - K072 para Calibração e o uso do Biocontrol N – K073 e do Biocontrol P - K074 para Controle de Qualidade.

Test Name	Code	Local Code	Channel	<input checked="" type="checkbox"/> Enable	Release	Modified on
COLESTEROL	COL	COL	00			17/01/2017

**Characteristics**

Sample Type

Number of Reagents

**Pre-dilution**

Diluent name

Factor

Incubation time (in cycles)

**Result**

Unit

Decimal Position

☐ Manual patient validation

**Reagent**

Reagent short name Reagent number

On board stability (days)

☐ Cassete

**Linearity (µmol/L)**

Low Limit High Limit

**Correlation**

Slope Intercept

**Delta Check**

Delta Check Validity Absolut Variation Relative Variation

**Reference Range (µLmol/L)**

	Man/Default	Woman	Child
<input checked="" type="checkbox"/> Low check	<input type="text" value="0"/>	<input type="text" value="0"/>	<input type="text"/>
<input checked="" type="checkbox"/> High check	<input type="text" value="200"/>	<input type="text" value="200"/>	<input type="text"/>

**Automatic Rerun**

☒ Post Dilution

☐ Post Concentration

**Rerun Range (µLmol/L)**

	Man/Default	Woman	Child
<input type="checkbox"/> Low check	<input type="text"/>	<input type="text"/>	<input type="text"/>
<input type="checkbox"/> High check	<input type="text"/>	<input type="text"/>	<input type="text"/>

Catálogo: K083

COLESTEROL MONOREAGENTE

Nº de testes: 1000

Versão: 17/01/2017

**Observações:** Para garantir o melhor desempenho nas análises, a Bioclin recomenda o uso do Biocal - K072 para Calibração e o uso do Biocontrol N – K073 e do Biocontrol P - K074 para Controle de Qualidade.

Test Name	Code	Local Code	Channel
<input type="text" value="COLESTEROL"/>	<input type="text" value="COL"/>	<input type="text" value="COL"/>	<input type="text" value="00"/>

**Checks**

☒ Reagent Limit Absorbance Check

Reagent Range Low

Reagent Range High

☒ Reagent Blank Limit Absorbance Check

Blank Range - Low limit

Blank Range - High limit

**Pre-dilution**

Type	Calibrador Diluent
<input type="text"/>	<input type="text"/>

Factor 1	Factor 2	Factor 3	Factor 4
<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
Factor 5	Factor 6	Factor 7	Factor 8
<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>

**Calibration**

Calibration mode

Level

Calibration Factor

Run(s)

**Validity**

☒ On request ☐ Time Validity

Interval	Time Unit
<input type="text" value="0"/>	<input type="text" value="Days"/>

**Factor calibration**

Low limit check

High limit check

Relative limit check

**Control required**

☒ Control 1

☒ Control 2

☐ Control 3

**Validity Backup**

☐ Backup time frame without calibration required

Calibrator Used

Interval	Time Unit
<input type="text" value="0"/>	<input type="text" value="Days"/>



Catálogo: K083

COLESTEROL MONOREAGENTE

Nº de testes: 1000

Versão: 17/01/2017

**Observações:** Para garantir o melhor desempenho nas análises, a Bioclin recomenda o uso do Biocal - K072 para Calibração e o uso do Biocontrol N – K073 e do Biocontrol P - K074 para Controle de Qualidade.

Test Name	Code	Local Code	Channel
<input type="text" value="COLESTEROL"/>	<input type="text" value="COL"/>	<input type="text" value="COL"/>	<input type="text" value="00"/>

**Correlation Factor**

Slope

Intercept

**Reaction Direction**

☒ Reaction Direction Check

Reaction Direction

**Sample Limit Check**

Sample Limit

Sample Limit Cycle

**Definition**

Calculation Type

**OD Deviation check**

☐ Linear regression

r2 restored  SD

☐ First point

First point threshold

☐ Last point

SD factor

**Antigen excess activation**

Antigen excess limit (%)

Antigen excess point

**Reaction Limit check**

Reaction limit absorbance

Cycle

**First Reading**

Cycle

**Last Reading**

Cycle