

Catálogo: K094

COLINESTERASE

Nº de Testes: 100

Versão: 05/04/2020

Cód.: 00

K094-1

Reagente Nº 1 – Tampão Reagente – 1 x 20 mL

Reagente Nº 2 – Substrato – 1 x 5 mL

Os reagentes são prontos para uso.

Chemistry Parameters¹

Method No.	#	Method Name	COLI	R1	Tampão	Reagent (µl)	200	Water (µl)	0
Print Name	COLINESTERASE			R2	Substrato	50	0		
Sample Type	Serum			Diluent	H ₂ O				
Unit	U/L			Sample Ppt Wash	Disable				
Assay Type	Rate								

Measuring Points	Start	End
1	29	38
2	Disable	0

Stirring Speed
R1 Middle R2 Middle

Wave Length
Pri 415 Sec -

Normal Range

Nº	Normal Range Name	Min	Máx
1	Male-G1	0.00	0.00
2	Male-G2	0.00	0.00
3	Male-G3	0.00	0.00
4	Female-G1	0.00	0.00

Normal	Sample Volume (µl)
Low	Normal
High	
Dilution	0.0 < 5 < 0.0
Rerun (High/Prozone)	
Dilution	0.0 < 5 < 0.0
Rerun (Low)	
Dilution	0.0 < 5 < 0.0

Technical Range

Conc.	50	-	25000
(mAbs./10)	-999999	-	999999
Prev. Result Comparison (%)	999999.00	0.00	%
Abnormal Range			
(Conc.)	0.00	-	999999.00
Panic Range			
(Conc.)	0.00	-	999999.00

Decimal Point	0	Profile SI	Disable
---------------	---	------------	---------

Para garantir o melhor desempenho das análises, a Bioclin recomenda o uso do Kit de [Biocal – K072](#) para Calibração e o uso do [Biocontrol N – K073](#) e [Biocontrol P – K074](#) para Controle Interno de Qualidade.